



Agilent Technologies

Setting Network Parameters 6890N Gas Chromatograph

A05294

This document is believed to be accurate and up-to-date. However, Agilent Technologies, Inc. cannot assume responsibility for the use of this material.
The information contained herein is intended for use by informed individuals who can and must determine its fitness for their purpose.

Configuring IP Addresses from the Keyboards

6890N GC Configuration

The Agilent 6890N GC communicates with a ChemStation via the LAN. It does not require the BootP Server software. This procedure applies only when adding a 6890N GC. On a multi-instrument system where BootP Server is already installed to provide IP addresses to instruments, BootP should remain installed and enabled. 6890N/ChemStation LAN communication operates whether BootP is present or not.

To enter the IP address into a 6890N GC:

- 1 On the 6890N keyboard, press Options. Scroll to Communication and press Enter. The following screen appears:**

```
COMMUNICATION SETPTS
---- LAN ----
IP: 000.000.000.000 <
GW: 000.000.000.000
SM: 000.000.000.000
Enable DHCP OFF
---- RS-232 ----
```

- 2 Enter the IP address for your 6890N GC. Enter the numbers separated by dots and press Enter. The GC displays a message instructing you to power cycle the instrument. Do NOT power cycle yet.**
- 3 Press Clear. Scroll to GW. Enter the Gateway number and press Enter. The GC displays a message instructing you to power cycle the instrument. Do NOT power cycle yet.**
- 4 Press Clear. Scroll to SM and press Mode/Type. Scroll to the appropriate subnet mask from the list of modes and press Enter. The GC displays a message instructing you to power cycle the instrument.**
- 5 Power cycle the instrument to apply the LAN setpoints to the card.**
- 6 Press Options. Scroll to Communications and press Enter. Confirm that the correct setpoints are present.**