

Configuring A GPIB (HPIB) GC Instrument Module Using The ChemStation Configuration Editor Windows® 95 / Windows NT 4.0®

Before an ONLINE instrument session can be utilized, the instrument module(s) must be configured for the ChemStation software via the Configuration Editor application. The following is a procedure for performing this process:

I. Start The HP Configuration Editor Application

1. Select the **Start** button on the Windows Task Bar
2. From the **Start** menu select **Programs**
3. From the **Programs** menu select **HP ChemStations**
4. From the **HP ChemStations** menu select **HP Configuration Editor**

II. Configure the Overall Instrument Type

1. Select the Instrument # (1,2,3,4) panel in the middle of the Configuration Editor panel
2. Pull down the **Configure** menu and select **Instruments**

*Note: If the instrument session is already configured and only module addition is required, skip this step (II) by immediately selecting the **OK** button on this panel*

3. From the instrument type list, select the instrument type you wish to configure

- A. 5890 GC
- B. 6890 GC
- C. 35900 A/D

- 35900C – External A/D
- 35900E – External A/D
- 35900D – Internal A/D Card

This document is believed to be accurate and up-to-date. However, Agilent Technologies, Inc. cannot assume responsibility for the use of this material. The information contained herein is intended for use by informed individuals who can and must determine its fitness for their purpose.

MS-DOS®, Windows® and Windows NT® are registered trademarks of Microsoft Corp.

4. Enter an instrument name into the **Instrument Name** field
5. Select the **OK** button

III. Adding / Configuring Modules On The GC System

1. 5890 GC System
 - A. Enter the GPIB (HPIB) address for the 5890 mainframe into the **Instrument HPIB** address field (factory default GPIB (HPIB) address is 15)
 - B. Select additional module from the **Options** list, enter the GPIB (HPIB) address for the module being added, and select the **Add** button (see Appendix for module GPIB (HPIB) address listing)
 - C. Repeat step "B" above for additional modules
 - D. Select the **OK** button
2. 6890 GC System
 - A. Enter the GPIB (HPIB) address for the 6890 mainframe into the **Link Options HPIB** address field (factory default GPIB (HPIB) address is 0)
 - B. Select additional module from the **Options** list, enter the GPIB (HPIB) address for the module being added, and select the **Add** button (see Appendix for module GPIB (HPIB) address listing)
 - C. Repeat step "B" above for additional modules
 - D. Select the **OK** button
3. 35900 A/D
 - A. 35900C / 35900E
 - 1) Enter the GPIB (HPIB) address for the 35900 into the **Interface** address field (factory default GPIB (HPIB) address is 13)
 - 2) Select acquisition channels **A and/or B** (must select at least 1)
 - 3) Select the **OK** button
 - B. 35900D
 1. Enter the Select Code for the 35900 into the **Interface** address field (factory default Select Code is 6)
 - 2) Select acquisition channels **A and/or B** (must select at least 1)
 - 3) Select the **OK** button

IV. Repeat Steps II and III Above For Additional Instrument

V. Save The Configuration And Exit The Application

1. Pull down the **File** menu and select **File**
2. Pull down the **File** menu and select **Exit**

Appendix

GC Module / Default (Factory) GPIB (HPIB) Address Listing

5890 GC

5890 GC Instrument (Mainframe)	15
Autosampler	8

6890 GC

6890 GC Instrument (Mainframe)	0
--------------------------------	---

35900 A/D

35900C	13
35900E	13
35900D	Select Code 6