



Agilent Technologies



Accurate quantitation of regulated mycotoxins by UHPLC-MS/MS



Dr. Thomas Glauner

Senior LCMS Applications Chemist
Agilent Technologies

In cooperation with:

Elisabeth Varga*

Rudolf Krska*

Rainer Schuhmacher*

Franz Berthiller*

*Center for Analytical Chemistry, Department of Agrobiotechnology (IFA-Tulln),
University of Natural Resources and Life Sciences Vienna, Austria



Agilent Technologies

Confidentiality Label
October 13, 2011

Outline

- Mycotoxins
- Stable isotope dilution assay (SIDA) for the mycotoxins regulated in the European Union
 - Reasons
 - Approaches
 - Sample Preparation and Method
 - Results
- Summary and conclusions



Mycotoxins

Background

myces (Greek) = fungus

toxicum (Latin) = toxic

} = Mycotoxin

- low molecular weight, toxic, secondary metabolites of fungi
- produced by e.g.:
 - Fusarium sp., Aspergillus sp., Penicillium sp.
- toxicity:
 - acute toxic, carcinogenic, mutagenic, teratogenic, estrogenic and immunotoxic effects



provided by Prof. Bürstmayr

Mycotoxins

Why are they an issue?


- >25% of all agricultural commodities are contaminated with mycotoxins
- annual losses of several hundred million tons of food worldwide
- annual economical losses: 1 billion USD (US only)
- 100+ countries have regulations for the control of mycotoxins in food and feed



Mycotoxins

Relevance for food control

- Notifications concerning mycotoxins
 - (RASFF-Annual reports 2002-2008)



Food contaminants online

FERA FC24 database access

- Free access to the RASFF alerts and notifications via FC24 database
- Registration through Agilent website
(<http://www.chem.agilent.com/en-US/solutions/foodtesting-agriculture/Pages/default.aspx>)

Agilent Technologies

Home Page | Search by Commodity | Log out

Search the Contaminants Database

Food type: ☒ Food ☐ Feed

Choose From
Commodity (ordered by group): **Bivalve molluscs**

OR
Commodity (ordered A - Z): (Select Commodity)

Food Contaminant: (All Food Contaminants)

Country of Origin: (All Countries)

RASFF Notification Type: ☒ Alert Notifications ☐ Information Notifications ☐ Border Rejections ☐ All Notifications

Tel: +44 (0)1904 462613 fc24@fera.gov.uk
Built by [KM@Fera](#)

[Legal Notices](#) | [Terms & Conditions](#)
[Acknowledgement](#) | [Disclaimer](#)
© Crown Copyright Fera 2010


Accurate quantitation of mycotoxins

Reasons

- European Commission Regulation (EC) No 1881/2006
 - set maximum residue levels (MRL) for mycotoxins
- Single target versus multi-target methods

BUT:


- Electrospray ionisation (ESI)
 - matrix effects hamper accurate mass spectrometric quantification
- Quantification of regulated mycotoxins at a very high degree of accuracy is required



Matrix effects in ESI-MS and quantitation

Approaches

- Dilution of the sample
 - method less sensitive
- Matrix matched calibration
 - tedious
 - differences within one commodity not compensated
- Standard addition to each sample
 - more runs
 - more costs (time and standards)
- Internal calibration
 - similar compounds (ZAN for ZEN)
 - deuterium or ^{13}C -labelled compounds
 - until this year: only single analyte or group analyte IS-addition
 - usually associated with rather high costs




Stable Isotope Dilution Assay (SIDA)

Aims

- Development of a method fulfilling:
 - covering all regulated mycotoxins in solid food matrices
 - providing best possible accuracy
 - easy to handle
 - cost effective
- Stable isotope dilution assay (SIDA) for LC-MS/MS
 - 11 mycotoxins
 - ^{13}C -labelled compounds as internal standards
 - validation of the method for maize




Multiple Reaction Monitoring Principles



Sample preparation

Universal extraction procedure



Stable Isotope Dilution Assay (SIDA)

HPLC method



Agilent 1290 Infinity LC system consisting of:

- binary pump
- wellplate sampler
- column compartment
- diode array detector (not used)

HPLC method

Separation column: ZORBAX Eclipse Plus C-18 RRHD column,
100 x 2.1 mm, 1.8 μ m @ 30°C

Mobile phase: A: 5 mM HCOONH₄ + 0.1% formic acid
B: methanol + 5 mM HCOONH₄ + 0.1% formic acid

Flow: 0.35 ml/min

Gradient:	0.00 min	30 % B
	0.50 min	30 % B
	8.00 min	100 % B
	9.50 min	100 % B
	9.60 min	30 % B

Inj.Vol.: 3 μ l

Stable Isotope Dilution Assay (SIDA)

MS method

Spray chamber conditions:

Gas temp.:	140°C
Dry gas:	16 l/min
Nebulizer:	25 psi
Sheath gas temp:	350°C
Sheath gas flow:	11 l/min



	Positive	Negative
CapVoltage:	4000 V	3000 V
Nozzle voltage	0 V	0 V


Automatic setup of MRM tables based on selected cycle time, retention times and retention time windows for the individual compounds

• Cycle time	400 ms
• Interscan delay	3.5 ms
• Total No. of MRMs	33
• Maximum No. Of concurrent MRMs	12
• Minimum Dwell time	39.8 ms
• Maximum Dwell time	196.5 ms



Dynamic MRM functionality


Comparison of MRM and DMRM



- 2 x shorter cycle times supports narrow chromatographic peaks, more analytes or longer dwell per analyte.

Dynamic MRM functionality

DMRM simulation




Concurrent
Compounds



Stable Isotope Dilution Assay (SIDA)

Chromatogram


- due to same MRM transitions baseline separation required for:
 - fumonisin B2 and B3
 - aflatoxin G1 and ^{13}C -aflatoxin B1
 - aflatoxin G2 and ^{13}C -aflatoxin B2




Internal calibration in solvent

Aflatoxin B1


- Challenging compound due to low MRLs
 - 0.1 µg/kg in processed cereal based baby food
 - 2 to 12 µg/kg in nuts and cereals




Overlay of 4 individual calibrations
acquired within 45 hour worklist.



0.0075 ng/ml



0.0225 ng/ml



0.075 ng/ml




Internal calibration in solvent


Ochratoxin A

- Challenging compound due to low MRLs
 - 0.5 µg/kg in processed cereal based baby food
 - 3.0 / 5.0 µg/kg in processed / unprocessed cereals
 - 10.0 µg/kg in dried vine fruit

Ochratoxin A - 9 Levels, 9 Levels Used, 36 Points, 36 Points Used, 0 QCs




ISTD-10000 [Ochratoxin A]




0.0077 ng/ml

ISTD-3000 [Ochratoxin A]



0.0230 ng/ml

ISTD-1000 [Ochratoxin A]



0.0765 ng/ml

Overlay of 4 individual calibrations
acquired within 45 hour worklist.



Agilent Technologies

Confidentiality Label
October 13, 2011

Validation of SIDA method

Experimental setup and results



- Linear range (external calibration in solvent)
 - 4 orders of magnitude for all toxins, 5 orders for Aflatoxins, T-2, and ZEN
- Costs
 - Additional price per IS per sample is between 0.01 to 1.40 €
 - Price for all 11 IS per sample < 2.00 €
- Full validation for maize
 - Maize: - matrix for which most mycotoxins are regulated
 - known for matrix effects and matrix interferences
 - more costs (time and standards)
 - Spiking with native mycotoxins before extraction
 - Six concentration levels with 3 replicates
 - Spiking with ¹³C-labelled mycotoxins before analysis to compensate matrix effects in ESI
 - No sample clean-up



Validation of SIDA method in maize

Extraction of spiked blank maize and reference materials

- Blank maize sample spiked with native mycotoxins before extraction
 - includes 10-fold dilution of matrix in the final extract due to extraction procedure



Validation of SIDA method

Results – Sample preparation

■ Extraction efficiency

- Determined by spiking of blank samples before extraction
- First extraction: 80% acetonitrile content (60 min)
 - recovery between 80 and 110% except for FB1 and FB2
- Second extraction: 20% acetonitrile content (30 min)
 - improved extraction recovery for FB1 and FB2 to approx. 90%

■ Matrix effects

- Signal suppression
 - 50 to 60% aflatoxins
 - 50% DON
- Signal enhancement
 - Fumonisin, HT-2, T2, OTA
- Effectively compensated by ISTD



Validation of SIDA method

Results for maize

Analyte	LOQ in µg/kg	R_A^{**} in % $\pm$ RSD in %	
Aflatoxin B ₁	0.04	105	6
Aflatoxin B ₂	0.04	100	4
Aflatoxin G ₁	0.05	101	5
Aflatoxin G ₂	0.24	101	8
Deoxynivalenol	2.5	99	9
HT-2 toxin	2.0	98	7
T-2 toxin	0.17	99	6
Ochratoxin A	0.23	93	7
Zearalenone	0.97	103	11
Fumonisin B ₁	2.5	101	10
Fumonisin B ₂	0.64	88	7


** average for triplicate samples and 6 spiking levels



DMRM database for mycotoxins

Customize your mycotoxin method

- Multi-mycotoxin method for 242 mycotoxins and other fungal metabolites has been developed
 - Validated for different nuts
 - Transitions are shortly available as DMRM database



Summary

- UHPLC-MS/MS method
 - Improved chromatographic resolution
- Multiple extraction steps
 - Enhancement of extraction efficiency especially for fumonisins
- Dynamic MRM with fast polarity switching
 - Most abundant ionization mode and maximized dwell times within a single run
- Addition of internal standards after extraction
 - Compensation for matrix effects
 - Minimized costs
- Apparent recoveries of 88 to 105% for all mycotoxins
 - Evaluated by extraction of spiked maize samples
 - Validated by correct quantitation of 18 reference materials covering all toxin groups
- Sensitivity suitable for MRLs
 - Improved sensitivity of G6490 allows to omit sample concentration resulting in easier handling and improved robustness



Questions



Validation of SIDA method

Method characteristics and requirements

Analytes	Linear range ng/mL	LOQs (maize) µg/kg	MRLs (EC Reg. No 1881/2006)	Commodities
Aflatoxin B ₁	0.0075 - 74.6	0.04	0.1	processed cereal-based baby food
			2.0-12 sum of aflatoxins: 4.0-15.0	nuts and cereals
Deoxynivalenol	0.23 - 225	2.5	200 500 - 1750	processed cereal-based baby food processed / unprocessed cereals, bread, pasta, breakfast cereals
Fumonisin B ₁	0.075 - 249	2.5	200 800 1000 / 4000	processed maize-based baby food maize-based breakfast cereals maize / unprocessed maize
Fumonisin B ₂	0.075 - 251	0.64		
HT-2 toxin	0.2 - 202	2.0	implementation of MRLs is expected in the near future	unprocessed cereals and cereal products
T-2 toxin	0.023 - 75.4	0.17		
Ochratoxin A	0.023 – 23.0	0.23	0.5 3.0 / 5.0 10.0 15 / 20 / 80	processed cereal-based baby food processed / unprocessed cereals dried vine fruit spices / liquorice root / extract
Zearalenone	0.076 - 252	1.0	20 50 75-350	processed cereal-based baby food bread, biscuits, breakfast cereals processed / unprocessed cereals

