

Factors Influencing the T_m of Nucleic Acid Hybrids

Factor	Influence on T_m
Ionic strength	T_m increases 16.6°C for each one-fold increase in monovalent cations, between 0.01 and 0.40 M NaCl.
Base composition	AT base pairs are less stable than GC base pairs in aqueous solutions containing NaCl. The difference is negligible in tetramethyl ammonium chloride.
Destabilizing agents	Each 1% formamide reduces the T_m by $\sim 0.6^{\circ}\text{C}$. 6M urea reduces the T_m by about 30°C .
Mismatched base pairs	T_m is reduced by 10°C for each 1% of mismatch.
Duplex length	Negligible effect with probes > 500 bp.