



WELCOME

This event will begin soon

Music is now
playing.

If you are unable to
hear the audio,
please take a
moment to test your
system by clicking
HELP in the
resource box.

Interested in conducting your own webinar? Email webinars@bnpmmedia.com

Determination of Multi-Pesticide Residues in Dried Tea Samples using an Optimized Extraction / Clean-up Regime and the Agilent 7000 Series Tandem Quadrupole GC/MS/MS System

**Prof. Jana Hajslova, Dr. Tomas Cajka, Dr. Jana Pulkrabova,
Veronika Bachanova, Lucie Drabova, Kamila Kalachova**

**Department of Food Analysis and Nutrition
Institute of Chemical Technology
Prague, Czech Republic**

in collaboration with
Chris Sandy
EMEA GC-MS Food Applications Chemist
Agilent Technologies UK



How to Use Your Console

RESOURCES

Enlarge Slides
Download Slides

Invite a Friend
Email BNP Media


HELP

-Enlarge Slides for better viewing
-HELP link
-Download & print slides
-Invite a friend to join
-Email webinars@bnpmedia.com

Submit a Question

Submit a Question

Disable Pop-Up Blocker



Archive

www.foodsafetysummit.com

Add to your favorites

Determination of Multi-Pesticide Residues in Dried Tea Samples using an Optimized Extraction / Clean-up Regime and the Agilent 7000 Series Tandem Quadrupole GC/MS/MS System

**Prof. Jana Hajslova, Dr. Tomas Cajka, Dr. Jana Pulkrabova,
Veronika Bachanova, Lucie Drabova, Kamila Kalachova**

**Department of Food Analysis and Nutrition
Institute of Chemical Technology
Prague, Czech Republic**

in collaboration with
Chris Sandy
EMEA GC-MS Food Applications Chemist
Agilent Technologies UK






Presentation Overview

1.



Introduction

2.




Optimization of sample prep / method validation

3.



GC/MS/MS Analysis

4.



Results

5.



Summary



1. Introduction

- Collaboration at ICT started March 2011
- Key objectives
 - Optimization of sample prep
 - Development of robust / rapid GC method > 100 target analytes
- All analytical development performed at ICT
- Analysis performed at ICT and in UK




Tea – *Camelia sinensis*




WW Tea Production by Country - 2008

Tea – *Camelia sinensis*




Production processes for Tea




The many varieties of Tea...

- White Tea: Wilted and un-oxidized
- Yellow Tea: Unwilted and un-oxidized, but allowed to yellow
- Green Tea: Unwilted and un-oxidized
- Oolong Tea: Wilted, bruised and partially oxidized
- Black Tea: Wilted, sometimes crushed and fully oxidized
- Post-fermented Tea: Green Tea that has been allowed to ferment / compost

- Flavoured Tea: Fruit, herbal etc




The many varieties of Tea...



Other teas are available

Common Pests in Tea Plantations



Major tea producing countries



Pest	Common Name
Insects	Thrips
	Slug caterpillars
	Bunch caterpillars
	Leaf roller
	Looper
	Tea mosquito
	Aphids
	Jassids
	Helopeltis
	Flushworm
Mites	Purple, pink and red
Fungi	Blister blight
	Black rot
	Red rust



Agro-chemicals used on Tea Plantations

Chemical Class	Example Chemicals	Pest / Disease
Insecticides	OCPs, OPs, Pyrethroids	Leafhoppers, plant bugs, aphids, thrips, leaf eating / leaf folding caterpillars, beetles etc
Acaricides	Amitraz, Azinphos ethyl, Diazinon, Naled	Mites, Ticks etc
Fungicides	Carbendazim, Hexaconazole, Propiconazole, Tridemorph	Fungal diseases – Blister blight, Grey blight etc



EU Pesticide MRL Legislation

REGULATION (EC) NO 396/2005 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL

of 23 February 2005

on maximum residue levels of pesticides in or on food and feed of plant and animal origin and amending Council Directive 91/414/EEC

(Text with EEA relevance)

Recital 5

(5) One of the most common methods of protecting plants and plant products from the effects of harmful organisms is the use of active substances in plant protection products. However, a possible consequence of their use may be the presence of residues in the treated products, in animals feeding on those products and in honey produced by bees exposed to those substances. According to Council Directive 91/414/EEC of 15 July 1991 concerning the placing of plant protection products on the market (1), **public health should be given priority over the interests of crop protection, thus it is necessary to ensure that such residues should not be present at levels presenting an unacceptable risk to humans and, where relevant, to animals.** MRLs should be set at the lowest achievable level consistent with good agricultural practice for each pesticide with a view to protecting vulnerable groups such as children and the unborn.

➤ Requirement for strict residue controls to safeguard public health

FC24 Food Commodities and Contaminants Database

FC24 Food Commodities and Contaminants Database. Agilent offers our customers **free access** to the FERA FC24 database, which pulls together relevant food contaminants and residues information from across the EU. In a matter of seconds, you can search the database and have the information you need at your fingertips.

<https://secure.fera.defra.gov.uk/foodcontaminants/agilent/index.cfm>



FC24 Food Commodities and Contaminants Database



Agilent Technologies



[Home Page](#)

[Search by Commodity](#)

[Log out](#)

EU & CODEX MRL Database (under Regulation 396/2005/EC)

Active Substance

Food Type

Food

Annex

Current / Proposed MRLs

Also Display Residue Definition: Yes ☐ No ☒

The absence of an active substance or food from these lists does not imply any error in the data as this dataset is complete. Absence of data indicate that the default MRL applies (0.01*).



FC24 Food Commodities and Contaminants Database



Agilent Technologies



[Home Page](#) | [Search by Commodity](#) | [Log out](#)

LIAISON EU & CODEX MRL Database (under Regulation 396/2005/EC) Search Results for Fenvalerate and Esfenvalerate (Sum of RR & SS isomers)

There is 1 request matching your criteria.

Active Substance	MRL mg/kg	Previous MRL (Date of Change)	Food type	Food	Comments	Codex MRL (derivation)	Annex
Fenvalerate and Esfenvalerate (Sum of RR & SS isomers)	0.05*	(1 September 2008)	Tea (dried leaves and stalks, fermented or otherwise of Camellia sinensis)	Tea		No Codex MRL has been set	II

* Level at or about the limit of determination (LOD)

(1) The code number is introduced by this Annex and is intended to set a classification under this and other related Annexes of Regulation (EC) No 396/2005.

(2) The scientific name of the items listed in the column "Examples of individual products within the groups to which the MRLs apply", where possible and relevant, is mentioned. As much as possible the International System of Nomenclature is followed.

(9) Part of the products to which the MRL apply.


MRLs with a green background indicate that the active has approval (in at least one product, either on-label or off-label) on the listed crop within the UK. Such approved uses at the LOD (i.e. those marked *) indicate that no detectable residue should occur under GAP

MRLs with no background colour indicate that no UK uses are currently approved. However, additional presence of a CODEX MRL would indicate approval somewhere in the world and residues above the MRL would not be expected under the appropriate GAP in those countries

[New Search](#)

Sources





2. Optimization of Sample prep / Method Validation




Tea matrix composition

Main constituents of tea leaves

Component	Content (% dry weight)
Polyphenols	25–35
flavanols	80% of total polyphenols
Sacharides	25
polysacharides	14–22
Proteins	15
Minerals	5
Free aminoacids	4
Chlorophyll	0.5
CAFFEINE	2.5–5.5

CAFFEINE



Content in cup of tea

Black Tea: 23–110 mg

Oolong Tea: 12–55 mg

Green Tea: 8–36 mg

White Tea: 6–25 mg



**40% dry matter soluble in water
(fermented tea)**



SAMPLE PREPARATION PROCEDURE



(1)

**Weigh 2 g of
homogenized
tea sample**



(2)

**Add 10 mL of
distilled water
Shake for 30 s**



(3)

**Matrix swelling
(hydration) for
30 min**



(4)

**Add 10 mL of
MeCN
Shake for 1 min**

SAMPLE PREPARATION PROCEDURE



(5)

Add 1 g NaCl



(6)

**Add 4 g MgSO₄
Shake for 1 min**



(7)

**Add internal
standard (TPP)**



(8)

**Centrifuge 5
min at 10,000
rpm**

SAMPLE PREPARATION PROCEDURE




(9)

Add 1 mL of MeCN extract into a 15-mL tube containing 1 mL of hexane and 5 mL of 20% (w/w) NaCl solution, shake for 1 min



(10)

Centrifuge 1 min at 10,000 rpm



(11)

Take a part of the hexane layer into a vial

Bond Elut – QuEChERS Consumables



(5)
Add 1 g NaCl




(6)
Add 4 g MgSO_4
Shake for 1 min




Agilent Bond Elut : P/No. 5982-5550
50x 50mL tubes + 50x salt packets
[4g MgSO_4 + 1g NaCl]




SAMPLE PREPARATION PROCEDURE



SAMPLE PREPARATION PROCEDURE



Easily prepare 6 samples / hour



■ Acetonitrile ■ Water ■ Hexane
■ Acetonitrile / 20% Aq. NaCl

GC/MS/MS Methodology (Summary)




GC–MS/MS analysis

- Injection: MMI, Cold splitless (2 μ L)
- Column: HP-5MSUI
(15 m $\times$ 0.25 mm ID $\times$ 0.25 μ m)
- Retention time locking
- Post-column, post-run back flush
- 2–4 MRM transitions per analyte

GC run: 20 min
 Backflush: 2 min
 MMI/GC cooling: 6 min
 Pre-injection steps : 1 min

**TOTAL : 29 min
per injection**



OPTIMIZATION OF SAMPLE PREPARATION

- **(1) Purification of crude MeCN extract**
 - dSPE, LLE
- **(2) Checking of sample preparation efficiency**
 - Gravimetric determination of co-extracts and determination of caffeine (dominating co-extract) using DART–MS
- **(3) LLE with hexane and NaCl solution**
 - Recovery of target analytes

OPTIMIZATION OF SAMPLE PREPARATION

- **(4) Comparison of various QuEChERS methods**
 - Original, buffered, citrate
- **(5) Matrix swelling prior to extraction**
 - QuEChERS (with water addition) vs. MeCN extraction (without water addition)
- **(6) Long term stability of analytes signal**
 - Repeated injections
- **(7) Method validation**
 - Recovery, repeatability, LOQ

(1) Purification of crude MeCN extract

(2) Checking of sample preparation efficiency

- Tested purification approaches
 - (A) dSPE followed by LLE
 - (B) dSPE without LLE
 - (C) only LLE (dSPE omitted)

(1) Purification of crude MeCN extract


(2) Checking of sample preparation efficiency

- Sorbents tested
 - 50 mg PSA + 150 mg MgSO_4 per 1 mL
 - 50 mg PSA + 50 mg CaCl_2 per 1 mL
 - 10 mg C18 + 50 mg PSA + 150 MgSO_4 per 1 mL
 - 10 mg C18 + 150 mg MgSO_4 per 1 mL
 - 10 mg C18 per 1 mL
- *Matrix content: 0.2 g/mL*

(1) Purification of crude MeCN extract

(2) Checking of sample preparation efficiency


- Content of caffeine in the examined extracts




(1) Purification of crude MeCN extract

(2) Checking of sample preparation efficiency

Content of caffeine in the examined extracts




- (1) Purification of crude MeCN extract
- (2) Checking of sample preparation efficiency
 - Gravimetrically determined co-extracts




(1) Purification of crude MeCN extract

(2) Checking of sample preparation efficiency

■ Gravimetrically determined co-extracts




- (1) Purification of crude MeCN extract**
- (2) Checking of sample preparation efficiency**
 - GC–MS analysis (MS¹ full scan) of the extracts




(1) Purification of crude MeCN extract

(2) Checking of sample preparation efficiency

- Influence of matrix co-extracts on RT, peak shape and signal intensity



100 ng/g
(20 ng/mL)




- (1) Purification of crude MeCN extract
- (2) Checking of sample preparation efficiency

- **CONCLUSIONS**

- **dSPE** (with $\text{MgSO}_4/\text{CaCl}_2$) is **not sufficient clean-up technique** for purification of tea MeCN extracts
- **LLE** (with hexane/20% NaCl (w/w) mixture) leads to **significant reduction of matrix co-extracts**
- *dSPE step with C18 sorbent (as proposed previously) can be omitted, thus, further increasing speed of sample preparation*

(3) LLE with hexane and NaCl solution

- 164 pesticides involved (MeCN extracts spiked before LLE, 1 mg/kg)
- MeCN extract : hexane : NaCl solution (1:1:5, v/v/v)
- In the case of 0% NaCl solution an emulsion formed; results obtained by LLE with MeCN extract : hexane (1:1, v/v)



(4) Comparison of various QuEChERS methods

Original QuEChERS

Buffered QuEChERS

Citrate QuEChERS

2 g sample + 10 mL distilled water

matrix swelling for 30 min

10 mL MeCN

10 mL 1% AcA MeCN

10 mL MeCN

Shake for 1 min

4 g MgSO₄, 1 g NaCl

**4 g MgSO₄,
1.7 g NaOAc · 3 H₂O**

**4 g MgSO₄, 1 g NaCl,
1 g Na₃Citr · 2 H₂O,
0.5 g Na₂HCitr · 1.5**

shake for 1 min, centrifuge 5 min at 10,000 rpm

1 mL of MeCN per 1 mL hexane and 5 mL 20% NaCl solution

shake for 1 min, centrifuge 1 min at 10,000 rpm


hexane layer into a vial for GC–MS/MS analysis

(4) Comparison of various QuEChERS methods

- **Content of caffeine in the examined extracts**
- **Gravimetrically determined co-extracts**
- **Recovery of target analytes**


(4) Comparison of various QuEChERS methods

■ Content of caffeine in the examined extracts




(4) Comparison of various QuEChERS methods

■ Gravimetrically determined co-extracts



(4) Comparison of various QuEChERS methods

Recovery of target analytes



pH of tea matrix: 5.1–5.2 (2 g + 10 mL distilled water) → no buffering required

(4) Comparison of various QuEChERS methods

■ CONCLUSIONS

- All three QuEChERS extraction techniques provide comparable results in terms of:
 - Recovery of target analytes
 - Amount of isolated co-extracts
- Original QuEChERS method used for validation experiments

(5) Matrix swelling prior the extraction

- Recently, studies dealing with the extraction of pesticide residues from tea using pure MeCN without previous addition of water for matrix hydration reported

G.F. Pang, *et al.*, *J. AOAC Int.* **94** (2011) 1253–1296

X.M. Xu, *et al.*, *J. Sep. Sci.* **34** (2011) 210–216

Z. Huang *et al.*, *J. Sep. Sci.* **32** (2009) 1294–1301

- **SANCO/12495/2011** strictly advises to use water for swelling the matrix before the extraction

‘To improve the extraction efficiency of low moisture containing commodities (cereals, dried fruits), it is recommended to add water to the samples before extraction is carried out.’

(5) Matrix swelling prior the extraction

- **In general, absence of water during the extraction of low moisture matrices results in cleaner extracts, but problem with lower recoveries of INCURRED pesticide residues**

M. E. Poulsen, H.B. Christensen, S.S. Herrmann, *Accred. Qual. Assur.* **14** (2009) 477–485

(5) Matrix swelling prior to extraction

Method A (with water)

2 g sample + 10 mL distilled water

matrix swelling for 30 min

10 mL MeCN

shake for 1 min

4 g MgSO₄, 1 g NaCl

shake for 1 min

centrifuge 5 min at 10,000 rpm

1 mL of MeCN per 1 mL hexane and 5 mL 20% NaCl solution,

shake for 1 min, centrifuge 1 min at 10,000 rpm

hexane layer into a vial for GC–MS/MS analysis

Method B (without water)

5 g sample + 15 mL MeCN

Turrax for 1 min

centrifuge 5 min at 10,000 rpm

MeCN into 25 mL vol. flask

repeat with 15 mL MeCN

combine MeCN extracts (25 mL)

(5) Matrix swelling prior the extraction

- **Content of caffeine in the examined extracts**
- **Gravimetrically determined co-extracts**
- **Comparison of recoveries: QuEChERS/LLE vs. extraction with pure MeCN**
- **Comparison of content of incurred residues obtained using both extraction techniques**

(5) Matrix swelling prior the extraction

■ Content of caffeine in the examined extracts

Method A (with water)

Caffeine: 3.8 mg/mL

Extract after LLE: 0.005 mg/mL

Method B (without water)

Crude extract: 0.13 mg/mL

Crude extract: 0.0001 mg/mL

■ Gravimetrically determined co-extracts

Method A (with water)

Crude extract: 38 mg/mL

Extract after LLE: 0.60 mg/mL


Method B (without water)

Crude extract: 1.3 mg/mL

Extract after LLE: 1.3 mg/mL

(5) Matrix swelling prior the extraction

■ Analysis of samples with incurred residues



Spikes (0.1 mg/kg)


Analyte	Recovery QuEChERS	Recovery Pure MeCN
Bifenthrin	78%	88%
Chlorpyrifos-Me	92%	85%
Cyhalothrin-λ	101%	86%
Endosulfan-α	77%	80%
Endosulfan-β	77%	80%
Endosulfan-SO ₄	90%	82%

(5) Matrix swelling prior the extraction

■ Analysis of samples with incurred residues

Incurred residues

Green tea



Spikes (0.1 mg/kg)

Analyte	Recovery QuEChERS	Recovery Pure MeCN
Bifenthrin	78%	88%
Buprofezin	87%	78%
Chlorfenapyr	121%	115%
Cyhalothrin-λ	101%	86%
Endosulfan-α	77%	80%
Endosulfan-β	77%	80%
Endosulfan-SO ₄	90%	82%
Fenvalerate	94%	89%
Triazophos	106%	84%


(5) Matrix swelling prior the extraction

■ CONCLUSIONS

- **Matrix swelling (hydration)** represents an **important step** during the analysis of pesticide residues in tea
- **Lower extraction efficiency of incurred residues** from tea samples if **pure MeCN** without matrix hydration is used although acceptable recoveries can be obtained by the analysis of spikes

(6) Long-term stability of analytes signal

- Repeated injections of 0.1 mg/kg spiked tea extracts




(7) Method validation

- **Recovery (%) & Repeatability (RSD, %)**
 - **Matrix: Black tea (levels: 0.01 , 0.1 , 1 mg/kg) $n=6$**
Green tea (levels: 0.01 , 0.1 , 1 mg/kg) $n=6$
- **LOQ**
 - **Expressed as the lowest calibration level (LCL)**


(7) Method validation

■ Recovery (%)




(7) Method validation


■ Repeatability (RSD, %)



■ LCLs For most of analytes ≤ 0.005 mg/kg (Corresponding to ≤ 0.25 ng/mL)

(7) Method validation






3. GC/MS/MS Analysis



7000 GC/MS/MS – Compact Design




7890 / 7000B GC/MS/MS System



Target analytes for Tea analysis by GC/MS/MS

1 Aldrin	26 Chlorpyrifos-methyl	51 Endosulfan-sulfate	76 Isufenphos	101 Phenothrin (sum)	126 Tetradifon
2 Ametryn	27 Chlozolinate	52 Endrin	77 Isufenphos-methyl	102 Phenthoate	127 Thiometon
3 Azinphos-ethyl	28 Cyfluthrin (sum)	53 Ethion	78 Kresoxim-methyl	103 Phenylphenol, <i>o</i> -	128 Tolclofos-methyl
4 Azinphos-methyl	29 Cyhalothrin-lambda	54 Ethoprophos	79 Lindane	104 Phosalone	129 Tolyfluanid
5 Azoxystrobin	30 Cypermethrin (sum)	55 Etrimfos	80 Malaoxon	105 Phosmet	130 Triadimefon
6 BHC-alpha	31 Cyprodinil	56 Fenamidone	81 Malathion	106 Pirimiphos-ethyl	131 Triadimenol
7 BHC-beta	32 DDD, <i>o</i> , <i>p</i> '-	57 Fenamiphos	82 Mecarbam	107 Pirimiphos-methyl	132 Triazophos
8 BHC-delta	33 DDD, <i>p</i> , <i>p</i> '- + DDT, <i>o</i> , <i>p</i> '-	58 Fenarimol	83 Metazachlor	108 Procymidone	133 Trifloxystrobin
9 Bifenthrin	34 DDE, <i>o</i> , <i>p</i> '-	59 Fenchlorphos	84 Methacrifos	109 Profenofos	134 Trifluralin
10 Bromophos-methyl	35 DDE, <i>p</i> , <i>p</i> '-	60 Fenitrothion	85 Methidathion	110 Propargite	135 Vinclozolin
11 Bromophos-ethyl	36 DDT, <i>p</i> , <i>p</i> '-	61 Fenoxycarb	86 Methiocarb	111 Propham	
12 Bromopropylate	37 Deltamethrin	62 Fenthion	87 Methoxychlor	112 Prothiofos	
13 Bupirimate	38 Diazinon	63 Fenvalerate (sum)	88 Myclobutanil	113 Pyrazophos	
14 Buprofezin	39 Dichlobenil	64 Flucythrinate (sum)	89 Naled	114 Pyridaben	
15 Cadusafos	40 Dichlofluanid	65 Fludioxonil	90 Nitrofen	115 Pyridaphenthion	
16 Carbophenothion	41 Dichlorvos	66 Fluvalinate (sum)	91 Nuarimol	116 Quinalphos	
17 Chinomethionat	42 Diclofop-methyl	67 Fonofos	92 Oxychlordane	117 Quintozene	
18 Chlordane, <i>cis</i> -	43 Dicloran	68 Haloxyfop-ethoxyethyl	93 Oxyfluorfen	118 Resmethrine (sum)	
19 Chlordane, <i>trans</i> -	44 Dieldrin	69 Haloxyfop-methyl	94 Paraoxon	119 Sulfotep	
20 Chlorfenapyr	45 Difenoconazol (sum)	70 Heptachlor	95 Parathion-ethyl	120 Tebuconazole	
21 Chlorfenvinphos	46 Diphenylamine	71 Heptachlor-epoxide (endo)	96 Parathion-methyl	121 Tecnazene	
22 Chlorobenzilate	47 Disulfoton	72 Heptachlor-epoxide (exo)	97 Penconazole	122 Tefluthrin, <i>cis</i> -	
23 Chlorothalonil	48 Disulfoton-sulfone	73 Heptenophos	98 Pencycuron	123 Terbufos	
24 Chlorpropham	49 Endosulfan-alpha	74 Hexachlorobenzene	99 Pendimethalin	124 Terbufos-sulfone	
25 Chlorpyrifos	50 Endosulfan-beta	75 Iprodione	100 Permethrin (sum)	125 Tetraconazole	

RTL GC/MS/MS Method – Set up




Load method, re-lock method on customer site and acquisition method is set-up


RTL 2x Pesticide Method + Post-Column Backflush

(1) Constant Pressure, Column 15 m X 0.25 mm id X 0.25 μ m HP-5MSUI

(2) Constant Pressure, Restrictor 0.50 m X 0.15 mm id x 0.15 μ m DB-5MSUI



Matrix Compounds may be left in the Column after each Injection



- These heavy matrix materials build up and travel further into the column with each injection.
- This build-up of heavy materials can cause retention time shifts, peak distortion, higher bleed, and loss of sensitivity

Back flushing After Each Injection



Back flush 20 sec

← Flow



Back flush 60 sec

← Flow



Back flush 90 sec

← Flow




Back flush 120 sec

← Flow




- Back flushing removes heavy materials after each injection, heaviest first

Post column, post run back flush




During the run, $P_1 > P_2$, Carrier gas flows towards MS

Post column, post run back flush




During the run, $P_1 > P_2$, Carrier gas flows towards MS

Post column, post run back flush



During post run back flush $P_1 < P_2$, **Carrier gas flow reversed in column**

Post-Column Back flush set-up – Step 1



Post-Column Back flush set-up – Step 2

Summary of Backflush Calculations

End Run and Use Post Run to Perform Backflush at:
20.003 min

Operating Conditions During Post Run Backflush

Oven Temperature: 280 °C
Transfer Line Temperature: 280 °C

Detector	Maximum Flow	Allowable Pressure	Flow at Chosen Pressure
Vacuum	75	94.318	35.212

Backflush Pressure: 60 psi
Inlet Pressure during Backflush: 1 psi
Void Volumes: 12.471


Backflush Flow to Inlet: 8.6566 mL/min
Void Time: 0.16038 min
Backflush Time: 2 min

OK Cancel Help


Annotations:

- Post-run oven temperature (points to Oven Temperature)
- Flow to MS during backflush (points to Flow at Chosen Pressure)
- Post-run pressure applied to Purged union (points to Backflush Pressure)
- Post-run pressure applied to column inlet (points to Inlet Pressure during Backflush)
- Reverse flow rate to inlet (points to Backflush Flow to Inlet)
- Post-run time (points to Backflush Time)

Post-Column Back flush set-up – Step 3




Back flush saves cycle time



Long-term stability of analytes signal



- Repeated injections of 0.1 mg/kg spiked tea extracts



Benefits of Backflush

- Consistent analyte retention times and responses
- Robust chromatography and consistent analyte chromatographic peak shapes
- Prevention of high boiling matrix from contaminating the MS ion source
- Extended column life-time and reduced cycle times by removing the need for high-temperature bake-out between runs







4. Results




TIC MRM :Tea Matrix Spiked at 60 ng/mL



Matrix matched Tea Standard - Dichlobenil





Ion ratio
+/- 20%





Early eluting Benzo-nitrile herbicide

Matrix matched Tea Standard – p,p'-DDT




Mid- eluting Organo-chlorine insecticide

Matrix matched Tea Standard - Deltamethrin



Late eluting Pyrethroid insecticide

Pesticide residues in Tea from Czech market: Results of pilot study ($n=45$)



Pesticide residues in tea samples ($n=45$)

Data from GC–MS/MS and LC–MS/MS analyses

Pesticide	Number of positive samples (45 Tea samples)	Maximum concentration (mg/kg)	EU MRL* (mg/kg)	Number of MRL exceedances
Chlorpyrifos-methyl	22	1.763	0.1	8
Imidacloprid	12	0.294	0.05	7
Acetamiprid	7	0.623	0.1	2
Myclobutanil	1	0.459	0.05	1
Pyraclostrobin	1	0.052	0.05	1
Thiophanate-Me	1	1.329	0.1	1
Acephate	1	0.091	0.05	1
Metolcarb	1	0.05	0.01	1
Monocrotophos	1	0.855	0.1	1
Triazophos	1	0.461	0.02	1
Paclobutrazol	1	0.056	0.02	1
Carbendazim	2	0.284	0.1	1
Buprofezin	2	0.09	0.05	1
Methomyl	2	0.303	0.1	1
Thiametoxam	2	0.169	0.1	1

* Regulation (EC) No 396/2005, Tea (dried leaves and stalks, fermented or otherwise, of *Camelia sinensis*)

Pesticide residues in tea samples ($n=45$)

Confirmed, quantified but compliant with EU MRL

Pesticide	Number of positive samples	Maximum concentration (mg/kg)	EU MRL* (mg/kg)
Endosulfan	16	0.107	30
Bifenthrin	13	0.862	5
Cyhalothrin-lambda	10	0.123	1
Propargite	6	0.119	5
Cypermethrin	5	0.312	0.5
Chlorpyrifos	3	0.081	0.1
DDT	3	0.03	0.2
Ethion	2	0.838	3
Metalaxyl	2	0.021	0.1
Pyridaben	2	0.048	0.05

* Regulation (EC) No 396/2005, Tea (dried leaves and stalks, fermented or otherwise, of *Camelia sinensis*)



Finding a Tea Matrix Blank - Analysis of Organic Tea Samples




Pesticide Residues in 4 Organic Tea Extracts by GC/MS/MS

Analyte	Ret Time (min)	Organic A	Organic G	Organic D	Organic W
		ppb	ppb	ppb	ppb
Biphenyl	3.55	0.33	0.67	0.24	0.41
o-Phenyl phenol	4.42	0.25	0.45	0.24	0.24
a-Endosulfan	11.32	0.89			
p,p'-DDE	12.02	0.18			
b-Endosulfan	12.59	1.15			
Endosulfan sulfate	13.40	3.31		0.60	
pp'-DDT	13.51	0.45		0.27	

Σ Endosulfans = 5.35 ppb, Equivalent to 0.027 mg/kg
 EU MRL (Tea) = 30 mg/kg



Endosulfan Sulfate Transitions in 4 Organic Tea Extracts




Endosulfan Sulfate (?) in Organic Tea Extract

– Using 2 MRM Transitions

Reagent blank

Matrix Matched
Std 5 ppb

Organic Tea A



Optimized Endosulfan Sulfate MRM Transitions in the G9250AA Database

Security Warning: Macros have been disabled. Options...


C6572 fx =IF(D6572=D6571,C6571,C6571+1)

	R	S	T	U	V	W	X	Y	Z	AA	AB	AC	AD	AE	AF
6558	35694-06-5	2,2',3,4,4',5-Hexachlorobiphenyl (E	13.44	false	289.9	LowRes	219.9	LowRes	10	35	0.1	1020	48%	Q5	
6559	35694-06-5	2,2',3,4,4',5-Hexachlorobiphenyl (E	13.44	false	359.9	LowRes	287.8	LowRes	10	25	0.1	970	46%	Q6	
6560	35694-06-5	2,2',3,4,4',5-Hexachlorobiphenyl (E	13.44	false	324.8	LowRes	289.8	LowRes	10	15	0.1	960	45%	Q7	
6561	35694-06-5	2,2',3,4,4',5-Hexachlorobiphenyl (E	13.44	false	361.9	LowRes	291.8	LowRes	10	25	0.1	870	41%	Q8	
6562	35694-06-5	2,2',3,4,4',5-Hexachlorobiphenyl (E	13.44	false	218.0	LowRes	183.0	LowRes	10	20	0.1	800	38%	Q9	
6563	35694-06-5	2,2',3,4,4',5-Hexachlorobiphenyl (E	13.44	false	361.9	LowRes	326.8	LowRes	10	15	0.1	750	35%	Q10	
6564	3735-33-9	Phosmet oxon	13.49	false	160.0	LowRes	77.0	LowRes	10	25	0.1	480	100%	Q0	
6565	3735-33-9	Phosmet oxon	13.49	false	160.0	LowRes	51.0	LowRes	10	40	0.1	230	48%	Q1	
6566	3735-33-9	Phosmet oxon	13.49	false	160.0	LowRes	133.0	LowRes	10	15	0.1	220	46%	Q2	
6567	3735-33-9	Phosmet oxon	13.49	false	133.0	LowRes	76.9	LowRes	10	15	0.1	40	8%	Q3	
6568	3735-33-9	Phosmet oxon	13.49	false	133.0	LowRes	50.9	LowRes	10	30	0.1	30	6%	Q4	
6569	3735-33-9	Phosmet oxon	13.49	false	133.0	LowRes	104.8	LowRes	10	5	0.1	20	4%	Q5	
6570	3735-33-9	Phosmet oxon	13.49	false	301.0	LowRes	191.8	LowRes	10	10	0.1	0	0%	Q6	
6571	3735-33-9	Phosmet oxon	13.49	false	172.8	LowRes	104.0	LowRes	10	15	0.1	0	0%	Q7	
6572	1031-07-8	Endosulfan sulfate	13.38	false	271.9	LowRes	237.0	LowRes	10	15	0.1	1480	100%	Q0	
6573	1031-07-8	Endosulfan sulfate	13.38	false	273.8	LowRes	238.9	LowRes	10	15	0.1	1080	73%	Q1	
6574	1031-07-8	Endosulfan sulfate	13.38	false	273.8	LowRes	236.9	LowRes	10	15	0.1	520	35%	Q2	
6575	1031-07-8	Endosulfan sulfate	13.38	false	271.9	LowRes	235.0	LowRes	10	15	0.1	300	20%	Q3	
6576	1031-07-8	Endosulfan sulfate	13.38	false	229.0	LowRes	194.0	LowRes	10	25	0.1	180	12%	Q4	
6577	1031-07-8	Endosulfan sulfate	13.38	false	386.6	LowRes	288.7	LowRes	10	5	0.1	80	5%	Q5	
6578	1031-07-8	Endosulfan sulfate	13.38	false	386.6	LowRes	252.7	LowRes	10	10	0.1	70	5%	Q6	
6579	52663-66-8	2,2',3,3',4,5'-Hexachlorobiphenyl (I	13.48	false	359.9	LowRes	289.9	LowRes	10	30	0.1	1300	100%	Q0	
6580	52663-66-8	2,2',3,3',4,5'-Hexachlorobiphenyl (I	13.48	false	287.9	LowRes	217.9	LowRes	10	35	0.1	1240	95%	Q1	
6581	52663-66-8	2,2',3,3',4,5'-Hexachlorobiphenyl (I	13.48	false	359.9	LowRes	324.9	LowRes	10	15	0.1	1190	92%	Q2	
6582	52663-66-8	2,2',3,3',4,5'-Hexachlorobiphenyl (I	13.48	false	324.8	LowRes	289.8	LowRes	10	20	0.1	1170	90%	Q3	
6583	52663-66-8	2,2',3,3',4,5'-Hexachlorobiphenyl (I	13.48	false	361.9	LowRes	289.9	LowRes	10	30	0.1	820	63%	Q4	
6584	52663-66-8	2,2',3,3',4,5'-Hexachlorobiphenyl (I	13.48	false	289.9	LowRes	219.9	LowRes	10	35	0.1	790	61%	Q5	
6585	52663-66-8	2,2',3,3',4,5'-Hexachlorobiphenyl (I	13.48	false	289.9	LowRes	217.9	LowRes	10	35	0.1	790	61%	Q6	
6586	52663-66-8	2,2',3,3',4,5'-Hexachlorobiphenyl (I	13.48	false	361.9	LowRes	326.8	LowRes	10	15	0.1	710	55%	Q7	
6587	52663-66-8	2,2',3,3',4,5'-Hexachlorobiphenyl (I	13.48	false	318.0	LowRes	183.0	LowRes	10	20	0.1	680	52%	Q8	

MRM Table DATABASE DB Compound List My Target Compound List CF, 40-min Method CP, 40-min Method


Ready Average: 52049.66963 Count: 231 Sum: 7286953.748 80%

Endosulfan Sulfate (✓) in Organic Tea Extract – Using 4 MRM Transitions




... demonstrates the value of access to multiple optimized transitions in the G9250AA MRM Database
-Not just to avoid matrix interferences but also for additional confirmation

Endosulfan α , β and Sulfate - 5975 EI Mode




Endosulfan α , β and Sulfate - 5975 NCI Methane Mode



Endosulfan α , β and Sulfate (✓) in Organic Tea Extract A

5975 NCI Methane SIM



Summary

- Simplified, rapid and optimized sample preparation for tea based on a modified QuEChERS regime with Liq/Liq extraction
- Sensitive detection – LCLs for most analytes ≤ 0.005 mg/kg
- Excellent recoveries (70–120%) down to 0.01 mg/kg level
- Excellent repeatability $\leq 20\%$
- RTL for ease of method set-up / maintenance
- Capillary Flow Technology: Back flush – GC method robustness
- MS/MS Selectivity – low matrix interferences
- Method suitable for detection and confirmation of pesticide residues well below EU MRLs (Regulation (EC) No.396/2005)



Acknowledgements

**Prof. Jana Hajslova, Dr. Tomas Cajka, Dr. Jana Pulkrabova,
Veronika Bachanova, Lucie Drabova, Kamila Kalachova**

**Department of Food Analysis and Nutrition
Faculty of Food and Biochemical Technology
Institute of Chemical Technology
Prague
Czech Republic**

■ jana.hajslova@vscht.cz ■

QUESTIONS?




Please Submit your questions by:

A screenshot of a web form for submitting questions. It features a large text input box with the instruction "Typing your question in the box and hitting submit." and a "Submit a Question" button at the bottom left. A vertical scrollbar is visible on the right side of the text box.

E:mail: chris_sandy@agilent.com

THANK YOU!



We apologize for this brief interruption but we are experiencing technical difficulties and will resume shortly.

