

GC Q-TOF Analysis of Fluorinated Alkyl Compounds in Bio-solid Matrix

Anthony Macherone, Ph.D.



Sr. Applications Chemist, Agilent Technologies
Visiting Scientist, Johns Hopkins School of Medicine

Outline

Brief Review: The human Genome and exposome

Measuring the exposome

Analysis of Fluorotelomer alcohols (FTOH)

Summary

Q&A

GC Q-TOF analysis of fluorinated alkyl compounds in bio-solid matrix

THE HUMAN GENOME AND EXPOSOME



Sequencing and mapping

Elucidation of gene expression and protein function

Identification biochemical pathways implicated in chronic diseases

- Cancer, diabetes, vascular and neurodegenerative

Opportunities for improved treatment and patient management

- Tailored drug treatments and gene therapy

Wild, CP. Cancer Epidemiol Biomarkers Prev 2005;14:1847-1850



Application of genetic testing

Identify individuals at greatest risk for chronic diseases

New possibilities for disease prevention.

- Moves away from individual observation (clinical) towards population observation (epidemiology)
- Reduced morbidity and mortality
 - Utilizing public health measures

Our inheritance, our future: realising the potential of genetics in the NHS. Retrieved April 15, 2013 from, http://webarchive.nationalarchives.gov.uk/+/www.dh.gov.uk/en/Publicationsandstatistics/Publications/PublicationsPolicyAndGuidance/DH_4006538



Exposure-Disease associations

Environmental exposures play an important role in many chronic diseases

Most exposure-disease associations remain ill defined

Genetic variants or in the human genome are of high prevalence but low penetrance

- Environmental exposures may contribute to population disease burden
- Questions if genetic knowledge has or will improve understanding of disease etiology at the population level

Vineis P, Schulte P, McMichael AJ. Lancet 2001;357:709 – 12

Peto J. Nature 2001;411:390 – 5

Cappuccio FP. Int J Epidemiol 2004;33:387 – 8

The Exposome and Exposomics

Exposome:


- CP Wild (2005) defined the exposome as the sum of all exogenous and endogenous exposures over a complete lifetime.
 - *Exogenous*: food, air, water and lifestyle choices
 - *Endogenous*: biosynthesized chemicals and metabolites

Exposomics:

- Measuring the exposome using “omics” tools: LC and GC mass spectrometry, NMR, bioinformatics



New paradigm: exposome and exposomics



S.M. Rappaport, J Expo Sci Environ Epidemiol, (2011) 21, pg. 6

Exposomics Funding and Studies

US:

- CDC, NIOSH, NIEHS <http://www.cdc.gov/niosh/topics/exposome/>
- EPA <http://epa.gov/ncct/edr/non-monotonic.html>

UK Biobank:

- Long-term follow-up of health (500,000 people ages 40 – 69)
<http://www.ukbiobank.ac.uk/2012/06/following-health/>

EC:

- €17.3 million, Four year study. Nature 491:647 (2012)

China:

- National Basic Research Program of China J. Trans. Med. 10:A41 (2012)

Agilent and Exposomics in 2013

ASMS

- Monitoring Steroidal analogues in clinical and environmental chemistry: one model for exposomics.
- The role of GC/QTOF in exposomics

Agilent Environmental Seminars

- Harmonized GC Triple Quadrupole analysis of steroidal analogues for clinical and environmental monitoring of the exposome

Genetic Engineering & Biotechnology News

The future of GC/Q-TOF in environmental analysis

- Advanced Techniques in Gas Chromatography-Mass Spectrometry for Environmental Chemistry

Agilent Academic Newsletter

- Top-Down Exposomics: Profiling Urinary Organic Acids

U Penn Multi-omics symposium

- Exposomics and Exposure: Biological and Environmental Monitoring of 17 β -Estradiol

GC Q-TOF analysis of fluorinated alkyl compounds in bio-solid matrix

MEASURING THE EXPOSOME



Measuring the exposome

Purpose: identify epidemiological risk factors

Multi-omics tools

EWAS: Exposome-wide Association Studies

Hybrid analytical methods: simultaneous targeted and non-targeted analysis


Mass spectrometry, sensor technology, imaging technology

Wild, CP. Int J Epidemiol 2012;41:24-32



Instrumentation

Estimated limits of detection (x-axis) and the typical number of compounds identified (y-axis) by NMR, GC/MS and LC/MS for typical metabolomics studies



Wishart D.S. Exploring the Human Metabolome by Nuclear Magnetic Resonance Spectroscopy and Mass Spectrometry. Lutz N, Sweedler J, Wevers R, (eds.) In: Methodologies for Metabolomics: Experimental Strategies and Techniques, New York, NY, 2013. Reprinted with the permission of Cambridge University Press and David S. Wishart



Agilent GC/MS portfolio in exposomics

Metabolomics, food safety, water analysis, toxicology, clin chem, epidemiology, disease correlation

7890B/

- 7200A, 5977A, 5975T, 240IT, 7000B



What to measure? What would harm us?

Sediments, debris, chemical spills, fuel leaks, fires, open air burning...

Salted farm lands

Huge amount of debris to be treated

- Particulate matters from temporary storage sites

Radioactive debris

Chemical pollution from all those sources

- Impact on ecosystem
- Impact on human health

Targeted and Untargeted Analysis

Targeted Analysis

- List of chemicals to be monitored
- Standards for quantification

Untargeted Analysis

- When targets are unknown
- No standard available

Need to combine targeted and untargeted analysis especially in times of trouble

- Understand fate and transport
- Predict / determine toxicity
- Correlate exposure to disease



Bottom-up vs. top-down exposomics

Bottom-up:

- Samples air, water, food, etc.
- Ties exposures directly to their immediate sources
- Misses essential features of the internal chemical environment


Top-down:

- Chemicals measured in blood, urine, tissue, etc.
- Identifies all potentially important exposures
- Provides no information about exogenous sources

Rappaport, SM. J Expos Sci Env Epidemiol 2011;21:5–9



Really need to look at the whole: Cyclical Exposomics



Biomarkers

Exposure vs. Disease

Not always one and the same

Organic acidurias

- Genetic: Maple syrup Urine disease
- Exogenous: exposure xylene isomers results in increased levels of methyl benzoate and methylhippurate
- Iatrogenic sources: administration of valproic acid results elevated levels of in 3-hydroxyisovalerate, 5-hydroxyhexanoate, 7-hydroxyoctanoate, p-hydroxyphenylpyruvate and dicarboxylic
- Gastrointestinal dysbiosis: increased urinary levels of citramalic acid, 5-hydroxy-methyl-furoic acid and tartaric acid arising respectfully from *Saccharomyces*, *Aspergillus* and *Candida* species overgrowth in the bowel
 - Tartaric aciduria from yeast overgrowth in the bowel, there has been at least one report of the manifestation of Autism like symptoms in twin brothers

Lord RS, Alexander **Bralley** J Alt Med Rev 2008;13:205-215

Kumps A, Duez P, Mardens Y. Clin Chem 2002;48:708–717

Shaw W, Kassen E, Chaves E. Clin Chem 1995;41:1094-1104



Meet in the middle exposomics

Exposure and disease

Identifying the overlap between exposure markers and predictive biomarkers of disease

Biomarkers of exposure elevated in certain disease states:

- Causal link of exposure and disease

Potentially new avenues of prevention by identifying environmental factors of disease

Chadeau-Hyam M, Athersuch TJ, Keun HC, et al. Biomarkers 2011;16(1):83-88

GC Q-TOF analysis of fluorinated alkyl compounds in bio-solid matrix

ANALYSIS OF FTOH



Fluorotelomer Alcohols (FTOH)

Source unknown, but probably intermediate degradation products of fluorinated polymers

Oxidize to form fluorinated carboxylic acid, some of which are toxic

Methods needed for studying their transport and fate in the environment

Have the form $\text{F}_3\text{C}(\text{CF}_2)_{\text{N}-1}(\text{CH}_2)_{\text{M}}\text{OH}$

N:M FTOH is the shorthand notation



Purpose of Study

Develop robust, sensitive and selective GC-MS/MS method in bio-solid matrix from waste water treatment plants

Study the fate and transport of fluorotelomers in the environment

Investigate degradation pathways of fluoropolymers

Gain understanding of human exposure and toxicity

- Evaluate gas chromatography-quadrupole time of flight mass spectrometry (GC Q-TOF MS)

FTOH Disease associations

Endocrine disruption

Estrogen-Like Properties of Fluorotelomer Alcohols as Revealed by MCF-7 Breast Cancer Cell Proliferation¹

- FTOH act as xenoestrogens *in vitro*
- Structural similarities of these compounds and 4-nonylphenol (reference xenoestrogen)
- Ligands for the estrogen receptor

Fluorotelomer alcohols induce hepatic vitellogenin through activation of the estrogen receptor in male medaka (Oryzias latipes)²

- Vitellogenin: biomarker for estrogen EDC exposure
- Dose-dependent interaction between medaka estrogen receptor alpha (ERalpha) and coactivator TIF2
- Estrogenic effects for medaka ERalpha descended in the order of estradiol >> 6:2 FTOH > 8:2 FTOH

¹Maras M, Vanparys C, Muylle Environ F, et al. Health Perspect 2006;114(1):100–105.

²Ishibashi H, Yamauchi R, Matsuoka M, et al. Chemosphere 2008;71(10):1853-9.

GC/Q-TOF Method

15 m HP-INNOWax column (19091N-131)

Oven Program

- 60 °C for 1 min
- 3 °C/min to 75 °C for 0 min
- 20 °C/min to 210 °C for 0 min

MMI: 2 µl cold, splitless injections

- 65 °C (0.01 min), 300 °C/min to 250 °C

Two minute post column Backflush

Positive CI, 20% CH₄


Capillary Flow Technology: Backflush

Advantages


- Reduces analysis time (more samples per day)
- Reduces contamination (and cleaning frequency) for the source
- Increases lifetime of analytical column
- Consistent GC retention times
- Consistent MS spectra through a sample sequence
 - Reduces chemical noise
 - Higher quality quantitation




Post-Column Example – Analysis




Post-Column Example - Backflush



No BF



With BF



Standards and Samples

Solvent standards, 100 ng/μl in MTBE

Two bio-solid extracts prepared in EtOAc

- One 1% and the other 5% lime treated
- Each extract split 50/50
 - One spiked with analytes and IS @ 50 ng/ml, the other not spiked


ID	Solvent	Vol (mL)	Description
1	Ethyl acetate	1	Extracted solvent, 1% lime treated biosolid
2	Ethyl acetate	1	Standards in extracted solvent
3	Ethyl acetate	1	Extracted solvent, 5% lime treated biosolid
4	Ethyl acetate	1	Standards in extracted solvent



Analytes and ISTDs


Acronym ▼	Type ▼	r.t. ▼	Exact Mass + H ▼
4:2 FTOH	Target	4.078	265.0269
5:1 FTOH	Target	5.221	301.0081
6:2 FTOH	Target	5.773	365.0206
6:1 FTOH	Target	6.188	351.0049
7:2 sFTOH	Target	6.237	415.0174
7:1 FTOH	Target	6.981	401.0017
8:2 MFTOH	ISTD	7.449	469.0334
8:2 FTOH	Target	7.471	465.0142
8:1 FTOH	Target	7.598	450.9985
9:1 FTOH	Target	8.098	500.9953
10:2 FTOH	Target	8.470	565.0078
10:1 FTOH	Target	8.523	550.9921
11:1 FTOH	Target	8.886	600.9889
d7-MeFOSE	ISTD	11.868	565.0466
MeFOSE	Target	11.889	558.0026
d9-EtFOSE	ISTD	11.897	581.0748
EtFOSE	Target	11.928	572.0183

Chromatogram




4:2 FTOH


+ TIC Scan Extracted Solvent 3.d




+ EIC(227.0087) Scan Extracted Solvent 3.d




227.0087 , 265.0253
Ratio = 77.0 (97.6 %)



4:2 FTOH - 1 Levels, 1 Levels Used, 1 Points, 1 Points Used, 2 QCs
 $y = 1896.210189 * x + 1.192093E-011$
 $R^2 = 1.00000000$




11:1 FTOH




7:2 sFTOH


+ TIC Scan Extracted Solvent 1.d



+ EIC(376.9986) Scan Extracted Solvent 1.d




7:2 sFTOH - 1 Levels, 1 Levels Used, 1 Points, 1 Points Used, 2 QCs




MeFOSE


+ TIC Scan Extracted Solvent 1.d



+ EIC (558.0026) Scan Extracted Solvent 1.d




MeFOSE - 1 Levels, 1 Levels Used, 1 Points, 1 Points Used, 2 QCs
 $y = 239.086133 * x + 0.000000E+000$
 $R^2 = 1.00000000$



Uncorrected Mass Accuracy (IRM off)

Acronym ▼	Formula ▼	Exact Mass + H ▼	Observed Mass ▼	Δppm ▼
4:2 FTOH	C ₆ H ₅ F ₉ O	265.0269	265.0270	-0.3773
6:2 FTOH	C ₈ H ₅ F ₁₃ O	365.0206	365.0206	0.0000
8:2 FTOH	C ₁₀ H ₅ F ₁₇ O	465.0142	465.0140	0.4301
10:2 FTOH	C ₁₂ H ₅ F ₂₁ O	565.0078	565.0078	0.0000
7:2 sFTOH	C ₉ H ₅ F ₁₅ O	415.0174	415.0190	-3.8553
5:1 FTOH	C ₆ H ₃ F ₁₁ O	301.0081	301.0079	0.6644
6:1 FTOH	C ₇ H ₃ F ₁₃ O	351.0049	351.0050	-0.2849
7:1 FTOH	C ₈ H ₃ F ₁₅ O	401.0017	401.0016	0.2494
8:1 FTOH	C ₉ H ₃ F ₁₇ O	450.9985	450.9985	0.0000
9:1 FTOH	C ₁₀ H ₃ F ₁₉ O	500.9953	500.9956	-0.5988
10:1 FTOH	C ₁₁ H ₃ F ₂₁ O	550.9921	550.9922	-0.1815
11:1 FTOH	C ₁₂ H ₃ F ₂₃ O	600.9889	600.9896	-1.1647
MeFOSE	C ₁₁ H ₈ F ₁₇ N O ₃ S	558.0026	558.0042	-2.8674
EtFOSE	C ₁₂ H ₁₀ F ₁₇ NO ₃ S	572.0183	572.0167	2.7971

Resolution (50 pg on column)



Loss Elucidation




Mass accuracy in facilitated structure elucidation of fragment ions

- Nominal loss of 38 m/z for n:2 FTOH (n=4,6,7,8,10)
 - Loss of F₂ ruled out by mass
 - Loss of H₂O and HF confirmed
 - MS/MS confirmed –H₂O, -HF losses
- Loss of HF determined for n:1 FTOH (n=5,6,7,8,9,10,11)
- Loss of H₂O determined for perfluorosulfonamidoethanols



Structure Elucidation

Nominal loss = 38 m/z





Acronym	Observed Base Peak m/z	Molecular ion $-F_2$ m/z	Δ ppm	Molecular ion $-H_2O$, $-HF$ m/z	Δ ppm
4:2 FTOH	227.0104	227.0301	86.7726	227.0102	-0.8810

Loss Elucidation & Mass Accuracy

Acronym	Formula	Loss	Exact Mass-Loss	Observed Mass 2	Δ ppm 2
4:2 FTOH	$C_6H_5F_9O$	-H ₂ O, -HF	227.0102	227.0105	-1.3215
6:2 FTOH	$C_8H_5F_{13}O$	-H ₂ O, -HF	327.0038	327.0040	-0.6116
8:2 FTOH	$C_{10}H_5F_{17}O$	-H ₂ O, -HF	426.9972	426.9975	-0.7026
10:2 FTOH	$C_{12}H_5F_{21}O$	-H ₂ O, -HF	526.9910	526.9915	-0.9488
7:2 sFTOH	$C_9H_5F_{15}O$	-H ₂ O, -HF	377.0006	377.0011	-1.3263
5:1 FTOH	$C_6H_3F_{11}O$	-HF	281.0017	281.0019	-0.7117
6:1 FTOH	$C_7H_3F_{13}O$	-HF	330.9988	330.9990	-0.6042
7:1 FTOH	$C_8H_3F_{15}O$	-HF	380.9954	380.9956	-0.5249
8:1 FTOH	$C_9H_3F_{17}O$	-HF	430.9923	430.9927	-0.9281
9:1 FTOH	$C_{10}H_3F_{19}O$	-HF	480.9894	480.9896	-0.4158
10:1 FTOH	$C_{11}H_3F_{21}O$	-HF	530.9860	530.9864	-0.7533
11:1 FTOH	$C_{12}H_3F_{23}O$	-HF	580.9834	580.9835	-0.1721
MeFOSE	$C_{11}H_8F_{17}NO_3S$	-H ₂ O	539.9920	539.9943	-4.2593
EtFOSE	$C_{12}H_{10}F_{17}NO_3S$	-H ₂ O	554.0077	554.0100	-4.1516

MS/MS



Summary of Study

Backflush required for retention time reproducibility

- Ranged from 2.3% - 6.9% RSD in the matrix

Detection limits (LOD) as determined by a S/N ratio of $\geq 10:1$

- Ranged from 2 fg – 100 fg on column

Excellent Mass Accuracy

- Typically < 2 ppm
- Range 0.1 ppm – 5 ppm

Mass Resolution $> 12,000$ for most samples (up to 20,000)

- Allowed differentiation from heavy chemical background

Summary of exposomics

Exposomics paradigm focuses on the environmental impact of exposure on disease

- New opportunities
 - Government funding and interest (academia and industry)
- Leverages GC/Q-TOF (GC/MS in general) power in environmental, food safety *and* biological monitoring
- Cyclical exposomics
 - Hybrid targeted / non-targeted panels focused on biologically relevant chemotypes
 - Multi disciplinary, multi-technique analyses

Software and Bioinformatics will be an integral component

- Multi-variant analysis
- Identify chemical fingerprints and biomarkers

GC/TOF and GC/Q-TOF in exposomics

Exposomics is a nascent field of study but GC/Q-TOF will be a primary tool

- Strong history in environmental analysis, current methods are easily ported
- Global monitoring of environmental contaminants
- Discovery new biomarkers for disease



Acknowledgements

The authors would like to gratefully acknowledge the United States Environmental Protection Agency and especially Shoji Nakayama, M.D., Ph.D. and Marc Mills, Ph.D. for providing standards and prepared samples for the study presented herein.



Disclaimer

The United States Environmental Protection Agency through its Office of Research and Development contributed to the research described here. Mention of trade names or commercial products does not constitute endorsement or recommendation for use by USEPA.



Bibliography

Cappuccio FP. Int J Epidemiol 2004;33:387-8

Chadeau-Hyam M, Athersuch TJ, Keun HC, et al. Biomarkers 2011;16(1):83-88

Ishibashi H, Yamauchi R, Matsuoka M, et al. Chemosphere 2008;71(10):1853-9

Kumps A, Duez P, Mardens Y. Clin Chem 2002;48:708–717

Lord RS, Alexander Bralley J Alt Med Rev 2008;13:205-215

Macherone, A, Nakayama SF. Dioxin 2011, Brussels, Belgium

Maras M, Vanparys C, MuylleEnviron F, et al. Health Perspect 2006;114(1):100–105

Nakayama SF, Macherone A, Kidus Tadele K and Mills MA. Dioxin 2010, San Antonio, TX

Our inheritance, our future: realising the potential of genetics in the NHS. Retrieved April 15, 2013 from, http://webarchive.nationalarchives.gov.uk/+www.dh.gov.uk/en/Publicationsandstatistics/Publications/PublicationsPolicyAndGuidance/DH_4006538

Peto J. Nature 2001;411:390 – 5

Rappaport, SM. J Expos Sci Env Epidemiol 2011;21:5–9

Shaw W, Kassen E, Chaves E. Clin Chem 1995;41:1094-1104

Vineis P, Schulte P, McMichael AJ. Lancet 2001;357:709 – 12

Wild, CP. Cancer Epidemiol Biomarkers Prev 2005;14:1847-1850

Wild, CP. Int J Epidemiol 2012;41:24-32

Thank you!

Questions?

