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
```
Method Report                    Date 08/19/2013  Time 12:35:06

                          *** Method Report ***

Method File (modified)

     Last Update              Date 8/19/2013  Time 12:34:26 PM
     Operator                 J
     Product                  Uv-visible ChemStation, Rev. B.04.03
     Number of Standards      5
     Auxiliary Spectra        0

Method Information:
     Default Method

Method Options:

     Pre measure macro:       none
     Post measure macro:      none
     Store Spectra to File:   none
     Auto Analyze:            off

     Single Analysis

Sampling System:         Sipper

     Pump Direction:          Clockwise
     Pump Time (s):           20
     Wash Time (s):           0
     Wait Time (s):           3
     Sample Return (%):       0
     Air Segment (s):         0

Data Analysis:

     Spectral Processing:

          Absorbance

     Use Wavelength(s):

          Single (nm):             430

          Evaluation:              SCA

          Calibrated at:           Date 8/19/2013  Time 10:22:56 AM
          Operator:                J

          Evaluation:              SCA

          Calibrated at:           Date 8/19/2013  Time 10:22:56 AM
          Operator:                J

          Weighting Method:        Least squares
          Calibration Curve:       C = k1 * A

          Analyte Name        Unit
          ------------------------
          Yellow              mg/l

Method Report Parameter:

     Information                        on
     Options/Info                       on
     Instruments/Acquisition            on
     Data Analysis                      on
     Report                             on
     Include Calibration Report         off

Calibration Report Parameter:

     Standard Spectra                   on
     Processed Standard Spectra         on
     Path Length Table                  on
     Data Analysis Parameters           on
     Coefficients                       on
     Used Wavelength Results            on
     Calibration Table(s) of Analytes   on
     Curve                              on
     Residual Spectra                   on
     Diagnostics                        on

Results Report Parameter:

     Sample Information                 off
     Sample Spectra                     off
     Processed Sample Spectra           off
     Used Wavelength Results            off
     Evaluation Results                 off
     Statistical Information            off
     Residual Spectra                   off
     Confirmation Results               off
     Include Method Report              off
     All Sample Spectra overlaid        on
     Summary                            on
     Evaluation Result Statistics       on

                          *** End Method Report ***
```

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