

AGILENT PREMIUM REFURBISHED INSTRUMENTS

THE TECHNOLOGY YOU WANT – AT A PRICE YOU CAN AFFORD

Product Overview

Dr. Fortin, SCF Pharma,
Quebec, Canada
gets excellent results
from his Agilent premium
refurbished LC

*"I already had
experience
with Agilent.
Agilent was
my first choice."*

Broad selection of premium refurbished instrumentation

Agilent Technologies' refurbished instrument program gives you the opportunity to purchase Agilent quality equipment at extremely affordable prices. Each instrument is fully refurbished to Agilent standards, including comprehensive electro-mechanical and operational performance validation.

We offer a broad selection of Agilent premium refurbished instruments. With all common analytical techniques readily available for delivery, we can also supply equipment for more complex applications.

The product overview represents what is typically available – for exact details of availability, contact your Agilent representative or visit our web site for contact information at www.agilent.com/chem/contactus



Product Category	Product Description	Available Series	Comments
Agilent HPLC and UHPLC Modules and Systems			
Solvent Delivery HPLC	Isocratic Pump refurbished	1100/1200	400 bar
	Binary Pump refurbished	1100/1200/1260	400 bar
	Quaternary Pump refurbished	1100/1200/1260	400 bar
	Capillary Pump refurbished	1100/1200/1260	400 bar
	Nanoflow Pump refurbished	1100/1200/1260	400 bar
	Preparative Pump refurbished	1200/1260	
Solvent Delivery UHPLC	Isocratic Pump refurbished	1260	600 bar
	Binary Pump refurbished	1200/1260/1290	600 bar
	Binary Pump refurbished	1290	1200 bar
	Quaternary Pump refurbished	1200/1260	600 bar
	Quaternary Pump refurbished	1290	1200 bar
Solvent Preparation	Standard Degassing Unit refurbished	1100/1200/1260	
	Micro Degassing Unit refurbished	1100/1200/1260	
	High Performance Degassing Unit refurbished	1260	
Sample Introduction HPLC	Autosampler refurbished	1100	400 bar; no thermostating capability; vials only
	Autosampler thermostatable refurbished	1100/1200	400 bar; vials only
	Well-plate autosampler refurbished	1100	400 bar
	Micro well-plate sampler refurbished	1100/1200	400 bar
	Preparative Autosampler refurbished	1100/1200/1260	400 bar
	Dual Loop Autosampler refurbished	1100/1200/1260	400 bar
	Thermostat refurbished	1100/1200/1290	
Sample Introduction UHPLC	High Performance Autosampler refurbished	1200	600 bar; vials and well-plates
	High Performance Autosampler refurbished	1260	600 bar; vials and well-plates
	High Performance Autosampler refurbished	1290	1200 bar; vials and well-plates
	Autosampler thermostatable refurbished	1260	600 bar
	Thermostat refurbished	1100/1200/1290	
Column Compartment	Thermostatted Column Compartment refurbished	1100/1200/1260	different valve options possible; max. temp 80°C
	Thermostatted Column Compartment refurbished	1290	different valve options possible; max. temp 80°C
UV Detection	Variable Wavelength Detector refurbished	1100/1200/1260/1290	different data rates from 14-160 Hz
	Multiple Wavelength Detector refurbished	1100/1200/1260	different data rates from 20-80 Hz
	Diode Array Detector refurbished	1100/1200/1260/1290	different data rates from 20-160 Hz
Special Detectors	Fluorescence Detector refurbished	1100/1200/1260	
	Refractive Index Detector refurbished	1100/1200/1260	
Valves	Valves refurbished	1100/1200/1260/1290	various types available
Fraction Collection	Analytical, preparative and micro scale refurbished	1100/1200/1260	
Handheld Controller	Instant Pilot refurbished	1200	
Integrated Systems	Isocratic and gradient LC Systems refurbished	1220	different configurations available

Product Category	Product Description	Available Series	Comments
Agilent LC/MS (Mass Spectrometry) Systems			
LC/MS instruments	Single Quadrupole refurbished	6100	
	Triple Quadrupole refurbished	6400	
	Time-of-Flight (TOF) refurbished	6200	
	Quadrupole Time-of-Flight (Q-TOF) refurbished	6500	
Ion Sources	Electrospray Ion Source (API) refurbished		
	Atmospheric Pressure Chemical Ionization (APCI) refurbished		
	Atmospheric Pressure Photo Ionization (APPI) refurbished		
	Multimode ESI/APCI refurbished		
	Dual-Spray Agilent Jet Stream (AJS) refurbished		
Agilent Capillary Electrophoresis			
CE Systems	Capillary Electrophoresis System refurbished	7100	
Agilent Gas Chromatography Systems			
GC Systems	GC Systems refurbished	6850/6890/7890	
	Various other accessories		different configurations available
Sample Injection	Autoinjector refurbished	7650	
	Autosampler refurbished	7693	
	Headspace Sampler refurbished	7694/7697	
Agilent GC/MS (Mass Spectrometry) Systems			
GC/MSD Systems	Quadrupole (GC/MSD) refurbished	5973N/5975/5977	different detectors and inlets available
	Triple Quadrupole (GC QQQ) refurbished	7000	different detectors and inlets available
	Quadrupole Time-of-Flight (GC Q-TOF) refurbished	7200	different detectors and inlets available
	Ion Trap refurbished	220/240	different detectors and inlets available

Please note:

This equipment is available on a first-come, first-served basis ONLY.

Learn more

www.agilent.com/chem/refurbished

Find an Agilent customer center

www.agilent.com/chem/contactus

USA and Canada

1-800-227-9770

agilent_inquiries@agilent.com

Europe

info_agilent@agilent.com

Asia Pacific

inquiry_lsca@agilent.com

This information is subject to change without notice.

© Agilent Technologies, Inc. 2013
Printed in the USA, August 1, 2013
5990-9570EN



Scan the QR code
with your smartphone
for more information



Agilent Technologies