# Nitrogen Generator NM32LA

# We have nothing to declare but our **Genius.**

Our new **Genius range** of Gas Generators is smoother, quieter, safer and more efficient than anything on the market. The **Genius range** offers superb technical performance in the lab and comes with our 'world class' after sales service as standard.

**PEAK SCIENTIFIC. GO WITH THE FLOW.** 



# Nitrogen Generator NM32LA

Features & Benefits:

Independent – No external compressor required

Intelligent - High duty and service indication

Quiet – Insulated compressor compartment for minimum disruption

Mobile – Flexibility to position generator where required, even under standard lab bench, the generator is supplied with caster wheels for easy mobility

Simple Installation – Generator designed as plug & play system

Economical – More cost effective than any other gas supply method

Convenient – Gas on demand, no health hazards, no need to worry about running out of gas

The NM32LA Nitrogen Generator is Peak Scientific's proven solution for customers requiring nitrogen for their LCMS application. Extensive testing of the generator by mass spec manufacturers throughout the development process guarantees the best nitrogen solution for your LCMS.

## **Technical Specifications**

Nitrogen	32 L/min @ 6.9 bar / 1.13 CFM @ 100 psi
Min/ Max Operating Temperature	5°C - 35°C / 41°F - 95°F
Max Relative Humidity	80% Non-Condensing
Max Altitude	2000 Metres
Particles	< 0.01µm
Gas Outlets	1x 1/4" BSP Female
Phthalates	None
Suspended Liquids	None
Noise Level	54 dB(A) @ 1m
Electrical Requirements	230v 50/60Hz 7.0A
Power Consumption	1610 watts
Dimensions (cm/ ins) WxDxH	60 x 75 x 71.2 / 23.6 x 29.5 x 28
Weight (kg/ lbs)	95 / 209

Ordering Information

Agilent Part Number:	G1953A
Peak Part Number:	10-6020
Annual Service Kit NM32LA	08-4780
Standard Maintenance Plan	09-3110
Complete Maintenance Plan	09-3000

#### For more Agilent information:

Find an Agilent customer center in your country: www.agilent.com/chem/contactus

**U.S. and Canada** 1-800-227-9770

agilent\_inquiries@agilent.com

Europe info\_agilent@agilent.com

Asia Pacific inquiry\_lsca@agilent.com



Agilent is a distributor of the Peak Nitrogen Generator

# Product Accreditations:











Peak Scientific UK

Tel: +44 (0)141 812 8100 Fax: +44 (0)141 812 8200

# Peak Scientific North America

Tel: +1 866 647 1649 Fax: +1 978 608 9503

## Peak Scientific India

Tel: +91 (0)40 2780 0663 Fax: +91 (0)40 2780 0663

## Peak Scientific China

Tel: +86 21 2890 9652 Fax: +86 21 2890 9999

#### Peak Scientific Germany

Tel: +49 (0)69 677 33 490 Fax: +49 (0)69 677 33 200

# Peak Scientific Africa

Tel: +27 (0)2185 16542 Fax: +27 (0)2185 16542

### Peak Scientific Mexico

Tel: +52 1 771 772 5757 Fax: +44 (0)141 812 8200