

EPA Method 504 by GC/MS

Column: DB-VRX

30 m x 0.25 mm I.D., 1.4 µm

J&W P/N: 122-1534

Oven: 40°C for 0.5 min

40-200°C at 10°/min

Injector: Splitless, 200°C

12 sec purge activation time

25 pg on-column of each component

Detector: MSD, 220°C transfer line, SIM

1. 1,1,1-Trichloroethane
 2. Carbon tetrachloride
 3. Trichloroethene
 4. Chloroform
 5. Tetrachloroethene
 6. Bromodichloromethane
 7. 1,2-Dibromopropane
 8. 1,2-Dibromoethane (EDB)
 9. Dibromochloromethane
 10. Bromoform
 11. 1,2-Dibromo-3-chloropropane (DBCP)
- Method 504 standards are found on pages 131

