

Less Stress. More Flexibility.

Agilent InfinityLab Flex Bench Family



Agilent
InfinityLab

The demand for flexibility in today's lab is high. The InfinityLab Flex Bench family provides easy ways for you to configure your instrument to suit your needs, anywhere in the lab. These three sturdy steel benches are all adjustable, enabling easy and safe access to instrument components, solvent bottles, pumps, columns, and accessories.

Agilent InfinityLab Flex Bench - the mobile solution

Move your instruments where you need them, when you need them. Whether you want to work next to your MS, or move analysis to new projects, the InfinityLab Flex Bench enables space optimization while taking the hassle out of instrument relocation.

Agilent InfinityLab Benchtop - the adjustable solution

Easily change modules and perform maintenance with your existing bench setup. The InfinityLab Benchtop fits into your existing lab setup, while providing an easy way to configure your LC system.

Agilent InfinityLab Flex Bench MS - the integrated solution

Optimize your whole lab setup. The InfinityLab Flex Bench MS allows you to keep your Agilent InfinityLab MSDs and Ultivo LC/MS systems in one stack with easy access to all areas. It improves lab environments with a gas cartridge holder, an integrated solution for waste management and a quiet cover to reduce pump noise.



Agilent InfinityLab Flex Bench



Agilent InfinityLab Benchtop



Agilent InfinityLab Flex Bench MS

Specifications and ordering Information

InfinityLab Flex Bench	Part Number
Flex Bench	5043-1252**
Flex Bench MS	G6015A***
Benchtop	5043-1711****
Flex Bench with power strip	5043-1759*
Flex Bench MS with power strip	G6015B*
Benchtop with power strip	5043-1740*
Accessories	
Shelf	5043-1750
Shelf (for a maximum 3 G1170A valve drives)	5043-1245
Drawer	5043-1735
Tray accessory kit	5043-1725
Waste bin for standard Flex Bench	5043-1278
Power cord	8121-2258
Flex Bench replacement hardware kit	5043-1289

* The power strip versions are only available in US, Canada, and Western Europe

** The Flex Bench includes the framework, four shelves, and one waste bin

*** The Flex Bench MS includes one shelf for InfinityLab MSD or Ultivo, three shelves for InfinityLab LC modules, one waste container with Stay Safe cap, one gas cartridge holder, plus the installation. Not included but strongly recommended is the Agilent G6011B Quiet Cover

**** The Benchtop includes the framework and three shelves

Agilent InfinityLab is an optimized portfolio of LC and LC/MS instruments, columns, and supplies designed to work together in perfect harmony.

Combined with Agilent OpenLab software and Agilent CrossLab services, it provides you with end-to-end solutions and support to make every day more productive.

To learn more about how to improve the efficiency of your LC workflow, visit: www.agilent.com/chem/infinitylab



Learn more at www.agilent.com/chem/flexbench

DE44304.4776736111

This information is subject to change without notice.

© Agilent Technologies, Inc. 2021
Published in the USA, April 19, 2021
5991-5163EN