

Benzodiazepines

Column: DB-35ms
30 m x 0.25 mm id, 0.25 μ m

P/N: 122-3832

Carrier: Helium at 31 cm/sec, measured at 50°C

Oven: 50°C for 0.5 min

150-340°C at 10°/min

340°C for 6 min

Injector: Splitless, 250°C

30 sec purge activation time

Detector: MSD, 280°C transfer line

full scan at m/z 40-400

1. Medazepam
2. Halazepam
3. Oxazepam
4. Lorazepam
5. Diazepam
6. Demoxepam
7. Desmethyldiazepam
8. Clobazam
9. Temazepam
10. Flunitrazepam
11. Delorazepam
12. Bromazepam
13. Prazepam
14. Flurazepam
15. Clonazepam
16. Alprazolam
17. Triazolam

