

# Constant EOF Independent of Buffer, pH from 2 to 11 (MD-General Purpose Kit)

## EOF of FC-SO<sub>3</sub> Surfactant Modified FC capillary

Column:  $\mu$ SIL-FC

50  $\mu$ m I.D. x 70 cm (50 cm to detector)

J&W P/N: 291-8111 (MD-Phosphate Kit)

Buffer: 10 mM phosphate, pH as indicated, +0.01% FC-SO<sub>3</sub> (J&W P/N: 590-4303)

Voltage: +20 kV

Injector: hydrostatic

Temperature: ambient

Detector: UV, 200 nm

Sample: 0.1% mesityl oxide

