

CE Separation of Globin Chains Without Urea

Column: μ SIL-FC
50 μ m I.D. x 75 cm (50 cm to detector)
J&W P/N: 291-8111 (MD-Phosphate Kit)
Buffer: 30 mM phosphate, pH 3.0
+ 0.05% of FC-PN (J&W P/N: 590-4301)
Voltage: 20 kV
Injector: hydrostatic
Temperature: ambient
Detector: UV, 200 nm
Sample: Diluted Hb standards (A, S, F, A₂) (Helena Lab)

