

## Organophosphorus Insecticides (EPA Method 1618)

**Column:** DB-1  
30 m x 0.53 mm I.D., 1.5 µm  
**J&W P/N:** 125-1032  
**Carrier:** Helium at 7 mL/min  
**Oven:** 100-250°C at 3°C/min  
250°C for 10 min  
**Injector:** Megabore Direct, 1 µL, 250°C  
**Detector:** FPD, 250°C, 526 nm filter

11. Ethyl parathion
12. Trichloronate
13. Merphos A
14. Tetrachlorvinphos
15. Fensulfothion
16. Sulprofos
17. EPN
18. Methyl azinphos
19. Ethyl azinphos
20. Coumaphos

1. Dichlorovos
2. Trichlorfon
3. Ethoprop
4. Phorate
5. Dimethoate
6. Terbufos
7. Diazinon
8. Methyl parathion
9. Ronnel
10. Malathion

