

## Anesthetics

**Column:** DB-1  
30 m x 0.53 mm I.D., 5.0  $\mu$ m  
**J&W P/N:** 125-1035  
**Carrier:** Helium at 6 mL/min  
**Oven:** 50°C for 3.75 min  
50-125°C at 30°/min  
**Injector:** Split 1:10, 5  $\mu$ L, 200°C  
**Detector:** FID, 250°C  
Nitrogen makeup gas at 30 mL/min

1. *iso*-Flurane
2. *en*-Flurane
3. Halothane
4. Methoxyflurane

