

## Clary Sage Oil

**Column:** DB-WAX  
30 m x 0.25 mm I.D., 0.25 µm  
**J&W P/N:** 122-7032  
**Carrier:** Helium at 30 cm/sec  
**Oven:** 70-180°C at 2°/min  
**Injector:** Split 1:100, 250°C  
**Detector:** FID, 300°C

1. Menthane
2. α-Pinene
3. Myrcene
4. β-Pinene
5. Limonene
6. γ-Terpinene
7. Camphor
8. Linalool
9. Linalyl acetate
10. Borneol
11. Carvone
12. Citronellol
13. Geraniol
14. Nerolidol

