

## CLP Pesticides

### Column: DB-XLB

30 m x 0.25 mm I.D., 0.25 µm  
**J&W P/N:** 122-1232  
 Carrier: Helium at 45 cm/sec, measured at 110°C  
 Oven: 110°C for 0.5 min  
 110-150°C at 25°/min  
 150-260°C at 12°/min  
 260-320°C at 15°/min  
 320°C for 2 min  
 Injector: Splitless, 250°C  
 30 sec purge activation time  
 50 pg composite 8081 standard,  
 and CLP-034 surrogate mix,  
 Accustandard Inc.  
 Detector: ECD, 350°C  
 Nitrogen makeup gas at 40 mL/min

- |                              |                             |
|------------------------------|-----------------------------|
| 1. Tetrachloro-m-xylene (SS) | 13. Dieldrin                |
| 2. α-BHC                     | 14. Endrin                  |
| 3. γ-BHC                     | 15. p, p'-DDD               |
| 4. β-BHC                     | 16. Endosulfan II           |
| 5. Heptachlor                | 17. p, p'-DDT               |
| 6. δ-BHC                     | 18. Endrin aldehyde         |
| 7. Aldrin                    | 19. Endosulfan sulfate      |
| 8. Heptachlor epoxide        | 20. Dibutylchlorendate (SS) |
| 9. γ-Chlordane               | 21. Methoxychlor            |
| 10. α-Chlordane              | 22. Endrin ketone           |
| 11. Endosulfan I             |                             |
| 12. p, p'-DDE                |                             |

SS - Surrogate Standard

Standards used were a composite of individual solutions supplied courtesy of Accustandard Inc., 25 Science Park, New Haven, CT 06511, 800-442-5290.

