



# Agilent Technologies

## **Agilent Bond Elut OMIX Tips**

A collection of citations to advance your research

### **Table of contents**

[Forensics and toxicology](#)

## Forensics and toxicology

### [Determination of Illicit Drugs in Urine and Plasma by Micro-SPE Followed by HPLC–MS/MS](#)

*Chromatographia*, **75**, 55-63 (2012)  
Sabino Napoletano *et al.*

**Tags**  
Bond Elut OMIX Tips, forensics and toxicology

#### **Abstract**

Agilent Bond Elut OMIX pipette tips were used in a micro-volume SPE procedure to extract drugs of abuse from urine and plasma. Published by Springer.

---

### [A \$\mu\$ -SPE procedure for the determination of cannabinoids and their metabolites in urine by LC–MS/MS](#)

*Journal of Pharmaceutical and Biomedical Analysis*, **91**, 169-175 (2014)  
Camilla Montesano *et al.*

**Tags**  
Bond Elut OMIX Tips, forensics and toxicology

#### **Abstract**

The authors report that extraction of cannabinoids required enzymatic hydrolysis for two hours using only 90  $\mu$ L of urine with Agilent Bond Elut OMIX C18 pipette tips, saving time, reducing enzyme requirements, and dramatically reducing costs. Published by Elsevier B.V.

---

[www.agilent.com/chem](http://www.agilent.com/chem)

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material. Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2014

Printed in the UK  
11 August, 2014

5991-4187EN