



# Agilent Seahorse XF Calibrant

p/n 100840-000 (500 mL), 103059-000 (100 mL)

## Data Sheet

XF Calibrant is optimized for XF sensor cartridge calibration, including XFe, XF, and XFp sensor cartridges. The calibrant is premixed and ready to use. XF Calibrant should be stored at room temperature. This material is not believed to be toxic, hazardous, or dangerous to health according to OSHA regulations.

Instructions for hydrating sensor cartridges with XF Calibrant can be found in the Basic Procedures section of our website: <http://www.agilent.com/en-us/products/cell-analysis-%28seahorse%29/basic-procedures-to-run-an-xf-assay>

### XF Calibrant Specifications

Volumes	500 mL, 100 mL
Storage	Room temperature
Expiration	2 years when properly stored

[www.agilent.com](http://www.agilent.com)

For Research Use Only. Not for use in diagnostic procedures.

This information is subject to change without notice.

© Agilent Technologies, Inc., 2017  
Published in the USA, March 21, 2017  
100840-503 Rev. B



**Agilent Technologies**