


High-quality DakoLink Cassettes for Tissue Processing



DakoLink Cassettes are high-density, reliable polymer cassettes which securely hold specimens during processing, embedding, and storage. They are specially designed to keep specimens safely submerged in liquids with flow-through slots that maximize fluid exchange and ensure proper drainage.

The unique design of the DakoLink Cassettes allows for easy, one-hand operation for processing tissue samples and biopsies. The one-piece integral lid eliminates the need for separate lids, can be opened and closed as often as necessary and relocks securely.

DakoLink Cassettes

- Highly resistant to solvents, stains, and harsh chemicals
- Special coating allows extremely durable, permanent, high-resolution text and barcodes
- Available in 3 styles (tissue, needle biopsy, and fine mesh) in an assortment of colors
- One-hand operation for ease of use
- One-piece lids lock securely

DakoLink Cassettes Specifications and Ordering Information

Dimensions: 1.11 in W x 1.58 in L x .23 in H

Materials: 45 degree coated face
Made from durable Delrin®

Quantity: 1440 cassettes to a box

Table 1. DakoLink Cassettes ordering information.

	Tissue Processing DL200	Needle Biopsy Processing DL201	Fine Mesh Processing DL202
White	DL200WE	DL201WE	DL202WE
Yellow	DL200YW	DL201YW	DL202YW
Pink	DL200PK	DL201PK	DL202PK
Aqua	DL200AA	DL201AA	DL202AA
Blue	DL200BE	DL201BE	DL202BE
Orange	DL200OE	DL201OE	DL202OE
Gray	DL200GY	DL201GY	DL202GY
Green	DL200GN	DL201GN	DL202GN
Lilac	DL200LC	DL201LC	DL202LC
Peach	DL200PH	DL201PH	DL202PH
Tan	DL200TN	DL201TN	DL202TN

www.agilent.com

This information is subject to change without notice.