

Nitrogen Generator NM32LA

We have nothing to declare but our **Genius.**

Our new **Genius range** of Gas Generators is smoother, quieter, safer and more efficient than anything on the market. The **Genius range** offers superb technical performance in the lab and comes with our 'world class' after sales service as standard.

PEAK SCIENTIFIC. GO WITH THE FLOW.



Nitrogen Generator NM32LA

Features & Benefits:

Independent – No external compressor required

Intelligent – High duty and service indication

Quiet – Insulated compressor compartment for minimum disruption

Mobile – Flexibility to position generator where required, even under standard lab bench, the generator is supplied with caster wheels for easy mobility

Simple Installation – Generator designed as plug & play system

Economical – More cost effective than any other gas supply method

Convenient – Gas on demand, no health hazards, no need to worry about running out of gas

The NM32LA Nitrogen Generator is Peak Scientific's proven solution for customers requiring nitrogen for their LCMS application. Extensive testing of the generator by mass spec manufacturers throughout the development process guarantees the best nitrogen solution for your LCMS.

Technical Specifications

Nitrogen	32 L/min @ 6.9 bar / 1.13 CFM @ 100 psi
Min/ Max Operating Temperature	5°C - 35°C / 41°F - 95°F
Max Relative Humidity	80% Non-Condensing
Max Altitude	2000 Metres
Particles	< 0.01µm
Gas Outlets	1x 1/4" BSP Female
Phthalates	None
Suspended Liquids	None
Noise Level	54 dB(A) @ 1m
Electrical Requirements	230v 50/60Hz 7.0A
Power Consumption	1610 watts
Dimensions (cm/ ins) WxDxH	60 x 75 x 71.2 / 23.6 x 29.5 x 28
Weight (kg/ lbs)	95 / 209

Ordering Information

Agilent Part Number:	G1953A
Peak Part Number:	10-6020
Annual Service Kit NM32LA	08-4780
Standard Maintenance Plan	09-3110
Complete Maintenance Plan	09-3000

Peak Scientific UK

Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America

Tel: +1 866 647 1649
Fax: +1 978 608 9503

Peak Scientific India

Tel: +91 (0)40 2780 0663
Fax: +91 (0)40 2780 0663

Peak Scientific China

Tel: +86 21 2890 9652
Fax: +86 21 2890 9999

Peak Scientific Germany

Tel: +49 (0)69 677 33 490
Fax: +49 (0)69 677 33 200

Peak Scientific Africa

Tel: +27 (0)2185 16542
Fax: +27 (0)2185 16542

Peak Scientific Mexico

Tel: +52 1 771 772 5757
Fax: +44 (0)141 812 8200

For more Agilent information:

Find an Agilent customer center in your country:
www.agilent.com/chem/contactus

U.S. and Canada	1-800-227-9770 agilent_inquiries@agilent.com
Europe	info_agilent@agilent.com
Asia Pacific	inquiry_lsca@agilent.com



Agilent Technologies

Agilent is a distributor of the Peak Nitrogen Generator

Product Accreditations:

