

The Table of NMR Frequencies* (By Mass)

Isotope	Spin	2.35 T	4.70 T	7.05 T	9.39 T	11.74 T	14.09 T	16.45T	17.61T	18.80T	21.14T	Natural abundance	Relative sensitivity	Absolute sensitivity
1 H	1/2	100.000	200.057	299.949	399.952	499.843	599.944	700.000	749.948	800.449	900.080	99.989	62.899	5877.739
2 H	1	15.351	30.710	46.044	61.395	76.729	92.095	107.454	115.122	122.874	138.168	0.015	0.228	0.003
3 H	1/2	106.664	213.389	319.938	426.605	533.152	639.924	746.648	799.924	853.791	960.061	0	73.330	0.000
3 He	-1/2	76.179	152.402	228.499	304.681	380.778	457.034	533.256	571.306	609.778	685.676	0.000	27.807	0.004
6 Li	1	14.716	29.441	44.141	58.857	73.557	88.288	103.013	110.363	117.795	132.457	7.59	0.200	1.422
7 Li	3/2	38.864	77.750	116.572	155.437	194.258	233.161	272.047	291.458	311.085	349.805	92.41	3.692	318.871
9 Be	-3/2	14.052	28.112	42.148	56.201	70.237	84.303	98.363	105.381	112.478	126.478	100.0	0.175	16.310
10 B	3	10.744	21.493	32.225	42.969	53.701	64.456	75.206	80.572	85.998	96.702	19.9	0.078	1.451
11 B	3/2	32.084	64.186	96.236	128.320	160.369	192.486	224.588	240.613	256.816	288.781	80.1	2.077	155.509
13 C	1/2	25.145	50.304	75.422	100.568	125.686	150.856	176.015	188.575	201.273	226.325	1.07	1.000	1.000
14 N	1	7.226	14.457	21.675	28.902	36.120	43.354	50.584	54.194	57.843	65.043	99.632	0.024	2.211
15 N	-1/2	10.137	20.279	30.405	40.542	50.668	60.815	70.957	76.020	81.140	91.239	0.368	0.066	0.023
17 O	-5/2	13.556	27.121	40.662	54.219	67.761	81.331	94.895	101.666	108.513	122.019	0.038	0.157	0.006
19 F	1/2	94.094	188.242	282.234	376.331	470.322	564.511	658.658	705.656	753.175	846.921	100.0	52.400	4897.181
21 Ne	-3/2	7.895	15.794	23.680	31.574	39.460	47.363	55.262	59.205	63.192	71.057	0.27	0.031	0.008
23 Na	3/2	26.452	52.919	79.342	105.795	132.218	158.697	185.163	198.376	211.734	238.088	100.0	1.164	108.801
25 Mg	-5/2	6.122	12.247	18.362	24.484	30.599	36.726	42.852	45.909	49.001	55.010	10	0.014	0.135
27 Al	5/2	26.057	52.129	78.157	104.215	130.244	156.327	182.398	195.413	208.572	234.533	100.0	1.113	103.999
29 Si	-1/2	19.867	39.746	59.591	79.459	99.305	119.192	139.070	148.994	159.027	178.821	4.683	0.493	2.159
31 P	1/2	40.481	80.985	121.422	161.904	202.340	242.862	283.365	303.585	324.028	364.359	100.0	4.172	389.947
33 S	3/2	7.676	15.356	23.024	30.700	38.368	46.052	53.732	57.566	61.443	69.090	0.76	0.028	0.020
35 Cl	3/2	9.798	19.601	29.389	39.187	48.974	58.782	68.586	73.479	78.427	88.189	75.78	0.059	4.190
37 Cl	3/2	8.156	16.316	24.463	32.619	40.766	48.930	57.090	61.164	65.283	73.408	24.22	0.034	0.772
39 K	3/2	4.666	9.336	13.997	18.663	23.325	27.996	32.665	34.996	37.352	42.002	93.258	0.006	0.557
40 K	-4	5.802	11.607	17.403	23.205	29.001	34.809	40.614	43.512	46.442	52.223	0.012	0.012	0.000
41 K	3/2	2.561	5.124	7.683	10.244	12.803	15.367	17.929	19.209	20.502	23.054	6.730	0.001	0.007
43 Ca	-7/2	6.730	13.464	20.187	26.917	33.639	40.376	47.110	50.471	53.870	60.575	0.135	0.019	0.002
45 Sc	7/2	24.292	48.597	72.863	97.155	121.420	145.737	170.042	182.175	194.443	218.645	100.0	0.902	84.263
47 Ti	-5/2	5.638	11.278	16.910	22.548	28.179	33.822	39.463	42.279	45.126	50.743	7.44	0.011	0.078
49 Ti	-7/2	5.639	11.281	16.914	22.554	28.187	33.831	39.474	42.290	45.138	50.756	5.41	0.011	0.057
50 V	6	9.970	19.946	29.906	39.876	49.836	59.816	69.792	74.772	79.807	89.741	0.25	0.062	0.015
51 V	7/2	26.303	52.621	78.895	105.199	131.474	157.803	184.121	197.259	210.542	236.748	99.75	1.145	106.705
53 Cr	-3/2	5.653	11.308	16.955	22.607	28.254	33.912	39.568	42.391	45.245	50.877	9.501	0.011	0.101
55 Mn	5/2	24.789	49.592	74.355	99.144	123.906	148.720	173.523	185.905	198.424	223.121	100.0	0.958	89.545
57 Fe	1/2	3.238	6.477	9.712	12.950	16.184	19.425	22.664	24.282	25.917	29.143	2.119	0.002	0.004
59 Co	7/2	23.727	47.468	71.169	94.897	118.598	142.349	166.089	177.941	189.923	213.563	100.0	0.840	78.522
61 Ni	-3/2	8.936	17.877	26.804	35.740	44.666	53.611	62.552	67.016	71.529	80.432	1.134	0.045	0.048
63 Cu	3/2	26.515	53.046	79.533	106.049	132.536	159.078	185.608	198.852	212.243	238.660	69.17	1.173	75.801
65 Cu	3/2	28.404	56.824	85.196	113.601	141.974	170.406	198.826	213.013	227.357	255.656	30.83	1.441	41.529
67 Zn	5/2	6.257	12.517	18.767	25.024	31.274	37.537	43.798	46.923	50.083	56.316	4.1	0.015	0.059
69 Ga	3/2	24.001	48.016	71.992	95.993	119.969	143.994	168.009	179.997	192.118	216.030	60.108	0.870	48.853
71 Ga	3/2	30.497	61.011	91.474	121.972	152.435	182.962	213.476	228.708	244.110	274.494	39.892	1.784	66.512
73 Ge	-9/2	3.488	6.979	10.463	13.952	17.436	20.928	24.418	26.161	27.922	31.398	7.73	0.003	0.019
75 As	3/2	17.123	34.255	51.359	68.483	85.587	102.727	119.859	128.411	137.059	154.118	100.0	0.316	29.510
77 Se	1/2	19.072	38.154	57.205	76.277	95.328	114.418	133.501	143.026	152.658	171.659	7.63	0.436	3.111
79 Br	3/2	25.054	50.123	75.151	100.206	125.233	150.313	175.381	187.895	200.548	225.510	50.69	0.989	46.864
81 Br	3/2	27.007	54.029	81.007	108.015	134.993	162.027	189.049	202.539	216.177	243.085	49.31	1.239	57.099
83 Kr	-9/2	3.848	7.697	11.541	15.389	19.232	23.083	26.933	28.855	30.798	34.631	11.49	0.004	0.038
85 Rb	5/2	9.655	19.316	28.961	38.616	48.261	57.926	67.586	72.409	77.285	86.904	72.17	0.057	3.818
87 Rb	3/2	32.721	65.461	98.147	130.869	163.555	196.309	229.049	245.392	261.917	294.517	27.83	2.204	57.314
87 Sr	-9/2	4.334	8.670	12.999	17.333	21.662	26.000	30.337	32.501	34.690	39.008	7	0.005	0.033
89 Y	-1/2	4.900	9.803	14.698	19.598	24.493	29.398	34.301	36.749	39.224	44.106	100.0	0.007	0.692
91 Zr	-5/2	9.296	18.598	27.884	37.181	46.467	55.773	65.074	69.717	74.412	83.674	11.22	0.050	0.527
93 Nb	9/2	24.476	48.966	73.416	97.893	122.343	146.843	171.333	183.559	195.919	220.305	100.0	0.922	86.197
95 Mo	5/2	6.157	12.317	18.468	24.625	30.775	36.938	43.099	46.174	49.283	55.417	15.92	0.017	0.259
97 Mo	-5/2	6.654	13.311	19.958	26.612	33.258	39.919	46.576	49.899	53.260	59.889	9.55	0.019	0.165
99 Ru	-3/2	4.605	9.213	13.813	18.418	23.019	27.628	32.236	34.536	36.862	41.450	12.76	0.006	0.073
101 Ru	-5/2	5.161	10.326	15.481	20.643	25.799	30.965	36.130	38.708	41.314	46.456	17.06	0.009	0.138
103 Rh	-1/2	3.186	6.375	9.558	12.744	15.927	19.117	22.305	23.897	25.506	28.681	100.0	0.002	0.190
105 Pd	-5/2	4.576	9.155	13.726	18.302	22.873	27.454	32.033	34.318	36.629	41.189	22.33	0.006	0.126
107 Ag	-1/2	4.048	8.098	12.141	16.189	20.233	24.285	28.335	30.357	32.401	36.434	51.839	0.004	0.202
109 Ag	-1/2	4.654	9.310	13.958	18.612	23.260	27.919	32.575	34.899	37.249	41.886	48.161	0.006	0.285
111 Cd	-1/2	21.215	42.443	63.636	84.852	106.044	127.281	148.508	159.105	169.819	190.956	12.8	0.601	7.185
113 Cd	-1/2	22.193	44.399	66.568	88.762	110.931	133.147	155.352	166.437	177.645	199.756	12.22	0.689	7.866
113 In	9/2	21.866	43.744	65.586	87.452	109.294	131.182	153.060	163.981	175.023	196.808	4.29	0.658	2.636

The Leader in Design. The Leader in Performance.

(Continued on page 2)



Agilent Technologies

The Table of NMR Frequencies* (By Mass) (Continued)

Isotope	Spin	2.35 T	4.70 T	7.05 T	9.39 T	11.74 T	14.09 T	16.45T	17.61T	18.80T	21.14T	Natural abundance	Relative sensitivity	Absolute sensitivity
115 In	9/2	21.913	43.838	65.726	87.640	109.528	131.463	153.388	164.333	175.399	197.230	95.71	0.662	59.196
115 Sn	-1/2	32.719	65.456	98.140	130.859	163.542	196.294	229.031	245.374	261.897	294.495	0.34	2.203	0.700
117 Sn	-1/2	35.632	71.285	106.879	142.512	178.105	213.774	249.426	267.223	285.218	320.719	7.68	2.846	20.424
119 Sn	-1/2	37.291	74.603	111.853	149.145	186.395	223.723	261.034	279.660	298.492	335.646	8.59	3.262	26.185
121 Sb	5/2	23.931	47.875	71.780	95.711	119.615	143.570	167.514	179.467	191.552	215.395	57.21	0.862	46.089
123 Sb	7/2	12.959	25.925	38.870	51.829	64.774	77.746	90.712	97.185	103.729	116.640	42.79	0.137	5.474
123 Te	-1/2	26.170	52.354	78.496	104.666	130.808	157.004	183.188	196.259	209.475	235.549	0.89	1.127	0.938
125 Te	-1/2	31.550	63.118	94.633	126.184	157.699	189.281	220.848	236.607	252.540	283.973	7.07	1.975	13.052
127 I	5/2	20.009	40.029	60.016	80.025	100.012	120.041	140.060	150.054	160.159	180.094	100.0	0.504	47.088
129 Xe	-1/2	27.810	55.636	83.416	111.227	139.007	166.846	194.671	208.562	222.606	250.314	26.44	1.353	33.433
131 Xe	3/2	8.245	16.494	24.730	32.976	41.211	49.465	57.714	61.832	65.996	74.211	21.18	0.035	0.698
133 Cs	7/2	13.116	26.240	39.342	52.459	65.560	78.690	91.813	98.365	104.989	118.056	100.0	0.142	13.264
135 Ba	3/2	9.934	19.874	29.798	39.733	49.656	59.601	69.541	74.503	79.520	89.418	6.592	0.062	0.380
137 Ba	3/2	11.113	22.232	33.333	44.446	55.547	66.671	77.790	83.341	88.953	100.025	11.23	0.086	0.906
138 La	5	13.194	26.396	39.576	52.771	65.951	79.158	92.360	98.950	105.613	118.759	0.09	0.144	0.012
139 La	7/2	14.126	28.259	42.370	56.496	70.606	84.746	98.879	105.935	113.068	127.142	99.91	0.177	16.553
141 Pr	5/2	30.620	61.257	91.844	122.465	153.052	183.703	214.340	229.634	245.097	275.604	100.0	1.806	168.763
143 Nd	-7/2	5.450	10.903	16.347	21.797	27.241	32.697	38.150	40.872	43.624	49.054	12.18	0.010	0.116
145 Nd	-7/2	3.360	6.722	10.078	13.438	16.795	20.158	23.520	25.198	26.895	30.243	8.3	0.002	0.019
147 Sm	-7/2	4.170	8.342	12.508	16.678	20.843	25.018	29.190	31.273	33.379	37.533	14.97	0.005	0.064
149 Sm	-7/2	3.440	6.882	10.318	13.758	17.195	20.638	24.080	25.798	27.535	30.963	13.83	0.003	0.033
151 Eu	5/2	24.860	49.734	74.567	99.428	124.261	149.146	174.020	186.437	198.992	223.760	47.82	0.966	43.189
153 Eu	5/2	10.980	21.966	32.934	43.915	54.883	65.874	76.860	82.344	87.889	98.829	52.18	0.083	4.060
155 Gd	-3/2	3.070	6.142	9.208	12.279	15.345	18.418	21.490	23.023	24.574	27.632	14.73	0.002	0.025
157 Gd	-3/2	4.030	8.062	12.088	16.118	20.144	24.178	28.210	30.223	32.258	36.273	15.65	0.004	0.060
159 Tb	3/2	24.040	48.094	72.108	96.148	120.162	144.227	168.280	180.287	192.428	216.379	100.0	0.874	81.670
161 Dy	-5/2	3.440	6.882	10.318	13.758	17.195	20.638	24.080	25.798	27.535	30.963	18.88	0.003	0.045
163 Dy	5/2	4.820	9.643	14.458	19.278	24.092	28.917	33.740	36.147	38.582	43.384	24.97	0.007	0.164
165 Ho	7/2	21.340	42.692	64.009	85.350	106.666	128.028	149.380	160.039	170.816	192.077	100.0	0.611	57.127
167 Er	-7/2	2.880	5.762	8.639	11.519	14.395	17.278	20.160	21.599	23.053	25.922	22.95	0.002	0.032
169 Tm	-1/2	8.290	16.585	24.866	33.156	41.437	49.735	58.030	62.171	66.357	74.617	100.0	0.036	3.350
171 Yb	1/2	17.499	35.009	52.489	69.989	87.469	104.986	122.495	131.236	140.073	157.508	14.31	0.337	4.508
173 Yb	-5/2	4.821	9.645	14.461	19.282	24.097	28.923	33.747	36.155	38.590	43.393	16.12	0.007	0.106
175 Lu	7/2	11.404	22.815	34.206	45.611	57.002	68.418	79.828	85.524	91.283	102.645	97.41	0.093	8.493
176 Lu	7	8.131	16.267	24.389	32.520	40.642	48.781	56.917	60.978	65.085	73.186	2.59	0.034	0.082
177 Hf	7/2	4.007	8.016	12.019	16.026	20.029	24.040	28.049	30.050	32.074	36.066	18.6	0.004	0.070
179 Hf	-9/2	2.517	5.035	7.550	10.067	12.581	15.101	17.619	18.876	20.147	22.655	13.62	0.001	0.013
181 Ta	7/2	11.990	23.986	35.963	47.953	59.929	71.931	83.927	89.916	95.971	107.916	99.988	0.108	10.131
183 W	1/2	4.166	8.335	12.497	16.664	20.825	24.996	29.165	31.246	33.350	37.501	14.31	0.005	0.061
185 Re	5/2	22.525	45.062	67.562	90.088	112.588	135.135	157.672	168.923	180.298	202.739	37.4	0.719	25.125
187 Re	5/2	22.752	45.516	68.243	90.995	113.722	136.497	159.261	170.625	182.115	204.783	62.6	0.741	43.338
187 Os	1/2	2.282	4.566	6.846	9.128	11.408	13.693	15.976	17.116	18.269	20.543	1.96	0.001	0.001
189 Os	3/2	7.765	15.535	23.292	31.058	38.815	46.588	54.358	58.236	62.158	69.895	16.15	0.029	0.445
191 Ir	3/2	1.718	3.437	5.153	6.871	8.587	10.307	12.026	12.884	13.752	15.463	37.3	0.000	0.011
193 Ir	3/2	1.871	3.743	5.612	7.483	9.352	11.225	13.097	14.032	14.976	16.840	62.7	0.000	0.024
195 Pt	1/2	21.497	43.006	64.479	85.977	107.450	128.969	150.477	161.215	172.071	193.488	33.832	0.625	19.756
197 Au	3/2	1.729	3.459	5.186	6.915	8.642	10.373	12.103	12.967	13.840	15.562	100.0	0.000	0.031
199 Hg	1/2	17.911	35.832	53.723	71.635	89.526	107.455	125.376	134.322	143.367	161.212	16.87	0.361	5.697
201 Hg	-3/2	6.611	13.227	19.831	26.442	33.047	39.665	46.280	49.582	52.921	59.508	13.18	0.018	0.224
203 Tl	1/2	57.123	114.279	171.340	228.465	285.526	342.707	399.862	428.394	457.242	514.154	29.524	11.724	323.499
205 Tl	1/2	57.684	115.401	173.022	230.708	288.329	346.071	403.787	432.599	461.730	519.201	70.476	12.073	795.177
207 Pb	1/2	20.921	41.853	62.751	83.672	104.570	125.512	146.444	156.894	167.459	188.302	22.1	0.576	11.895
209 Bi	9/2	16.069	32.148	48.200	64.269	80.321	96.407	112.485	120.511	128.626	144.636	100.0	0.261	24.393
235 U	7/2	1.841	3.684	5.523	7.365	9.204	11.047	12.890	13.810	14.739	16.574	0.72	0.000	0.000

*Based on "Provisional Recommendations for NMR Nomenclature: Nuclear Spin Properties and Conventions Chemical Shifts", Commission I.5, Molecular Structure and Spectroscopy, IUPAC, Nov. 2001

The Leader in Design. The Leader in Performance.

The Table of NMR Frequencies* (By Frequency)

Isotope	Spin	2.35 T	4.70 T	7.05 T	9.39 T	11.74 T	14.09 T	16.45T	17.61T	18.80T	21.14T	Natural abundance	Relative sensitivity	Absolute sensitivity
3 H	1/2	106.664	213.389	319.938	426.605	533.152	639.924	746.648	799.924	853.791	960.061	0	73.330	0.000
1 H	1/2	100.000	200.057	299.949	399.952	499.843	599.944	700.000	749.948	800.449	900.080	99.989	62.899	5877.739
19 F	1/2	94.094	188.242	282.234	376.331	470.322	564.511	658.658	705.656	753.175	846.921	100.0	52.400	4897.181
3 He	-1/2	76.179	152.402	228.499	304.681	380.778	457.034	533.256	571.306	609.778	685.676	0.000	27.807	0.004
205 Tl	1/2	57.684	115.401	173.022	230.708	288.329	346.071	403.787	432.599	461.730	519.201	70.476	12.073	795.177
203 Tl	1/2	57.123	114.279	171.340	228.465	285.526	342.707	399.862	428.394	457.242	514.154	29.524	11.724	323.499
31 P	1/2	40.481	80.985	121.422	161.904	202.340	242.862	283.365	303.585	324.028	364.359	100.0	4.172	389.947
7 Li	3/2	38.864	77.750	116.572	155.437	194.258	233.161	272.047	291.458	311.085	349.805	92.41	3.692	318.871
119 Sn	-1/2	37.291	74.603	111.853	149.145	186.395	223.723	261.034	279.660	298.492	335.646	8.59	3.262	26.185
117 Sn	-1/2	35.632	71.285	106.879	142.512	178.105	213.774	249.426	267.223	285.218	320.719	7.68	2.846	20.424
87 Rb	3/2	32.721	65.461	98.147	130.869	163.555	196.309	229.049	245.392	261.917	294.517	27.83	2.204	57.314
115 Sn	-1/2	32.719	65.456	98.140	130.859	163.542	196.294	229.031	245.374	261.897	294.495	0.34	2.203	0.700
11 B	3/2	32.084	64.186	96.236	128.320	160.369	192.486	224.588	240.613	256.816	288.781	80.1	2.077	155.509
125 Te	-1/2	31.550	63.118	94.633	126.184	157.699	189.281	220.848	236.607	252.540	283.973	7.07	1.975	13.052
141 Pr	5/2	30.620	61.257	91.844	122.465	153.052	183.703	214.340	229.634	245.097	275.604	100.0	1.806	168.763
71 Ga	3/2	30.497	61.011	91.474	121.972	152.435	182.962	213.476	228.708	244.110	274.494	39.892	1.784	66.512
65 Cu	3/2	28.404	56.824	85.196	113.601	141.974	170.406	198.826	213.013	227.357	255.656	30.83	1.441	41.529
129 Xe	-1/2	27.810	55.636	83.416	111.227	139.007	166.846	194.671	208.562	222.606	250.314	26.44	1.353	33.433
81 Br	3/2	27.007	54.029	81.007	108.015	134.993	162.027	189.049	202.539	216.177	243.085	49.31	1.239	57.099
63 Cu	3/2	26.515	53.046	79.533	106.049	132.536	159.078	185.608	198.852	212.243	238.660	69.17	1.173	75.801
23 Na	3/2	26.452	52.919	79.342	105.795	132.218	158.697	185.163	198.376	211.734	238.088	100.0	1.164	108.801
51 V	7/2	26.303	52.621	78.895	105.199	131.474	157.803	184.121	197.259	210.542	236.748	99.75	1.145	106.705
123 Te	-1/2	26.170	52.354	78.496	104.666	130.808	157.004	183.188	196.259	209.475	235.549	0.89	1.127	0.938
27 Al	5/2	26.057	52.129	78.157	104.215	130.244	156.327	182.398	195.413	208.572	234.533	100.0	1.113	103.999
13 C	1/2	25.145	50.304	75.422	100.568	125.686	150.856	176.015	188.575	201.273	226.325	1.07	1.000	1.000
79 Br	3/2	25.054	50.123	75.151	100.206	125.233	150.313	175.381	187.895	200.548	225.510	50.69	0.989	46.864
151 Eu	5/2	24.860	49.734	74.567	99.428	124.261	149.146	174.020	186.437	198.992	223.760	47.82	0.966	43.189
55 Mn	5/2	24.789	49.592	74.355	99.144	123.906	148.720	173.523	185.905	198.424	223.121	100.0	0.958	89.545
93 Nb	9/2	24.476	48.966	73.416	97.893	122.343	146.843	171.333	183.559	195.919	220.305	100.0	0.922	86.197
45 Sc	7/2	24.292	48.597	72.863	97.155	121.420	145.737	170.042	182.175	194.443	218.645	100.0	0.902	84.263
159 Tb	3/2	24.040	48.094	72.108	96.148	120.162	144.227	168.280	180.287	192.428	216.379	100.0	0.874	81.670
69 Ga	3/2	24.001	48.016	71.992	95.993	119.969	143.994	168.009	179.997	192.118	216.030	60.108	0.870	48.853
121 Sb	5/2	23.931	47.875	71.780	95.711	119.615	143.570	167.514	179.467	191.552	215.395	57.21	0.862	46.089
59 Co	7/2	23.727	47.468	71.169	94.897	118.598	142.349	166.089	177.941	189.923	213.563	100.0	0.840	78.522
187 Re	5/2	22.752	45.516	68.243	90.995	113.722	136.497	159.261	170.625	182.115	204.783	62.6	0.741	43.338
185 Re	5/2	22.525	45.062	67.562	90.088	112.588	135.135	157.672	168.923	180.298	202.739	37.4	0.719	25.125
113 Cd	-1/2	22.193	44.399	66.568	88.762	110.931	133.147	155.352	166.437	177.645	199.756	12.22	0.689	7.866
115 In	9/2	21.913	43.838	65.726	87.640	109.528	131.463	153.388	164.333	175.399	197.230	95.71	0.662	59.196
113 In	9/2	21.866	43.744	65.586	87.452	109.294	131.182	153.060	163.981	175.023	196.808	4.29	0.658	2.636
195 Pt	1/2	21.497	43.006	64.479	85.977	107.450	128.969	150.477	161.215	172.071	193.488	33.832	0.625	19.756
165 Ho	7/2	21.340	42.692	64.009	85.350	106.666	128.028	149.380	160.039	170.816	192.077	100.0	0.611	57.127
111 Cd	-1/2	21.215	42.443	63.636	84.852	106.044	127.281	148.508	159.105	169.819	190.956	12.8	0.601	7.185
207 Pb	1/2	20.921	41.853	62.751	83.672	104.570	125.512	146.444	156.894	167.459	188.302	22.1	0.576	11.895
127 I	5/2	20.009	40.029	60.016	80.025	100.012	120.041	140.060	150.054	160.159	180.094	100.0	0.504	47.088
29 Si	-1/2	19.867	39.746	59.591	79.459	99.305	119.192	139.070	148.994	159.027	178.821	4.683	0.493	2.159
77 Se	1/2	19.072	38.154	57.205	76.277	95.328	114.418	133.501	143.026	152.658	171.659	7.63	0.436	3.111
199 Hg	1/2	17.911	35.832	53.723	71.635	89.526	107.455	125.376	134.322	143.367	161.212	16.87	0.361	5.697
171 Yb	1/2	17.499	35.009	52.489	69.989	87.469	104.986	122.495	131.236	140.073	157.508	14.31	0.337	4.508
75 As	3/2	17.123	34.255	51.359	68.483	85.587	102.727	119.859	128.411	137.059	154.118	100.0	0.316	29.510
209 Bi	9/2	16.069	32.148	48.200	64.269	80.321	96.407	112.485	120.511	128.626	144.636	100.0	0.261	24.393
2 H	1	15.351	30.710	46.044	61.395	76.729	92.095	107.454	115.122	122.874	138.168	0.015	0.228	0.003
6 Li	1	14.716	29.441	44.141	58.857	73.557	88.288	103.013	110.363	117.795	132.457	7.59	0.200	1.422
139 La	7/2	14.126	28.259	42.370	56.496	70.606	84.746	98.879	105.935	113.068	127.142	99.91	0.177	16.553
9 Be	-3/2	14.052	28.112	42.148	56.201	70.237	84.303	98.363	105.381	112.478	126.478	100.0	0.175	16.310
17 O	-5/2	13.556	27.121	40.662	54.219	67.761	81.331	94.895	101.666	108.513	122.019	0.038	0.157	0.006
138 La	5	13.194	26.396	39.576	52.771	65.951	79.158	92.360	98.950	105.613	118.759	0.09	0.144	0.012
133 Cs	7/2	13.116	26.240	39.342	52.459	65.560	78.690	91.813	98.365	104.989	118.056	100.0	0.142	13.264
123 Sb	7/2	12.959	25.925	38.870	51.829	64.774	77.746	90.712	97.185	103.729	116.640	42.79	0.137	5.474
181 Ta	7/2	11.990	23.986	35.963	47.953	59.929	71.931	83.927	89.916	95.971	107.916	99.988	0.108	10.131
175 Lu	7/2	11.404	22.815	34.206	45.611	57.002	68.418	79.828	85.524	91.283	102.645	97.41	0.093	8.493
137 Ba	3/2	11.113	22.232	33.333	44.446	55.547	66.671	77.790	83.341	88.953	100.025	11.23	0.086	0.906
153 Eu	5/2	10.980	21.966	32.934	43.915	54.883	65.874	76.860	82.344	87.889	98.829	52.18	0.083	4.060
10 B	3	10.744	21.493	32.225	42.969	53.701	64.456	75.206	80.572	85.998	96.702	19.9	0.078	1.451
15 N	-1/2	10.137	20.279	30.405	40.542	50.668	60.815	70.957	76.020	81.140	91.239	0.368	0.066	0.023
50 V	6	9.970	19.946	29.906	39.876	49.836	59.816	69.792	74.772	79.807	89.741	0.25	0.062	0.015

The Leader in Design. The Leader in Performance.

(Continued on page 4)

The Table of NMR Frequencies* (By Frequency) (Continued)

Isotope	Spin	2.35 T	4.70 T	7.05 T	9.39 T	11.74 T	14.09 T	16.45T	17.61T	18.80T	21.14T	Natural abundance	Relative sensitivity	Absolute sensitivity
135 Ba	3/2	9.934	19.874	29.798	39.733	49.656	59.601	69.541	74.503	79.520	89.418	6.592	0.062	0.380
35 Cl	3/2	9.798	19.601	29.389	39.187	48.974	58.782	68.586	73.479	78.427	88.189	75.78	0.059	4.190
85 Rb	5/2	9.655	19.316	28.961	38.616	48.261	57.926	67.586	72.409	77.285	86.904	72.17	0.057	3.818
91 Zr	-5/2	9.296	18.598	27.884	37.181	46.467	55.773	65.074	69.717	74.412	83.674	11.22	0.050	0.527
61 Ni	-3/2	8.936	17.877	26.804	35.740	44.666	53.611	62.552	67.016	71.529	80.432	1.134	0.045	0.048
169 Tm	-1/2	8.290	16.585	24.866	33.156	41.437	49.735	58.030	62.171	66.357	74.617	100.0	0.036	3.350
131 Xe	3/2	8.245	16.494	24.730	32.976	41.211	49.465	57.714	61.832	65.996	74.211	21.18	0.035	0.698
37 Cl	3/2	8.156	16.316	24.463	32.619	40.766	48.930	57.090	61.164	65.283	73.408	24.22	0.034	0.772
176 Lu	7	8.131	16.267	24.389	32.520	40.642	48.781	56.917	60.978	65.085	73.186	2.59	0.034	0.082
21 Ne	-3/2	7.895	15.794	23.680	31.574	39.460	47.363	55.262	59.205	63.192	71.057	0.27	0.031	0.008
189 Os	3/2	7.765	15.535	23.292	31.058	38.815	46.588	54.358	58.236	62.158	69.895	16.15	0.029	0.445
33 S	3/2	7.676	15.356	23.024	30.700	38.368	46.052	53.732	57.566	61.443	69.090	0.76	0.028	0.020
14 N	1	7.226	14.457	21.675	28.902	36.120	43.354	50.584	54.194	57.843	65.043	99.632	0.024	2.211
43 Ca	-7/2	6.730	13.464	20.187	26.917	33.639	40.376	47.110	50.471	53.870	60.575	0.135	0.019	0.002
97 Mo	-5/2	6.654	13.311	19.958	26.612	33.258	39.919	46.576	49.899	53.260	59.889	9.55	0.019	0.165
201 Hg	-3/2	6.611	13.227	19.831	26.442	33.047	39.665	46.280	49.582	52.921	59.508	13.18	0.018	0.224
67 Zn	5/2	6.257	12.517	18.767	25.024	31.274	37.537	43.798	46.923	50.083	56.316	4.1	0.015	0.059
95 Mo	5/2	6.157	12.317	18.468	24.625	30.775	36.938	43.099	46.174	49.283	55.417	15.92	0.017	0.259
25 Mg	-5/2	6.122	12.247	18.362	24.484	30.599	36.726	42.852	45.909	49.001	55.010	10	0.014	0.135
40 K	-4	5.802	11.607	17.403	23.205	29.001	34.809	40.614	43.512	46.442	52.223	0.012	0.012	0.000
53 Cr	-3/2	5.653	11.308	16.955	22.607	28.254	33.912	39.568	42.391	45.245	50.877	9.501	0.011	0.101
49 Ti	-7/2	5.639	11.281	16.914	22.554	28.187	33.831	39.474	42.290	45.138	50.756	5.41	0.011	0.057
47 Ti	-5/2	5.638	11.278	16.910	22.548	28.179	33.822	39.463	42.279	45.126	50.743	7.44	0.011	0.078
143 Nd	-7/2	5.450	10.903	16.347	21.797	27.241	32.697	38.150	40.872	43.624	49.054	12.18	0.010	0.116
101 Ru	-5/2	5.161	10.326	15.481	20.643	25.799	30.965	36.130	38.708	41.314	46.456	17.06	0.009	0.138
89 Y	-1/2	4.900	9.803	14.698	19.598	24.493	29.398	34.301	36.749	39.224	44.106	100.0	0.007	0.692
173 Yb	-5/2	4.821	9.645	14.461	19.282	24.097	28.923	33.747	36.155	38.590	43.393	16.12	0.007	0.106
163 Dy	5/2	4.820	9.643	14.458	19.278	24.092	28.917	33.740	36.147	38.582	43.384	24.97	0.007	0.164
39 K	3/2	4.666	9.336	13.997	18.663	23.325	27.996	32.665	34.996	37.352	42.002	93.258	0.006	0.557
109 Ag	-1/2	4.654	9.310	13.958	18.612	23.260	27.919	32.575	34.899	37.249	41.886	48.161	0.006	0.285
99 Ru	-3/2	4.605	9.213	13.813	18.418	23.019	27.628	32.236	34.536	36.862	41.450	12.76	0.006	0.073
105 Pd	-5/2	4.576	9.155	13.726	18.302	22.873	27.454	32.033	34.318	36.629	41.189	22.33	0.006	0.126
87 Sr	-9/2	4.334	8.670	12.999	17.333	21.662	26.000	30.337	32.501	34.690	39.008	7	0.005	0.033
147 Sm	-7/2	4.170	8.342	12.508	16.678	20.843	25.018	29.190	31.273	33.379	37.533	14.97	0.005	0.064
183 W	1/2	4.166	8.335	12.497	16.664	20.825	24.996	29.165	31.246	33.350	37.501	14.31	0.005	0.061
107 Ag	-1/2	4.048	8.098	12.141	16.189	20.233	24.285	28.335	30.357	32.401	36.434	51.839	0.004	0.202
157 Gd	-3/2	4.030	8.062	12.088	16.118	20.144	24.178	28.210	30.223	32.258	36.273	15.65	0.004	0.060
177 Hf	7/2	4.007	8.016	12.019	16.026	20.029	24.040	28.049	30.050	32.074	36.066	18.6	0.004	0.070
83 Kr	-9/2	3.848	7.697	11.541	15.389	19.232	23.083	26.933	28.855	30.798	34.631	11.49	0.004	0.038
73 Ge	-9/2	3.488	6.979	10.463	13.952	17.436	20.928	24.418	26.161	27.922	31.398	7.73	0.003	0.019
149 Sm	-7/2	3.440	6.882	10.318	13.758	17.195	20.638	24.080	25.798	27.535	30.963	13.83	0.003	0.033
161 Dy	-5/2	3.440	6.882	10.318	13.758	17.195	20.638	24.080	25.798	27.535	30.963	18.88	0.003	0.045
145 Nd	-7/2	3.360	6.722	10.078	13.438	16.795	20.158	23.520	25.198	26.895	30.243	8.3	0.002	0.019
57 Fe	1/2	3.238	6.477	9.712	12.950	16.184	19.425	22.664	24.282	25.917	29.143	2.119	0.002	0.004
103 Rh	-1/2	3.186	6.375	9.558	12.744	15.927	19.117	22.305	23.897	25.506	28.681	100.0	0.002	0.190
155 Gd	-3/2	3.070	6.142	9.208	12.279	15.345	18.418	21.490	23.023	24.574	27.632	14.73	0.002	0.025
167 Er	-7/2	2.880	5.762	8.639	11.519	14.395	17.278	20.160	21.599	23.053	25.922	22.95	0.002	0.032
41 K	3/2	2.561	5.124	7.683	10.244	12.803	15.367	17.929	19.209	20.502	23.054	6.730	0.001	0.000
179 Hf	-9/2	2.517	5.035	7.550	10.067	12.581	15.101	17.619	18.876	20.147	22.655	13.62	0.001	0.013
187 Os	1/2	2.282	4.566	6.846	9.128	11.408	13.693	15.976	17.116	18.269	20.543	1.96	0.001	0.001
193 Ir	3/2	1.871	3.743	5.612	7.483	9.352	11.225	13.097	14.032	14.976	16.840	62.7	0.000	0.024
235 U	7/2	1.841	3.684	5.523	7.365	9.204	11.047	12.890	13.810	14.739	16.574	0.72	0.000	0.000
197 Au	3/2	1.729	3.459	5.186	6.915	8.642	10.373	12.103	12.967	13.840	15.562	100.0	0.000	0.031
191 Ir	3/2	1.718	3.437	5.153	6.871	8.587	10.307	12.026	12.884	13.752	15.463	37.3	0.000	0.011

*Based on "Provisional Recommendations for NMR Nomenclature:Nuclear Spin Properties and Conventions Chemical Shifts", Commission I.5, Molecular Structure and Spectroscopy, IUPAC, Nov. 2001

www.agilent.com

This information is subject to change without notice

© Agilent Technologies, Inc., 2011

Printed in USA, March 22, 2011

SI-0820

The Leader in Design. The Leader in Performance.



Agilent Technologies