


Agilent VWorks Automation Control Software 14.1



Which version do I need?

Compare the major features of VWorks 14.1 to previous versions.

VWorks Automation Control Software version features.

Features	VWorks 12.3	VWorks 13.1	VWorks 14.x Standard	VWorks 14.x Plus
Windows 10 (64 bit)	Win 7 (64 bit)	Yes	Yes	Yes
3rd party instrument support	Legacy	BioTek	BioTek*	BioTek*
Legacy Agilent instrument (Non-RoHS) support	Yes	Yes	Yes	Yes
Graphical interface forms authoring	Yes	Yes	Yes	Yes
Non-Windows admin access	No	No**	Yes	Yes
Integrated with Windows Active Directory	No	No	Yes	Yes
Records and data integrity	No	No	No	Yes
User accounts with permission control	Local	Local	Yes	Yes
Audit trails	Logs	Logs	Logs	Yes
Audit reports	No	No	No	Yes
E-signatures	No	No	No	Yes
Part number	Contact tech support		G5263AA	G5264AA

* Version 14.1 and later

** Workaround available

Support matrix

Different VWorks versions support Microsoft Windows and Agilent automation instruments.

VWorks	Features	Windows		Bravo		BenchCel		Lab Minihub		BenchBot	DDR	BioTek	
		10	7	XP (SP3)	Non-RoHS	RoHS	Non-RoHS	RoHS	Non-RoHS				RoHS
14.1	VWorks Standard version includes: – Integration with Agilent OpenLab Shared Services – Windows Active Directory login integration – Forms experiment database feature added – Additional compliance enabling features available in VWorks Plus	Yes (64 bit)	No	No	Yes	Yes	Yes	Yes	Yes	Yes	No	No	Yes
13.1	– Support of RoHS instruments – Removed 3rd party devices – Bravo NGS protocols validated for RoHS on 13.1.0.1366 – Last version supported on Windows 7	Yes (64 bit)	Yes (64 bit)	No	Yes	Yes Win 10 Only	Yes	Yes	Yes	Yes Win 10 Only	Yes	Yes	Yes 13.1.6 and later
13.0	– Windows 10 support – Limited support of 3rd party devices: Biotek washers, plate readers, thermo multidrop, liconic incubators	Yes (64 bit)	Yes (64 bit)	No	Yes	No	Yes	Yes	Yes	No	Yes	Yes	Yes
12.3	– Ability to run VWorks as Windows Non-Admin (Bravo only, not Biocels) – Biocels Support including 3rd party instruments		Yes (64 bit)	No	Yes	No	Yes	Yes	Yes	No	Yes	Yes	Yes
12.2	– Last version supported on WinXP (not SP1) – Last version to support 32 bit Windows (7 and XP)	Yes (64 bit)	Yes (not SP1)		Yes	No	Yes	Yes	Yes	No	Yes	Yes	Yes

www.agilent.com/chem/vworks-software-standard
www.agilent.com/chem/vworks-software-plus

DE78530146

This information is subject to change without notice.

© Agilent Technologies, Inc. 2021, 2022
 Printed in the USA, January 28, 2022
 5994-2893EN