

# Mass Spec/tacular!

The widest range of GC/MS solutions—  
from the name you can trust.

Another spectacular photograph from one of our employees: **Lee San Chung, Penang, Malaysia**

Whether you're looking to increase the productivity of your lab's routine screening, lower your detection limits, quantitate ultra-trace contaminants in very complex samples—or all of the above—be sure to look at Agilent's Mass Spec/tacular GC/MS lineup.

Our industry-leading GC/MSD and new Triple Quadrupole GC/MS solutions give you the analytical capabilities you need to keep pace with the most stringent new methods, plus day-after-day reliability to keep up with the most demanding sample loads.



#### **Agilent 5975C GC/MSD with Triple-Axis Detector**

Optimized performance from injection to final report—plus the industry's lowest LODs and LOQs

- Exceptional data – higher spectral quality and lower detection limits
- Synchronous SIM/Scan mode lets you capture SIM and full scan data in the same acquisition
- Fast, confident identification and quantitation with powerful Deconvolution Reporting Software



#### **Agilent 7000A Triple Quadrupole GC/MS**

Lowest limits of detection plus advanced, high-speed GC/MS/MS quantitation

- Routine femtogram-level sensitivity and superior selectivity—even in the most complex matrices
- Acquisition speeds of up to 500 MRM transitions per second let you confirm and quantitate more targets in a single run
- Easy-to-use MassHunter software simplifies data analysis, review and reporting

It's easy to choose the Agilent GC/MS or GC/MS/MS solution that exactly fits your lab's needs. Just call your Agilent representative or visit [www.agilent.com/chem/spectacularGCMS](http://www.agilent.com/chem/spectacularGCMS)

Download the complete information kit now.  
Visit [www.agilent.com/chem/spectacularGCMS](http://www.agilent.com/chem/spectacularGCMS)



5975C GC/MSD



7000A Triple Quadrupole GC/MS



7890A GC

Our measure is your success.



**Agilent Technologies**

## GC/MS and GC/MS/MS!





**Agilent Technologies**

Address, City, State, Zip

**Unprecedented sensitivity and performance.  
Day-after-day reliability**

Name  
Company  
Address  
City, State, Zip

