



Agilent Technologies

Agilent Technologies.
679 Springvale Road
Mulgrave, VIC, 3170
Australia

Ph: +61 (3) 9560 7133
Fax: +61 (3) 9560 7950
www.agilent.com

July 26, 2016

Dear Customer,

Thank you for purchasing Agilent Cary WinUV Pharma software that supports the Cary 60 UV-Vis spectrophotometers. We trust they fully met your expectations.

Agilent is committed to providing the best solutions for customers operating in regulated environments. As such we are contacting you to inform you that in two specific Cary WinUV Pharma releases (v. 5.0.0.1005 and v. 5.0.0.1008), the Configuration Manager username and password are stored in a text file in the Program directory of the PC where the Configuration Manager is installed. Our records indicate that the version of Cary WinUV 5.0 Pharma that you purchased belongs to one of these two releases.

To help your system maintain the highest levels of regulatory compliance we have created a new Configuration Manager executable file that replaces the existing file and removes the text files. We strongly recommend that you take this opportunity to update your software using the following link, which contains instructions and an updated Configuration Manager application:

<http://www.chem.agilent.com/Library/software/Public/ConfigManagerReplacement.exe>

Your continued satisfaction with Agilent products is always our priority and we sincerely regret any inconvenience that this may cause you. Please contact us if you have further questions.

Thank you again for choosing Agilent,

Dr. Kevin Grant
Molecular Spectroscopy Product Manager
Agilent Technologies