

These documents in this set show the compatibility of 1100/1200 series modules' firmware with

- previous module revisions,
- the Instant Pilot,
- the Control Module (up to B.04.02),
- the ChemStation,
- Certy,
- Picard/Tof software.

Document Name	For FW revisions (set)
FW_Compatibility_A061x_B061x	A.06.1x / B.06.1x/B.06.2x and B.02.xx (IP)
FW_Compatibility_A.0606_B0106	A.06.02/06 / B.01.02/06 and B.01.02/B.02.05 (IP)
FW_Compatibility_A0601	A.06.01 / B.01.01
FW_Compatibility_A511	A.05.11/A.05.13
FW_Compatibility_A509	A.05.09/A.05.10
FW_Compatibility_A506	A.05.06/A.05.08
FW_Compatibility_A505	A.05.05
FW_Compatibility_A504	A.05.01/A.05.04