

5975T LTM GC/MSD system consists of two major parts:

	5975 Series MSD	LTM GC
SW Tool for FW update	MSD Productivity Chemstation (starting from SW ver. E.02.02)	Agilent Instrument Utilities (starting from SW ver. B.01.05)
Latest FW revision	X.XX.XX	X.XX.XX.XXX

5975T.APCSetupInstaller.X.XX is the tool to install or remove AUX EPC on 5975T system. For detailed instruction please refer to the Aux EPC installation guide (shipped with 75T Aux EPC Kit accessory-G4367A or available at 75Tsuppot portal)