

**Solvent Compressibility Tables for  
 BINARY Pumps G4220A/B and G7120A  
 QUATERNARY Pumps G4204A and G7104A/C**

**Caution: Max length for solvent names is 30 characters!**

**NOTE:**

**The Quaternary pumps G4204A, G7104A and G7104C are only using a subset of the available solvent types, see table below.**

Pure Solvents or Solvent Mixtures	File Name
Water (BIN/QUAT)	Water V.03.XML
Methanol (BIN/QUAT)	Methanol V.03.XML
Acetonitrile (BIN/QUAT)	Acetonitrile V.03.XML
Ethanol (BIN/QUAT)	Ethanol V.03.XML
Methanol/ACN 1:1 (BIN/QUAT)	Methanol-ACN V.03.XML
Isopropanol (BIN/QUAT)	Isopropanol V.03.XML
Organic(BIN/QUAT)	Organic V.03.XML
Aqueous (BIN/QUAT)	Aqueous V.03.XML
Ethyl Acetate (BIN)	Ethyl Acetate V.01.XML
Heptane (BIN)	Heptane V.01.XML
Hexane (BIN)	Hexane V.01.XML
Isooctane (BIN)	Isooctane V.01.XML

Selectable Percentage of Solvent A in Solvent B	File Name
% Acetonitrile in Methanol (BIN)	Acetonitrile in Methanol V.01.XML
% Acetonitrile in Water (BIN/QUAT)	ACN in Water V.02.XML
% Methanol in Water (BIN)	Methanol in Water V.02.XML
% Isopropanol in Water (BIN/QUAT)	Isopropanol in Water V.01.XML
% Ethanol in Water (BIN/QUAT)	Ethanol in Water V.01.XML
% (Methanol/ACN 1:1) in Water (BIN/QUAT)	Methanol-ACN in Water V.01.XML
% Organic in Aqueous (BIN/QUAT)	Organic in Aqueous V.03.XML