




# Agilent Top Seller CE Kit

## Agilent Extended Light Path Capillaries


### Maximize your CE system's performance and experience unmatched sensitivity with Agilent Extended Light Path Capillaries!

Use Agilent Technologies Extended Light Path Capillaries ("bubble" cell capillaries) to improve your CE sensitivity 3-to 5-fold over standard capillaries.

With Extended Light Path Capillaries the inner diameter is increased only at the detection window, offering the sensitivity of a wide inner diameter capillary and the low current generation of a narrow one. Resolution is not sacrificed when used with Agilent Technologies' matching optical alignment interfaces.



CZE analysis of a tryptic digest of carbonic anhydrase using a standard capillary compared with the Agilent Extended Light Path Capillary



### The Agilent Top Seller CE kit (order number 5067-4621) includes:

- capillaries with 50 µm id and 150 µm path length, 2 each with 40, 56 and 72 cm effective length
- one alignment interface
- one free cybertool

In addition you receive a 16% discount.



Receive a **FREE Victorinox Cybertool** with the purchase of an Agilent Top Seller CE Kit of Extended Light Path Capillaries.

For more information contact your local sales representative or order online at:  
[www.agilent.com/chem/store](http://www.agilent.com/chem/store)

