

AGILENT ASSAYMAP SOLUTION ASSAYMAP SYRINGE TEST KIT

Ensure Optimal AssayMAP
Bravo Performance

Highlight

- Simple and rapid test to verify AssayMAP syringe performance



Confirm AssayMAP syringe performance

The Agilent AssayMAP Syringe Test Kit is designed to simultaneously test the integrity of all 96 syringes in the AM Bravo Head. The kit is simple to use and rapidly generates an Excel-based report that identifies any syringe that should be replaced. The AssayMAP Syringe Test Kit contains sufficient reagents to perform 20 tests.

Product description	Part number	Kit components
Agilent AssayMAP Syringe Test Kit	G5496-60050	Rack of resin free AssayMAP Cartridges AssayMAP Syringe Test Solution 2x 96 well Greiner flat bottom plates 2x 12 reservoir Seahorse plates

Required by not included in the kit	Part number	Description
Syringe Test Software		Included in AssayMAP Protein Sample Prep Workbench version 1.0.0.1301 or later.
Spectrophotometer		A plate reader capable of measuring absorbance at 425 nm and generating text or Microsoft Excel files.
96-Well Plates	Greiner 655101 Nunc 269620 Corning 9017	Labware options are listed. Two 96-well plates required per test if the labware is not reused. Labware sufficient for a single test is included in the kit.
Reservoir Plates	Seahorse 201254-100 Axygen RES-SW96-LP Agilent G5498B#049	Labware options are listed. Two 96-well plates required per test if the labware is not reused. Labware sufficient for a single test is included in the kit.

www.agilent.com/lifesciences/AssayMAP

This item is intended for Research Use Only.
Not for use in diagnostic procedures. This
information is subject to change without notice.

© Agilent Technologies, Inc. 2015
Published in the USA, August 14, 2015
5991-4587EN

