

THE POWER OF XENON: LIGHT YEARS AHEAD IN FLUORESCENCE MEASUREMENT TECHNOLOGY

The Cary Eclipse—a guaranteed instrument for accurate measurements

Agilent is proud to offer a 10 year guarantee on the unique Agilent Xenon flashlamp in every new Cary Eclipse fluorescence spectrophotometer. The Cary Eclipse is a robust, versatile spectrophotometer that costs virtually nothing to operate year to year. It has been purposely designed and built so that you get the accurate results that you are looking for in the quickest possible time—and all without instrument downtime due to maintenance or failure.

AGILENT
10 YEAR
LAMP LIFE
GUARANTEE



Did You Know?

The unique Xenon flashlamp in the Agilent Cary Eclipse delivers roomlight immunity, protects your samples from photodegradation, and now comes with a 10 year warranty. This high-brightness, high flash-frequency source preserves the integrity of your precious biological samples and increases lamp life by only flashing when a measurement is taken. It never needs warm up time and is central to the unique measurement technologies that ensure that you and your instrument can quickly adapt to changing analytical challenges.

With the Cary Eclipse you can minimize your instrument downtime, reduce your ongoing costs to virtually zero, and improve laboratory safety.

To learn more about the Agilent Cary Eclipse visit:
www.agilent.com/chem/cary-eclipse



Agilent Technologies

Agilent Cary Eclipse

Superior sensitivity where you need it most, renowned reliability, and enhanced sampling capabilities.



Save time and money

- 10 year lamp life guarantee
- Reduce ongoing lab costs
- Start work immediately with the zero warm up flashlamp



Simplify your analysis

- Streamline workflows
- Operate with covers open
- Easily access samples and accessories



Improve lab safety

Eliminate costly, time-consuming, and dangerous lamp replacements.



Flexible configurations

Use software to switch seamlessly between

- Fluorescence
- Phosphorescence
- Chemi/bio-luminescence
- Time-resolved modes



Uncompromised performance

- Proven sensitivity where you need it
- Rapid and precise kinetics measurements
- Accurate, stable, and reliable thermal control accessories



Save

With the power of Xenon you can save money and help care for the environment.



Want to learn more?

Ask for a demo and see how the proven choice for applications in academia, life sciences, biotechnology, chemicals, energy, and food testing analyses can exceed your analytical expectations.

To learn more about the Agilent Cary Eclipse visit:
www.agilent.com/chem/cary-eclipse

Information, descriptions and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2016
Published in USA, March 11, 2016
5991-6737EN