

ACCELERATE YOUR SUCCESS WHEN YOU UPGRADE YOUR MASSHUNTER SOFTWARE

A Faster Route to Insight

MassHunter Quantitative Analysis B.08.00

- Customize your direct-to-PDF reports with the new PDF Report Builder.
- Examine more samples with increased linear dynamic range for GC/Q-TOF data using SureMass.
- Measure the quality of your separation using metrics from the United States Pharmacopeia, Chinese Pharmacopeia, Japanese Pharmacopeia, and the European Pharmacopeia and implement Review-by-Exception on these and other Outliers.
- Extend, Customize and Integrate MassHunter Quantitative Analysis with our new Software Development Kit.
- Get support for Windows 10 in addition to Windows 7.



Learn more:
www.agilent.com/chem/masshunterinsights

MassHunter Qualitative Analysis B.08.00

- Increase productivity with a simplified user interface.
- Meet SANTE guidelines with GC/Q-TOF All Ions feature enhancements.
- Get support for Windows 10 in addition to Windows 7.

GC/MS Acquisition B.07.05

- Transform your GC experience with support for the new Intuvo 9000 GC.
- Speed up method creation and optimize instrument duty cycle using dMRM for GC/TQ.
- Expand your lab's capacity by seamlessly loading a GC/SQ method on a GC/TQ.
- Clean less and analyze more by adding JetClean to your existing GC/MSD or GC/TQ.
- Get support for Windows 10 in addition to Windows 7.

MassHunter Acquisition B.08.02 for LC/TQ and B.08.00 for LC/(Q-)TOF

- Capitalize on Agilent's trusted 1260 Infinity II LC and highly flexible 1290 Infinity II LC systems.
- Improve instrument performance with SWARM tuning for Agilent 6500 Series LC/Q-TOF systems.
- Get support for Windows 10 in addition to Windows 7.

MassHunter BioConfirm B.08.00

- Increase productivity with a simplified user interface.
- Easily create disulfide mapping and scrambling workflows.

MassHunter Walkup C.02.02

- Increase efficiency and reduce errors by creating your own sample information import workflow.
- Ensure reliable results and maximum uptime with elegant equilibration when switching chromatographic conditions.
- Put data precisely where needed with automated data file transfers and copying.

Faster, easier, and more productive
MS analysis from Agilent MassHunter
www.agilent.com/chem/masshunterinsights

Agilent Products are for Research Use Only.
Not for use in diagnostic procedures.
Information, descriptions and specifications in this
publication are subject to change without notice.

© Agilent Technologies, Inc. 2016
Published in USA, October 31, 2016
5991-7491EN