

Strengthen Your Regulatory Audit Process

Agilent OpenLab Data Management and MassHunter Networked Workstation 11.0



Compliant-ready software solution for LC/TOF and LC/Q-TOF systems

Preparing for regulatory audits can be challenging and stressful. Improve your audit experience by combining MassHunter with OpenLab data management and be confident that your lab is doing everything it can to satisfy requirements for US FDA Title 21 CFR Part 11 and related global regulations.

By pairing MassHunter with OpenLab Server/ECM XT for secure storage, your lab will benefit from extensive technical controls including:

- **Secure on-premises or cloud storage** of active records including data files, methods, worklists, report templates, and results
- **Audit trails** that clearly show the file history with details of who did what, when, and why
- **eSignatures** to ensure the appropriate staff members sign reports in the proper order
- Assign **access privileges** based on the function of the user and the data each can access

Make OpenLab part of your data integrity strategy

Foolproof disaster recovery

Enhance your disaster recovery plan and minimize the effects of operational disruption. Support business continuity by restoring and recovering important information in a matter of minutes.

Prevent fraud and identify unauthorized changes or deletions

Store data and files in a secure location with password-protected access. The software also tracks and saves every version, so you can recover previous versions and review the changes.

Ensure that all important data is securely stored and easy to find

Secure data from your other analytical systems too. Content import tools make certain that all important instrument data is automatically secured at the appropriate time. Flexible search functions ensure that your data will be easy to find.

Agilent
OpenLab



Secure storage for your lab

Together MassHunter with OpenLab Server/ECM XT provides the technical controls to securely acquire, process, report, and store data from LC/TOF and LC/Q-TOF systems.

Contact your Agilent representative, or visit:
www.agilent.com/chem/openlabecmxt

Learn more about our data management offerings:

www.agilent.com/chem/datamanagement

Download the OpenLab Server/ECM XT compliance white paper:

<https://www.agilent.com/cs/library/whitepaper/public/wp-cfr-part-11-compliance-openlab-server-ecm-xt-5994-3228en-agilent.pdf>

Download the MassHunter compliance white paper:

<https://www.agilent.com/cs/library/whitepaper/public/white-paper-regulatory-compliance-MassHunter-Workstation-5994-2902EN.pdf>

DE44407.7043518519

This information is subject to change without notice.

© Agilent Technologies, Inc. 2021
Published in the USA, August 1, 2021
5994-3807EN