### Agilent InfinityLab Quick Connect Fittings

Less stress, more reliable connections



# Save time and trouble with these key features



Truly finger-tight connection



Spring-loaded design



Reliable and reusable





Fast, easy connections

### Three simple steps, one perfect connection


Badly connected fittings are one of the top reasons for broad or tailing peaks, loss of resolution, and overall poor chromatography. With Agilent InfinityLab Quick Connect fittings, you can be sure you're getting a perfect column connection, every time.






The proprietary Agilent design features a spring-loaded mechanism for zero dead volume and a sure, tight connection to 1300 bar.

# Significantly reduce your costs and increase economic value



#### Agilent InfinityLab Quick Connect fittings can

deliver superior chromatographic performance and save you money:



Learn more about how Agilent InfinityLab Quick Connect fittings can significantly reduce your costs and improve performance: www.agilent.com/chem/infinitylab

This information is subject to change without notice.

© Agilent Technologies, Inc. 2019 Published in the USA, February 28, 2019 5994-0749EN

