

Performance Example




Figure 1: Chromatograms of THC, 11-OH-THC and THC-COOH at 0.25, 0.5 and 0.5 ng/ml respectively demonstrating rapid and sensitive analysis in seven minutes




Figure 2: Calibration lines for THC, 11-OH-THC and THC-COOH demonstrating excellent sensitivity and linearity

Key Benefits

- Rapid, high-throughput analysis of THC and its metabolites within five minutes.
- Sensitive and accurate quantification down to 0.25µg/l for THC and 0.5µg/l for THC-COOH and 11OH-THC.
- Two transitions per analyte provides confident identification and quantification.
- Easy and fast results interpretation through advanced software tools.

Learn more:

www.agilent.com/chem

Email:

info_agilent@agilent.com

Find a customer centre in your country:

www.agilent.com/chem/contactus

For Forensic Use.

Information is subject to change without notice.

© Agilent Technologies, Inc. 2010

April 21, 2010

5990-5784EE

