
Get More Done. Get Cary WinUV 5.3.

Upgrade your Agilent Cary 60 UV-Vis with the latest version of Cary WinUV





The <u>Agilent Cary 60 UV-Vis spectrophotometer</u> is a flexible, powerful, and reliable UV-Vis system, ideal for labs that want to get their work done. Add even more power to your Cary 60 UV-Vis by upgrading to the latest version of WinUV.

Agilent Cary WinUV software

The easy-to-use <u>Agilent Cary WinUV software</u> is built around the needs of routine analysts using the Agilent Cary 60 UV-Vis spectrophotometer.

What's new in Cary WinUV 5.3?

- Support for Microsoft Windows 11 operating systems.
- Up to 10 replicates and averages facilitate validation/verification as per ICH Q2 (R1).
- Automatically calculates measurement statistics streamlines method development.
- Updated 'Cary Help' and 'Cary Tutorial' modules speed up learning and simplify operating the Cary 60 UV-Vis.

Agilent Cary WinUV Pharma software

Cary WinUV Pharma Software supports the requirements of pharmaceutical laboratories operating in regulated environments:

- Secured data storage and user/access privilege management.
- Can be deployed on individual workstations, or in a client/server configuration.
- Data is stored in a secure Microsoft SQL Server Express database...

What's new in Cary WinUV Pharma 5.3?

- Support for Microsoft Windows 11 operating systems.
- Use of Microsoft SQL Server Express 2022 database with Microsoft support until 2028.

Special promotional discounts for Agilent Cary WinUV customers

Take advantage of this limited time offer:

Save up to 70% on software updates, and

30%

on software qualification services

Contact your local Agilent Representative for more information about the special upgrade promotion agilent.com/en/contact-us/page



This information is subject to change without notice.



