

Agilent InteroCyte Shadow-Free Imaging Microplates

Cell culture and imaging microplates for whole well label-free cell imaging assays




Shine the light on your label-free imaging assays


Easily eliminate meniscus shadows that can negatively affect the data from your whole well imaging assays

- **Illuminate your imaging assays:** unique well design eliminates meniscus shadows for more robust whole well imaging data
- **Versatile assays options:** an enhanced 96-well plate with standard and meniscus-free imaging capabilities for adherent and suspension cell applications
- **Setting the bar for precise cell analysis:** superior cell counting at the well's edge, perfect for uneven or low cell densities
- **Extra features, same price:** more functionality without extra cost

A representative image taken from the InteroCyte versus conventional microplate



4x image montage from InteroCyte plate



4x image montage from commonly used plate with column shaped well

Ordering information

| Description | Part Number |
|---|-------------|
| InteroCyte shadow-free imaging microplate, 96-well, tissue-culture treated, irradiated, with low evaporation lid, 9 per case | 204973-100 |
| InteroCyte shadow-free imaging microplate, 96-well, tissue-culture treated, irradiated, with low evaporation lid, 36 per case | 205070-100 |

| Dimensions (mm) | |
|-----------------------|---------------------------------|
| Plate dimensions | 127.76 x 85.47 (length x width) |
| A1 row offset | 14.38 |
| A1 column offset | 11.24 |
| Well spacing | 9.00 |
| Well top diameter | 8.10 |
| Well bottom diameter | 6.00 |
| Well bottom thickness | 1.02 |
| Plate height | 17.4 |

To learn more, or to place an order through the Agilent eStore, visit the following website or scan the QR code.
www.agilent.com/lifesciences/microplates



Related products

Agilent BioTek Cytation 5



Agilent BioTek Cytation 1



Agilent BioTek Lionheart FX



Agilent xCELLigence RTCA eSight



DE-007592

This information is subject to change without notice.