



Top Recommended Supplies for MP-AES

Easy-fit One Piece Torch

- The one-piece design offers the best plasma stability and ease of use
- The 1.2 mm ID injector and ball joint base ensures easy set-up
- Plug-and-play design automatically aligns torch and connects the gases for fast start-up and reproducible performance

G8003-70000 Torch, Easy-fit, for Agilent 4200/4210 MP-AES

Wavelength Calibration Standard

- Calibrates and verifies that your system is performing at the highest levels
- Ready-to-use without further dilution
- With the assurance of NIST traceability, this ensures accurate, reliable calibration with great value

6610030100 Bottle, ICP-OES wavelength calibration solution; 500 mL



Agilent OneNeb Series 2 Nebulizer

- Offers superior performance with MP-AES instruments
- Utilizing flow blurring nebulization, it provides improved sensitivity and greater tolerance to dissolved solids
- Inert construction makes it robust and compatible with most samples, including HF acid

G8010-60293 OneNeb series 2 inert concentric nebulizer

Peristaltic Pump Tubing

- Standard PVC tubing suitable for most aqueous/acidic samples
- Tubing in other materials also available for compatibility with high acids and organic solvents
- Tubing should be replaced regularly (every 1–2 weeks)

3710068100 Peristaltic pump tubing, PVC Solva Flex, orange/green tabs with flared ends, 15/pk
or
 3710067900 Peristaltic pump tubing, PVC Solva Flex, blue/blue, 15/pk

