


Limited Time Upgrade Offer

ICP-MS MassHunter Workstation Plus for 21 CFR Part 11 compliance



Now available

The Agilent ICP-MS MassHunter Workstation Plus (G7216AA) is now available. This standalone compliance system for ICP-MS operates on a single PC while facilitating compliance with FDA 21 CFR Part 11, EU Annex 11, and similar regulations in other countries. It is the perfect solution for laboratories seeking cost-effective compliance.

What is the difference between Workstation Plus and SDA?

The Agilent ICP-MS MassHunter Workstation Plus is the successor to the Spectroscopy Database Administrator (SDA) for ICP-MS, with SDA reaching the end of support period in 2027. It offers several significant advantages over the older SDA for ICP-MS system. The table on the right shows a comparison of features for data integrity.

Limited-time upgrade offer

If you wish to upgrade your current SDA for ICP-MS to the new Agilent ICP-MS MassHunter Workstation Plus, a limited-time upgrade offer is available.

Smooth transition

Agilent provides a migration tool and services to help you upgrade from SDA for ICP-MS to ICP-MS MassHunter Workstation Plus with minimal disruption. This ensures that your lab can continue operations using MassHunter Workstation Plus while old data on SDA is being migrated.

Please contact your [local Agilent office/representative](#) for further details.

For the details of ICP-MS MassHunter Workstation Plus, refer to the whitepaper "[Support for 21 CFR Part 11 and Annex 11 Compliance: ICP-MS MassHunter Workstation Plus Software](#)".

Table 1. Feature comparison between SDA for ICP-MS and ICP-MS MassHunter Workstation Plus

	Spectroscopy Database Administrator (SDA) for ICP-MS	ICP-MS MassHunter Workstation Plus
Audit Trail	●	●
System Activity Log	●	●
Access Control	●	●
E-Signature in MassHunter		●
Result Lock		●
Project Management (grouping configuration and data)		●
Secure Data Storage	●	●
Data Integrity Check (checksum verification)	●	●
Revision Control	●	●
Access to Past Revisions		●
Capacity	10 GB per database	Unlimited
Data Backup Location	Local	Local or on a Windows share
Supported on Windows 11		●

DE-002493

This information is subject to change without notice.

© Agilent Technologies, Inc. 2024
Published in the USA, December 18, 2024
5994-7956EN