


Shop for Agilent xCELLigence Products in Real Time


Purchase online through the Agilent Store*



Agilent is committed to making your life easier and your day more productive. No matter how you order products for your lab, the Agilent online store helps you save time and focus on more important tasks.

The Agilent xCELLigence Real-Time Cell Analyzer (RTCA) S16 is now available for purchase through the Agilent online store.*

Agilent RTCA instruments provide unparalleled insight into cellular function through impedance-based technology that uses gold biosensors to continuously monitor cell proliferation, morphological changes, and adherence in a label-free, real-time manner.



Time Minutes Hours Days

- Morphological changes
- Receptor activation
- Cell barrier function disruption
- Changes in cell numbers:
 - Proliferation
 - Cytotoxicity

Monitor live cells in real time. RTCA instruments are sensitive enough to capture cellular changes within minutes and robust enough to collect data for several days.

Key benefits

- **Simplified purchasing journey:** Buy any time, at your convenience.
- **Extra savings:** No Shipping and Handling charges for online orders over \$2,000, leading to savings up to 3% of the order value (\$3,000 maximum).

Learn more

www.agilent.com/chem/s16

Worldwide technical support
cellanalysis.support@agilent.com

USA and Canada

1-800-227-9770

Select option 3, then 8

* **Terms and conditions:** Promo only available with purchase online at www.agilent.com/chem/s16. Promotion limited to online purchase in US, Canada, UK, Ireland, Spain, France, and Singapore. Promotion cannot be combined with other offers.

The xCELLigence RTCA S16 is our most compact, entry-level instrument—and now it can be yours with a click of a button at a discounted rate.

- 16-well format suited for assay development and optimization
- Simply seed cells and start monitoring
- Operates in a standard CO₂ incubator for physiological conditions
- User-friendly software, including EC₅₀/IC₅₀ and % cytotoxicity plotting
- Highly sensitive, reproducible, and direct assays for immune cell killing, viral CPE, cytotoxicity, adhesion, cell signaling, barrier function, and more

Ordering information

Part Number	Product Description
380601310	xCELLigence RTCA S16 bundle including instrument, control unit, and software
Related Products	
5469830001	E-Plate 16 (1 × 6 plates)
5469813001	E-Plate 16 (16 × 6 plates)
6324738001	E-Plate VIEW 16 (1 × 6 plates)
6324746001	E-Plate VIEW 16 (6 × 6 plates)
300600890	E-Plate 16 PET (1 × 6 plates)
300600880	E-Plate 16 PET (6 × 6 plates)
6465382001	E-Plate insert S16 for co-culture (1 × 6 inserts)

Navigate your Agilent online experience

Get what you need on your terms on agilent.com, intuitively, with many time-saving features. From search to quote to purchase and beyond, we've got the roadmap to get you back to work.

My Account	Maintain your profile and view history and updates
My Quotes	Create your own quote with the "Create a Quote" option
Share	Forward your quote to colleagues within your organization
My Orders	Stay up to date on order status, tracking, and reorder information
My Service	Get service plans for instruments and renew expiring service plans or warranties
Product Availability and Promotion Acceptance	View live inventory updates and enable promotions at checkout
My Favorites	Save lists of frequently ordered supplies and parts

[Click here to order your RTCA S16 today.](#)

www.agilent.com/chem/s16

For Research Use Only. Not for use in diagnostic procedures.

RA44922.7411921296

This information is subject to change without notice.

© Agilent Technologies, Inc. 2023

Published in the USA, January 27, 2023

5994-5640EN