

安捷伦自动化解决方案移液吸头

选择指南

描述	货号	移液范围 (μL)	Bravo 移液头	适用的微孔板
 10 μL 吸头, 384 孔板架, 50 个/盒	10734-202	0.3-10	96ST, 384ST	96, 384, 1536
 10 μL 吸头, 无菌, 384 孔板架, 50 个/盒	10734-212	0.3-10	96ST, 384ST	96, 384, 1536
 10 μL 吸头, 导电, 384 孔板架, 50 个/盒	10734-302	0.3-10	96ST, 384ST	96, 384, 1536
 10 μL 吸头, 嵌套式, 导电, 384 孔板架, 90 个/盒	21740-202	0.3-10	96ST, 384ST	96, 384, 1536
 10 μL 吸头, 无菌, 导电, 384 孔板架, 50 个/盒	10734-312	0.3-10	96ST, 384ST	96, 384, 1536
 10 μL 吸头, 嵌套式, 无菌, 导电, 384 孔板架, 90 个/盒	21740-212	0.3-10	96ST, 384ST	96, 384, 1536
 30 μL 吸头, 384 孔板架, 50 个/盒	11484-202	0.5-30	96ST, 384ST	96, 384, 1536
 30 μL 吸头, 无菌, 384 孔板架, 50 个/盒	11484-212	0.5-30	96ST, 384ST	96, 384, 1536
 30 μL 吸头, 导电, 384 孔板架, 50 个/盒	11484-302	0.5-30	96ST, 384ST	96, 384, 1536
 30 μL 吸头, 嵌套式, 导电, 384 孔板架, 90 个/盒	21750-202	0.5-30	96ST, 384ST	96, 384, 1536
 30 μL 吸头, 无菌, 导电, 384 孔板架, 50 个/盒	11484-312	0.5-30	96ST, 384ST	96, 384, 1536
 30 μL 吸头, 嵌套式, 无菌, 导电, 384 孔板架, 90 个/盒	21750-212	0.5-30	96ST, 384ST	96, 384, 1536
 30 μL 吸头, 无菌, 带滤芯, 导电, 384 孔板架, 50 个/盒	11484-322	0.5-15	96ST, 384ST	96, 384, 1536
 70 μL 吸头, 384 孔板架, 50 个/盒	19133-102	0.75-70	96ST, 384ST	96, 384
 70 μL 吸头, 无菌, 384 孔板架, 50 个/盒	19133-112	0.75-70	96ST, 384ST	96, 384
 70 μL 吸头, 导电, 384 孔板架, 50 个/盒	19133-202	0.75-70	96ST, 384ST	96, 384
 70 μL 吸头, 无菌, 导电, 384 孔板架, 50 个/盒	19133-212	0.75-70	96ST, 384ST	96, 384
 70 μL 吸头, 无菌, 带滤芯, 导电, 384 孔板架, 50 个/盒	19133-142	0.75-50	96ST, 384ST	96, 384
 70 μL 吸头, 大孔, 无菌, 导电, 384 孔板架, 50 个/盒	19134-012	2.0-70	96ST, 384ST	仅 96
 70 μL 吸头, 大孔, 导电, 384 孔板架, 50 个/盒	19134-002	2.0-70	96ST, 384ST	仅 96
 70 μL 吸头, 大孔, 带滤芯, 无菌, 导电, 384 孔板架, 50 个/盒	19134-022	2.0-50	96ST, 384ST	仅 96
 250 μL 吸头, 96 孔板架, 50 个/盒	19477-002	2.0-250	96LT	96, 384
 250 μL 吸头, 无菌, 96 孔板架, 50 个/盒	19477-012	2.0-250	96LT	96, 384
 250 μL 吸头, 带滤芯, 无菌, 96 孔板架, 50 个/盒	19477-022	2.0-180	96LT	96, 384
 250 μL 吸头, 大孔, 96 孔板架, 50 个/盒	19477-032	5.0-250	96LT	仅 96
 250 μL 吸头, 大孔, 无菌, 96 孔板架, 50 个/盒	19477-072	5.0-250	96LT	仅 96
 250 μL 吸头, 大孔, 带滤芯, 无菌, 96 孔板架, 50 个/盒	19477-082	5.0-180	96LT	仅 96

查找当地的安捷伦客户服务中心

www.agilent.com/chem/contactus-cn

免费专线:

800-820-3278

400-820-3278 (手机用户)

联系我们:

LSCA-China_800@agilent.com

在线询价:

www.agilent.com/chem/erfq-cn

如需了解更多信息, 请访问:

www.agilent.com/chem/pipettetips

www.agilent.com

DE38642296

本文中的信息、说明和指标如有变更, 恕不另行通知。

© 安捷伦科技 (中国) 有限公司, 2016, 2020, 2022

2022 年 12 月 21 日, 中国出版

5991-1902ZHCN