



## Thermostatted Rectangular Cell Holder

**Part number:** 0010072200

Installation category I

Pollution degree 2

Safety class 1 (EN61010-1)

### Introduction

The Thermostatted Rectangular Cell Holder enables you to perform temperature controlled measurements for samples of low concentration. It requires, but does not include, a circulating water bath.

The instructions for installing the accessory are provided online with the Cary WinUV software. To view these instructions, start the Cary WinUV software, choose 'Help Topics' from the 'Help' menu and select 'Accessories' from the Cary WinUV Help home page. Click the 'Thermostatted rectangular cell holder' option from the list of Cary 50 accessories.

### Safety information

Before using this accessory, you must read the Safety Practices and Hazards section in your Cary operating instructions.

If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

### Specifications

For environmental and power supply requirements, refer to your Cary operating instructions.

#### Weights

Packed: 0.38 kg (0.84 lb)

Unpacked: 0.23 kg (0.5 lb)

#### Dimensions (W x D x H)

Packed: 185 x 190 x 95 mm (7.3 x 7.5 x 3.7 in)

Unpacked: 110 x 60 x 60 mm (4.3 x 2.4 x 2.4 in)



**Agilent Technologies**

## Thermostatted Rectangular Cell Holder

This information is subject to change without notice.



8510161500

Part Number: 8510161500

Edition 01/11

Issue 2

Printed in Australia

© Agilent Technologies, Inc. 1999, 2011

Agilent Technologies

679 Springvale Road

Mulgrave, VIC 3170