



Part number: 0210163190

Installation category I
Pollution degree 2
Safety class 3 (EN61010-1)

Introduction

The Extended Sample Compartment increases the size of the Cary sample compartment, enabling you to use accessories such as the 18 Cell Holder.

The instructions for installing the accessory are provided online with the Cary WinUV software. To view these instructions, start the Cary WinUV software, choose 'Help Topics' from the 'Help' menu and select 'Accessories' from the Cary WinUV Help home page. Click the 'Extended Sample Compartment' option from the list of Cary 50 accessories.

Safety information

Before using this accessory, you must read the Safety Practices and Hazards section in your Cary operating instructions.

If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

Specifications

For environmental requirements, refer to your Cary operating instructions.

Weights

Packed: 1.2 kg (2.6 lb)
Unpacked: 0.8 kg (1.8 lb)

Dimensions (W x D x H)

Packed: 290 x 200 x 190 mm (11.4 x 7.9 x 7.5 in)
Unpacked: 260 x 190 x 170 mm (10.2 x 7.5 x 6.7 in)

This information is subject to change without notice.



8510164900

Part Number: 8510164900

Edition 01/11

Issue 2

Printed in Australia

© Agilent Technologies, Inc. 1999, 2011

Agilent Technologies

679 Springvale Road

Mulgrave, VIC 3170