

AriaMx Real-Time PCR System

TOTAL CONFIDENCE qPCR



SPEED

AGILITY

PRECISION

Know More 

The following fixes have been addressed in Aria Software Version 1.3:

- Installer Unable to Install AriaDx. The Config didn't updated after install Dx.
- No Image on Splash Screen Defect: 353959
- Install Shield: CFR should auto selected in Dx mode
- Install Shield: Dx Mode Install Folder Path Should be 'Agilent Aria' include its contents
- Install Shield: On Dx mode, mode labeling didn't updated
- Help: AriaMx Exist on AriaDx SW Help Window
- AriaMx/Dx: Version Number didn't updated
- IQ Test: Dx Mode Default Save Path contain Mx
- Instrument Filtering: Mx/Dx SW able to view Mx/Dx instrument
- Import/Export: Mx & Dx Application able to open cheating experiment file
- Mx SW: Application always prompt for license when launched
- Generate Report - Mx / Dx Standard Mode / CFR Mode
- Mx/Dx Generate report: Experiment Name truncated
- Install Shield: Unable to install on New Empty Machine
- CFR Mode: Application able to view, open, run, generate report for Mx and Dx experiment which is on same database
- Mx/Dx Standard Mode: My Templates Tab refer to wrong Template Path Mx/Dx Standard Mode: Menu File -> Print (Instrument name & Instrument Type column label missing)
- Dx CFR: Create Archive Database Window unable to connect Database
- Dx CFR: Menu -> Print Throw Exception Message
- Dx CFR Mode: Add Database Feature able to add Mx database instance as Primary database
- Dx: Instrument File Explorer Display Mx Files But Not Dx Files
- Mx/Dx: Instrument File Explorer Default Copy Location is Wrong
- Mx/Dx CFR Mode: Menu File -> Print (Database Host column label missing)