

# Agilent 700 시리즈 축방향 ICP-OES 소모품

## 애질런트 기기를 위한 애질런트 소모품

애질런트는 실험실의 생산성을 최적화하기 위해 최선을 다하고 있습니다. 이 참조 안내서는 700 시리즈 축방(수평 방향) ICP-OES 시스템에 대해 가장 흔히 주문하는 소모품 및 부품을 보여줍니다. 이 리스트를 보유하면 필요한 소모품을 빠르게 찾아볼 수 있고 기기 가동 중단 시간을 최소화할 수 있습니다.



펌프 및 연결 튜빙	부품 번호
Tubing, PVC, white/white, 12/pk	3710034400
Tubing, PVC, blue/blue, 12/pk	3710034600
Tubing, PVC, black/black, 12/pk	3710027200
Tubing, PVC, orange/white, 12/pk	3710046900
Y-piece connector, for online addition of internal standard/ionization buffer	1610132400
Replacement tubing, for ICP-OES aqueous samples, spray chamber to torch, per m	3710033400
Replacement tubing, for ICP-OES aqueous samples, spray chamber to peristaltic pump waste line (glass, cyclonic spray chamber), per m	3710024600
Replacement tubing, spray chamber to peristaltic pump waste line (Sturman-Masters spray chamber), per m	3710033900

스프레이 챔버	부품 번호
Spray chamber, single-pass, glass cyclonic	2010081700
Spray chamber, double-pass, glass cyclonic	7910043700
Sturman-Masters spray chamber, inert, white, for HF samples or when using the V-groove nebulizer for fusions or slurries. Requires mounting bracket p/n 410328600	110593190
Transfer tube, FEP, for axial ICP	7910051590
O-ring, 1/4 in id x 3/8 in od x 1/16 in, nitrile, for nebulizer inlet	6910001200
O-ring kit, for cyclonic glass spray chamber	9910081600

예비 튜빙 키트	부품 번호
Sample compartment/torch box tube connection kit, includes all tubing for spray chamber area and spare nipples for drain (excluding pump tubing)	9910057100
Tubing and connectors kit, includes 2/pk standard peristaltic pump tubing, snap connectors, other connectors used in sample compartment, and torch box area, ideal for spare startup kit	9910097500
Internal standard kit 1, with black/black tubing, includes pump tubing, Ezy-Fit connectors, 5 Y-pieces, and connecting tubing	9910124000
Internal standard kit 2, with orange/white tubing, includes pump tubing, Ezy-Fit connectors, 5 Y-pieces and, connecting tubing	9910124100
Tubing kit, aqueous, for SeaSpray/K-style nebulizers and glass, cyclonic spray chamber	9910123600
Tubing kit, aqueous, for V-groove nebulizer and Sturman-Masters spray chamber	9910123800

\* Each aqueous tubing kit includes standard PVC pump tubing, FEP transfer tube, drain tubing, waste line, EzyFit connectors, nebulizer capillary tubing, connectors, and connecting tubing

Tubing kit, organic, for SeaSpray/K-style nebulizers and glass, cyclonic spray chamber	9910124200
Tubing kit, organic, for V-groove nebulizer and Sturman-Masters spray chamber.	9910124400

\* Each organic tubing kit includes PVC SolvaFlex pump tubing for use with low volatile organic solvents, FEP transfer tube, drain tubing, waste line, EzyFit connectors, nebulizer capillary tubing, connectors, and connecting tubing

Nebulizer 및 소모품	부품 번호
Concentric glass nebulizer, SeaSpray/Ezylok, includes sample and gas inlet fittings	2010096400
Concentric glass nebulizer, K-style (Conikal for organic solvents)	2010106800
OneNeb series 2 inert concentric type nebulizer. For use with HF and high TDS solutions.	G8010-60293
V-groove nebulizer kit, inert, for Sturman-Masters spray chamber	9910057400
Sample inlet fitting, EzyFit, for concentric glass nebulizer	9910107700
Gas inlet connector, Ezylok, for concentric glass nebulizer	9910127800
Drain tubing, per m, 1/8 in id x 1/4 in od, connects argon gas supply to nebulizer and drain tubing for SPS 3	3710031500

토치 및 소모품	부품 번호
Torch, low flow, quartz, with 2.4 mm id injector (standard torch)	2010090400
Torch, low flow, quartz, with 1.4 mm id injector (for organics)	2010104800
Torch, low flow, quartz with 0.8 mm id injector (for volatile organics)	2010104700
Torch, high solids, with tapered 2.4 mm id injector (for high TDS samples)	2010094800
GazFit fittings for one piece and semidemountable ICP-OES torches, 4/pk	9910107100
Silicon tubing, 3/16 in id x 5/16 in od, per m, connects plasma and auxiliary gas lines to the argon outlet	2410023800

Pre-optics 콘	부품 번호
Pre-optics cone for Agilent Vista/700 series axial ICP-OES	3610007300
Estele cleaning powder, for cleaning interface cone	8110002800

SPS 3 자동 시료 주입기	부품 번호
Inert probe, 0.8 mm id, PTFE-sleeved	9910111900
Inert probe, 1.0 mm, PTFE-sleeved, used with diluter	9910112000
Inert probe, 1.3 mm id, PTFE-sleeved	9910130900
Peristaltic pump tubing, 3-bridged, gray/gray, 12/pk	3710049000
Silicon tubing, 1 mm id, per m, to connect to nebulizer	3710026400
Drain tubing, per m, 1/8 in id x 1/4 in od, connects argon gas supply to nebulizer and drain tubing for SPS 3	3710031500
Polypropylene test tubes, 17 mm od, 16.5 mL, 1,000/pk	00000001600L
Polypropylene tubes, 16 mm od, 1,000/pk	190049700
Digestion tubes, polypropylene, 5 mL graduations to 50 mL, 30 mm od x 115 mm, 500/pk 500/pk	190047900
Rack overlay kit, 40-position sample rack for Agilent SPS 3 autosampler, allows 40 x 18 mm od tubes to be used with this rack	6610026000
Rack overlay kit, 60-position sample rack for Agilent SPS 3 autosampler, allows 60 x 13 mm od tubes to be used with this rack	6610026100

SVS 1 생산성 패키징용 예비 부품	부품 번호
Valve, 4-port, for SVS. Replace every four to six months, depending on use.	3610010100
SVS, EzyFit, for concentric glass nebulize, 10/pk	1610136400
Flangeless fitting, 1/pk	1610136500
Ferrule, 1/pk	1610136600
FEP capillary tubing, 0.0625 in od x 0.02 in id	3710052400

VGA 77	부품 번호
Sodium tetrahydroborate reagent, 500 g, hazardous	8210029100
Peristaltic pump tubing, PVC, for pumping reagents, black/black, 12/pk	3710027200
Peristaltic pump tubing, for pumping reagents, purple/black, 12/pk	3710027300
Tubing and connector kit	9910061900
ICP-OES hydride module	9910062200
Gas-liquid separator	9910062000
Pump tube beds kit	9910050700

소모품 키트	부품 번호
Matrix Valu-Pack for axial ICP. Includes high solids and semidismountable torch, extra quartz injectors, 2 spray chambers (glass cyclonic and inert Sturman-Masters), V-groove nebulizer, O-ring kit, tubing/connectors, and EzyFit/GazFit fittings.	9910112600
Extenda Valu-Pack for axial ICP. Includes 7 torches, 3 glass cyclonic spray chambers (1 double-pass and 2 single-pass), SeaSpray nebulizer, 2 O-ring kits, tubing and connectors kit, and EzyFit/GazFit fittings.	9910112800

다원소 검량 표준물질	부품 번호
Calibration mix 1, 125 mL. Contains 100 mg/L Sb, Mo, Sn, and Tl in 2% HNO <sub>3</sub> plus 0.5% HF.	6610030500
Calibration mix 2, 125 mL. Contains 100 mg/L Ag, Al, As, Ba, Be, Cd, Co, Cr, Cu, Mn, Ni, Pb, Se, Tl, Th, U, V, and Zn in 5% HNO <sub>3</sub> .	6610030600
Calibration mix majors, 125 mL. Contains 500 mg/L Ca, Fe, K, Mg, and Na in 5% HNO <sub>3</sub> .	6610030700
Calibration mix 3, organic, 100 mL. Contains 500 ppm Ag, Al, B, Ba, Ca, Cd, Cr, Cu, Fe, Mg, Mn, Mo, Na, Ni, P, Pb, Si, Sn, Ti, V, and Zn in mineral oil.	6610030200

파장 검량 용액	부품 번호
ICP-OES wavelength calibration concentrate, 500 mL. Contains 50 mg/L Al, As, Ba, Cd, Co, Cr, Cu, Mn, Mo, Ni, Pb, Se, Sr, Zn, and 500 mg/L K in 5% HNO <sub>3</sub> . Dilute x10 before use.	6610030000
ICP-OES wavelength calibration solution, 500 mL. Contains 5 mg/L Al, As, Ba, Cd, Co, Cr, Cu, Mn, Mo, Ni, Pb, Se, Sr, Zn, and 50 mg/L K in 5% HNO <sub>3</sub> . Ready to use.	6610030100
Calibration blank solution for ICP-OES, 500 mL, 0 ppm. Contains pure ASTM Type 1A water with 5% HNO <sub>3</sub> for use as a calibration blank or dilution of the wavelength calibration concentrate. Ready to use.	5190-7001

내부 표준물질/이온화 억제제 500mL 병으로 제공	부품 번호
ICP internal standard, 125 mL bottle, 100 mg/L <sup>6</sup> Li, Sc, Y, In, Tb, Bi in 5% HNO <sub>3</sub>	6610030400
Cesium buffer solution Contains 1% Cs (from carbonate) in 5% nitric acid, 500 mL bottle	5190-8343
Yttrium solution – 1% 10,000 mg/L in 5% nitric acid. Also used for instrument optimization. 500 mL bottle	5190-8233
Multi-element internal standard (mix 4) Contains 50 µg/mL <sup>6</sup> Li, Sc, 25 µg/mL Ge, Te, 10 µg/mL Bi, In, Tb in 2% HNO <sub>3</sub> with trace HF. 100 mL bottle	5190-8593

애질런트는 원자 분광기용 단일 원소 및 다원소 표준물질, 이온화 버퍼, 유기 금속 및 바이오디젤 표준물질, 용매 및 환경 표준물질을 제공합니다. NIST 추적성을 보장한 동시에, 애질런트의 인증 표준물질 및 표준 혼합물은 정확하고 신뢰할 수 있는 높은 가치의 검량을 제공합니다.

## 애질런트 분광기 소모품- 귀하가 원하는 결과를 얻을 수 있도록 도와드립니다

보다 확실한 분광기 결과를 위해 고품질의 애질런트 소모품으로 귀하의 전체 시스템을 최적화할 수 있습니다.

당사의 소모품, 키트, 액세서리 등은 ISO 9001의 인증을 받은 품질 관리 시스템에서 애질런트의 엄격한 기준에 맞춰 설계, 생산, 테스트됩니다.

왜 분석 결과에 영향주는 위험을 감수하며 정품 애질런트 소모품보다 품질이 낮은 제품을 선택하십니까?

중공 음극관 램프와 AA, ICP-OES, ICP-MS, MP-AES 제품용 소모품을 포함한 애질런트 원자 분광기 소모품의 전체 목록을 살펴보시려면 다음 사이트를 방문하시기 바랍니다. [www.agilent.com/chem/supplies](http://www.agilent.com/chem/supplies)

다음 사이트에서 주문하십시오.

[www.agilent.com/chem/wheretobuy](http://www.agilent.com/chem/wheretobuy)

이 정보는 사전 고지 없이 변경될 수 있습니다.

© Agilent Technologies, Inc. 2019  
2019년 5월 1일, 한국에서 인쇄  
5990-9475KO

서울시 용산구 한남대로 98, 일신빌딩 4층 우)04418  
한국애질런트테크놀로지스(주) 생명과학/화학분석 사업부  
고객지원센터 080-004-5090 [www.agilent.co.kr](http://www.agilent.co.kr)