

Agilent ChemStation Compatibility Matrix

	A.03.0x	A.04.0x	A.05.0x	A.06.0x	A.07.0x	A.08.0x	A.09.0x	A.10.0x	B.01.01	B.01.03	B.01.03 SR1	B.02.01
Win 3.1 / WFW 3.11	Yes	Yes	No	No	No	No	No	No	No	No	No	No
Win 95	No	Yes	Yes	Yes	No	No	No	No	No	No	No	No
Win NT 4.0	No	No	Yes	Yes	Yes	Yes	Yes	No	No	No	No	No
SP3	No	No	Yes	Yes	No	No	No	No	No	No	No	No
SP4	No	No	No	Yes	Yes	No	No	No	No	No	No	No
SP5	No	No	No	No	No	No	No	No	No	No	No	No
SP6a	No	No	No	No	Yes	Yes	Yes	No	No	No	No	No
Win 2000 Pro	No	No	No	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes
SP2	No	No	No	No	No	No	Yes	Yes	No	No	No	No
SP3	No	No	No	No	No	No	A.09.03 Only	Yes	No	No	No	No
SP4	No	No	No	No	No	No	A.09.03 Only	Yes	Yes	Yes	Yes	Yes
Win XP Pro SP1A	No	No	No	No	No	No	No	Yes	Yes	No	No	No
Win XP Pro SP2	No	No	No	No	No	No	No	Yes	Yes	Yes	Yes	Yes
LAN Communication	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
82335 GPIB Card	Yes	Yes	Yes	Yes	No	No	No	No	No	No	No	No
82341 GPIB Card	No	No	Yes	Yes	Yes	Yes	Yes	Yes	No	No	No	No
82350A PCI GPIB Card	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
82350B PCI GPIB Card	No	No	No	No	No	No	A.09.03 Only	Yes	Yes	Yes	Yes	Yes

	A.03.0x	A.04.0x	A.05.0x	A.06.0x	A.07.0x	A.08.0x	A.09.0x	A.10.0x	B.01.01	B.01.03	B.01.03 SR1	B.02.01
82357A USB GPIB Interface	No	No	No	No	No	No	A.09.03 Only	Yes	No	Yes - No 1100, A/D or GC support	Yes - No 1100, A/D or GC support	Yes - No 1100, A/D or GC support
GPIB I/O Libraries Rev: F.01.02	No	No	No	Yes	No	No	No	No	No	No	No	No
GPIB I/O Libraries Rev: F.01.01	No	No	No	Yes	No	No	No	No	No	No	No	No
GPIB I/O Libraries Rev: G.02.02.01	No	No	No	No	Yes	No	No	No	No	No	No	No
GPIB I/O Libraries Rev: H.01.02.00	No	No	No	Yes	Yes	Yes	Yes	Yes	No	No	No	No
GPIB I/O Libraries Rev: L.02.01	No	No	No	No	No	No	W2K Only	Yes	W2K Only	Yes	No	No
GPIB I/O Libraries Rev: M.01.01	No	No	No	No	No	No	No	No	Yes, GC only tested	Yes	Yes	Yes
35900C External A/D	Yes	Yes	Yes	Yes	Yes	Yes	No	No	No	No	No	No
35900D PC A/D Card	Yes	Yes	Yes	Yes	Yes	Yes	No	No	No	No	No	No
35900E External A/D	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes, LAN only	Yes, LAN only	Yes, LAN only	Yes, LAN only
5890 Series A	Yes	Yes	No	No	No	No	No	No	No	No	No	No

	A.03.0x	A.04.0x	A.05.0x	A.06.0x	A.07.0x	A.08.0x	A.09.0x	A.10.0x	B.01.01	B.01.03	B.01.03 SR1	B.02.01
5890 Series II and (+)	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
6850 (LAN Only)	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
6890	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
6890 Enhancements	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
6890N	No	No	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes
7673A (18594A)	Yes	Yes	No	No	No	No	No	No	No	No	No	No
7673B (18594B)	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	No	No	No
7673C (G1512A)	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
7683	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
7683 (variable speed)	No	No	No	Yes	Yes	Yes	Yes	Yes	NO	Yes	Yes	Yes
19405A (SECM)	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	NO	NO	NO	NO
19405B (SECM)	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	NO	NO	NO	NO
1090 LC System	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
1046 FLD	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes	Yes
1049 ECD	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes	Yes
1050 LC Modular System	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	No	No	No
system	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes, LAN only	Yes, LAN only	Yes, LAN only	Yes, LAN only

	A.03.0x	A.04.0x	A.05.0x	A.06.0x	A.07.0x	A.08.0x	A.09.0x	A.10.0x	B.01.01	B.01.03	B.01.03 SR1	B.02.01
1200 LC modular	see SN 1200-012: 1200 series vs 1100 series Overviw	see SN 1200-012: 1200 series vs 1100 series Overviw	see SN 1200-012: 1200 series vs 1100 series Overviw	see SN 1200-012: 1200 series vs 1100 series Overviw	see SN 1200-012: 1200 series vs 1100 series Overviw	see SN 1200-012: 1200 series vs 1100 series Overviw	see SN 1200-012: 1200 series vs 1100 series Overviw	see SN 1200-012: 1200 series vs 1100 series Overviw	see SN 1200-012: 1200 series vs 1100 series Overviw	see SN 1200-012: 1200 series vs 1100 series Overviw	Yes, requires FW 6.2	Yes
G1601A CE System	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes	Yes	Yes
G2201 CE/MS									No	Yes	Yes	Yes
G1946A LC/MS									No	No	No	No
G1946B LC/MS									Yes	Yes	Yes	Yes
G1946C LC/MS									Yes	Yes	Yes	Yes
G1946D LC/MS									Yes	Yes	Yes	Yes
G1956A LC/MS									Yes	Yes	Yes	Yes
G1956A LC/MS									Yes	Yes	Yes	Yes
G2182AA GPC	No	No	No	Yes	Yes	Yes	Yes	Yes	No	No	No	Yes
G1115AA UV-VIS	No	No	Yes	Yes	Yes	Yes	No	No	No	No	No	No
G1116AA UV-VIS	No	No	Yes	Yes	Yes	Yes	No	No	No	No	No	No
G1117AA UV-VIS	No	No	Yes	Yes	Yes	Yes	No	No	No	No	No	No
G1118AA UV-VIS	No	No	Yes	Yes	Yes	Yes	No	No	No	No	No	No