

Support \ Compatibility Matrix

GC/LC/CE/AD/LCMS ChemStation Revisions A.03.0x→A.10.0x

	A.03.0x	A.04.0x	A.05.0x	A.06.0x	A.07.0x	A.08.0x	A.09.0x	A.10.0x
Windows® 3.1 / WFW® 3.11	Yes	Yes	No	No	No	No	No	No
Windows®95	No	Yes	Yes	Yes	No	No	No	No
Windows® NT Workstation 4.0	No	No	Yes	Yes	Yes	Yes	Yes	No
Service Pack 3	No	No	Yes	Yes	No	No	No	No
Service Pack 4	No	No	No	Yes	Yes	No	No	No
Service Pack 5	No	No	No	No	No	No	No	No
Service Pack 6a	No	No	No	No	Yes	Yes	Yes	No
Windows® 98	No	No	No	No	No	No	No	No
Windows® ME	No	No	No	No	No	No	No	No
Windows®2000 Professional	No	No	No	No	No	Yes(12)	Yes	Yes
Service Pack 2	No	No	No	No	No	Yes(12)	Yes(A.09.01)	Yes
Service Pack 3	No	No	No	No	No	No	Yes(A.09.03)	Yes
Service Pack 4	No	No	No	No	No	No	No	Yes
Windows® XP	No	No	No	No	No	No	No	Yes
Service Pack 1a	No	No	No	No	No	No	No	Yes
LAN Communication	No	No	No	Yes	Yes	Yes	Yes	Yes
82335 GPIB Card (Windows® 3.1x and Windows® 95 only)	Yes	Yes	Yes	Yes	No	No	No	No
82341 GPIB Card (Not on Windows® 3.1x)	No	No	Yes	Yes	Yes	Yes	Yes	No
82350A PCI GPIB Card (Only on Windows® NT/2000/XP)	No	No	No	Yes(11)	Yes	Yes H.01.02	Yes J.02.01	Yes L.02.01
82350B PCI GPIB Card	No	No	No	No	No	No	Yes L.02.01	Yes L.02.01
GPIB I/O Libraries Rev. F.01.02	No	No	Yes	Yes(8)	No	No	No	No
GPIB I/O Libraries Rev. G.02.02	No	No	No	No	Yes	No	No	No
GPIB I/O Libraries Rev. H.01.02	No	No	No	Yes(9)	Yes(9)	Yes(9)	No	No
GPIB I/O Libraries Rev. J.02.01	No	No	No	No	No	Yes(14)	Yes	No
GPIB I/O Libraries Rev. L.02.01	No	No	No	No	No	No	Yes	Yes
1946A LC/MSD (GPIB only)	No	No	Yes	Yes	Yes	Yes(15)	No	No
1946B/C/D LC/MSD (LAN only)	No	No	No	No	Yes	Yes	Yes	Yes(16)
1956A/B LC/MSD (LAN only)	No	No	No	No	No	No	Yes	Yes
35900C External A/D	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No
35900D PC A/D Card	Yes	Yes	Yes(1)	Yes	Yes	No	No	No
35900E External A/D	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes

This document is believed to be accurate and up-to-date. However, Agilent Technologies, Inc. cannot assume responsibility for the use of this material. The information contained herein is intended for use by informed individuals who can and must determine its fitness for their purpose. Windows®, Windows® 95, Windows® NT, Windows® 2000 Pro & ME, Windows® XP are registered trademarks of Microsoft Corp.

	A.03.0x	A.04.0x	A.05.0x	A.06.0x	A.07.0x	A.08.0x	A.09.0x	A.10.0x
5890 Series A	Yes	Yes	No	No	No	No	No	No
5890 Series II and (+)	Yes	Yes	Yes	Yes	Yes	Yes	Yes(17)	Yes(17)
6850N (LAN only)	No	No	No	Yes(10)	Yes	Yes	Yes	Yes
6890A	Yes(7)	Yes	Yes	Yes	Yes	Yes	Yes(18)	Yes(18)
6890A Enhancements	No	Yes(2)	Yes	Yes	Yes	Yes	Yes(18)	Yes(18)
6890N (LAN only)	No	No	No	No	No	Yes	Yes	Yes
7673A (18594A)	Yes	Yes	No	No	No	No	No	No
7673B (18594B /G1512A)	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
7683	No	No	Yes(1)	Yes	Yes	Yes	Yes	Yes
7683 (variable speed)	No	No	No	Yes	Yes	Yes	Yes	Yes
19405A (SECM)	Yes	Yes(3)	Yes	Yes	Yes	Yes	Yes	Yes
19405B (SECM)	No	Yes(2)	Yes	Yes	Yes	Yes	Yes	Yes
1090 LC System	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
1050 LC Modular System	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
1100 LC Modular System	No	Yes(4)	Yes(4)	Yes(4)	Yes(4)	Yes(4)	Yes(4)	Yes(4)
G1601A CE System	Yes	Yes(5)	Yes	Yes	Yes	Yes	Yes	Yes
G2182AA GPC	No	No	No	Yes(6)	Yes	Yes	Yes	Yes

- (1) Supported on revisions A.05.02 and above
- (2) Supported on revisions A.04.02 and above
- (3) 19405A (ONLY) Supported on revision A.04.01
19405B (ONLY) Supported on revision A.04.02
- (4) G1321A FLD supported on rev $\geq$ A.06.01
G1362A RID supported on rev $\geq$ A.06.04
G1365A MWD supported on rev $\geq$ A.06.04
G1315B DAD supported on rev $\geq$ A.08.01
G1367A/77A/78A WPS supported on rev $\geq$ A.08.03
G1361/91A PREP PUMP supported on rev $\geq$ A.08.03
G1376A CAP PUMP supported on rev $\geq$ A.08.03
G1389A THM MICRO ALS supported on rev $\geq$ A.08.03
G2260A PREP ALS supported on rev $\geq$ A.09.01
G1364A FC supported on rev $\geq$ A.09.01
G1156A VALVES supported on rev $\geq$ A.09.03
G2226 NANO PUMP supported on rev $\geq$ A.09.03
G2257A SAMPLE CAPACITY supported on rev $\geq$ A.09.03
G1364B\C\D FC supported on rev $\geq$ A.10.01
- (5) CEC-mode supported on revisions A.04.0x and above
- (6) Supported on revisions A.06.03 and above
- (7) Supported on revisions A.03.02 and above
- (8) Supported on revisions A.06.01 /A.06.02 (ChemStation CD)
- (9) Requires H.01.03.71 Beta Patch for PCI style 82350A card (on ChemStation CD in \Support\GPIB folder)
- (10) Supported on revisions A.06.03 and above.
- (11) 82350A Requires GPIB drivers from the A.07.0x or A.08.0X CD-ROM (\GPIB folder). The 82350B tested using GPIB driver revision L.02.01 only.

Bonnie

This document is believed to be accurate and up-to-date. However, Agilent Technologies, Inc. cannot assume responsibility for the use of this material. The information contained herein is intended for use by informed individuals who can and must determine its fitness for their purpose. Windows®, Windows® 95, Windows® NT, Windows® 2000 Pro & ME, Windows® XP are registered trademarks of Microsoft Corp.

- (12) Microsoft® Windows 2000 Professional supported on ChemStation revision A.08.03 and higher.
(14) Use J.02.00 only with A.08.03 or A.08.04 on Windows® 2000 Professional, for NT use H.01.02 see (10)
(15) LC/MSD did not support any HPIB/GPIB instruments after revision A.08.01
(16) Only the 1946C and 1946D LAN MSD's are supported on revision A.10.0x
*Special Note pertains to all LC/MSD systems – When Agilent ChemStation is controlling any LC/MSD, no additional instruments on the PC are supported.
(17) The 5890A Series II GPIB Instrument requires latest firmware upgrade p/n 05890-80325 - Eprom
(18) The 6890A GPIB Instrument requires latest firmware upgrade p/n G1530-61706 - ROM Rev. A.03.08

*Note on older IO Libraries – Visit <http://adn.tm.agilent.com/> to obtain most of the older revisions.

This document is believed to be accurate and up-to-date. However, Agilent Technologies, Inc. cannot assume responsibility for the use of this material.
The information contained herein is intended for use by informed individuals who can and must determine its fitness for their purpose.
Windows®, Windows® 95, Windows® NT, Windows® 2000 Pro & ME, Windows® XP are registered trademarks of Microsoft Corp.