

## Learn convenient ChromPilot instrument control (P/N 899-0032)

### Description

This software training course is designed for present users of the WinGPC UniChrom ChromPilot and for future clients interested in implementing a convenient and cost effective infrastructure. It offers lectures and practical PC sessions given by an experienced polymer chemist with extensive knowledge about chromatography and WinGPC UniChrom.

Participants should be familiar with the basic WinGPC features and functions. Previous attendance of the WinGPC UniChrom software training course (P/N 899-0029) is recommended, but not required.

After successful participation each attendee should be able

- to administrate the ChromPilot
- to program sample sequences efficiently
- to set intelligent and money saving end actions



### Program

<b>08.50</b>	<b>Login and Communication's Check</b>
<b>09.00</b>	<b>Welcome, speaker introduction, explanation of training tools</b>
<b>09.15</b>	<b>Introduction to WinGPC UniChrom with ChromPilot</b> WinGPC UniChrom philosophy Comparison controlled/uncontrolled systems Supported systems
<b>09.30</b>	<b>ChromPilot installation and setup</b> PC requirements and communication ports ChromPilot Connection Manager Add and manage systems Verify communication
<b>10.00</b>	<b>Advanced device control</b> ChromPilot Instrument Manager Manage instruments and single modules, status information Recommended module settings Module end actions
<b>11.15</b>	<b>Practical Session part I</b>
<b>12.00</b>	<b>Lunch</b>
<b>13.00</b>	<b>ChromPilot sample sequences</b> Sample and sequence manager Import/export options, priority samples Start/pause/stop, stop with end actions
<b>14.30</b>	<b>Practical Session part II</b>

## Learn convenient ChromPilot instrument control (P/N 899-0032)

---

**15.15 Documentation**  
Instrument log books and audit trails  
Print sequence lists  
Additional module information from WinGPC with ReportDesigner

**15.45 Practical Session part III: Troubleshooting**

**Approx. 16.00 Course review and end of training course**

### Registration and organization

PSS Polymer Standards Service GmbH  
Training Academy  
P.O. Box 3368  
55023 Mainz, Germany  
Phone: 0049- 6131-96239-30; Fax: 0049- 6131-96239-11  
[info@pss-polymer.com](mailto:info@pss-polymer.com)  
<https://pss-polymer.com/training-events/event-list.html>

### How does it work?

- After receiving your registration, we will send you the registration confirmation.
- We will contact you to get to know you and to identify your course goals. We would be happy to help you check the suitability of your PC / laptop.
- Approx. one week before the course, you will receive a package from PSS with the printed course documents.
- One day before the event you will receive an e-mail with the access link.
- The training itself is interactive, you have time and opportunity to ask questions during the presentations and practical sessions. Please make sure that you have a microphone so that you can talk to us.
- You will receive your certificate of attendance one day after the training.

### Conditions of participation

To participate you need a PC/laptop with keyboard and mouse, internet access and audio functions. A headset is recommended. We recommend a quiet room for attendance. Attendance from home office is possible.

The course is conducted using a software training tool.

[Test the suitability of your system here](#)

A separate, additional access to a WinGPC installation for the practical sessions is provided.

PSS TeamViewer licenses are used here.

[Test the suitability of your system here](#)

**Participation fee** EURO 725,-

Includes participation (lectures and practical session), training documents and technical setup. Training documents are sent in printed form in advance. Your certificate of participation will be sent electronically.

### Cancellation policy

We ask for your understanding that if you cancel up to 2 weeks before the start of the course, a cancellation fee of 50% of the participation fee will be due. If you cancel at a later date, the participation fee can no longer be reimbursed, but a substitute participant can be provided.

DE44693546  
5994-6198EN

2

