



Fast Track to Lab Automation

Increasing the pace of discovery

Join Agilent and HighRes Biosolutions for an exclusive event at the F1 Arcade in Boston Seaport

Experience an exciting evening of discovery and high-speed fun as we explore the latest advances accelerating scientific workflows. Network with automation experts, enjoy immersive F1 racing experiences, and connect with peers driving the future of discovery. Drinks, bites, and fast-paced inspiration included!

Space is limited. Register today!

© Agilent Technologies, Inc. 2025

Event Details

Date: Tuesday February 10th, 2026

Time: 06:00 p.m. - 09:00 p.m.

Location:

F1 Arcade Boston Seaport

Conveniently located near the Convention Center
87 Pier 4 Blvd.

Boston, MA 02210

